



100 Majestic Way, Bangor, PA 18013 / www.biospectra.us

Effective Date:	23-Jun-2020	23-Jun-2023	: Date of Next Review																				
Initiated By:	Yencho, Amy	3.1	: Supersedes																				
Reason for Revision:	to add additional tests for GT5221 code																						
Approval:	<table> <thead> <tr> <th>Approvers</th><th>Time</th><th>Group</th><th>Name</th></tr> </thead> <tbody> <tr> <td>Date</td><td></td><td></td><td></td></tr> <tr> <td>23-Jun-2020</td><td>09:01:11 AM</td><td>SNR MGMT</td><td>Yencho, Amy</td></tr> <tr> <td>23-Jun-2020</td><td>09:04:08 AM</td><td>EDITOR</td><td>Saam, Amy</td></tr> <tr> <td>23-Jun-2020</td><td>09:50:34 AM</td><td>QUALITY</td><td>Miller, Jenna</td></tr> </tbody> </table>			Approvers	Time	Group	Name	Date				23-Jun-2020	09:01:11 AM	SNR MGMT	Yencho, Amy	23-Jun-2020	09:04:08 AM	EDITOR	Saam, Amy	23-Jun-2020	09:50:34 AM	QUALITY	Miller, Jenna
Approvers	Time	Group	Name																				
Date																							
23-Jun-2020	09:01:11 AM	SNR MGMT	Yencho, Amy																				
23-Jun-2020	09:04:08 AM	EDITOR	Saam, Amy																				
23-Jun-2020	09:50:34 AM	QUALITY	Miller, Jenna																				

GUANIDINE THIOCYANATE TESTING METHODS

VIEW ONLY

Printed On:	19-Oct-2020 10:49:54 AM	Ledergerber, Alissa	: Printed By
Print Expiration:	This copy is VIEW ONLY. It is not to be printed or distributed in any manner.		
Notice:	The information contained herein is the property of BioSpectra. The recipient is responsible for its safe-keeping and the prevention of unauthorized appropriation, use, disclosure and copying.		

TABLE OF CONTENTS

1. PURPOSE:.....	3
2. SCOPE:.....	3
3. RESPONSIBILITIES:	3
4. REFERENCES:	3
5. EQUIPMENT:	4
6. ANALYTICAL PROCEDURES:	4

VIEW ONLY

Printed On:	19-Oct-2020 10:49:54 AM	Ledergerber, Alissa	: Printed By
Print Expiration:	This copy is VIEW ONLY. It is not to be printed or distributed in any manner.		
Notice:	The information contained herein is the property of BioSpectra. The recipient is responsible for its safe-keeping and the prevention of unauthorized appropriation, use, disclosure and copying.		

1. PURPOSE:

- 1.1. To provide the Quality Control (QC) Laboratory personnel with a procedure for analyzing Guanidine Thiocyanate Raw Material, In-Process, Finished Goods, and Stability in the QC Laboratory at all BioSpectra facilities.

2. SCOPE:

- 2.1. Applies to the testing of Guanidine Thiocyanate Raw Material, In-Process, Finished Goods, and stability in the QC Laboratory at all BioSpectra facilities. Methods include testing for all types of Guanidine Thiocyanate sold by BioSpectra; only the specific tests required for the requested type must be tested.

3. RESPONSIBILITIES:

- 3.1. The Executive Director of Quality Control is responsible for the implementation, control, training and maintenance of this procedure.
- 3.2. The QC Analysts are responsible for compliance with the terms of this procedure. This includes notifying the appropriate personnel if any analyses fail to meet their respective specifications.
- 3.3. All QC Laboratory personnel are responsible for reviewing the appropriate SDS's prior to handling any chemicals used in this procedure.

4. REFERENCES:

- 4.1. *ACS, Reagent Chemicals*, current edition
- 4.2. [Balance SOP](#)
- 4.3. [Blue M Convection Oven Operation and Calibration SOP](#)
- 4.4. *Current EP/BP*
- 4.5. *Current JP*
- 4.6. *Current USP*
- 4.7. [DNase \(Endonuclease\) Assay](#)
- 4.8. [DNase \(Exonuclease\) Assay](#)
- 4.9. [Laboratory Chemicals](#)
- 4.10. [Laboratory Notebooks](#)
- 4.11. [Lambda 25 UV/VIS Operation and Calibration](#)
- 4.12. [Metrohm Titrando 907 Auto-Titrator SOP](#)
- 4.13. [MP50 Melting Range Operation and Calibration SOP](#)
- 4.14. [NexION 350X ICP-MS SOP](#)
- 4.15. [Protease Assay](#)
- 4.16. [Pipette SOP](#)
- 4.17. [Result Reporting](#)
- 4.18. [RNase \(Ribonuclease\) Assay](#)
- 4.19. [Spectrum Two UATR SOP](#)
- 4.20. [Standardization of Titrants](#)
- 4.21. [XL200 pH/mV/Conductivity Meter SOP](#)
- 4.22. [VWR Gravity Convection Oven and Calibration \(Model Number 414005-106\)](#)
- 4.23. Analytical Procedure for Gel Assays
- 4.24. Analytical Procedure for Protease Assay

Printed On:	19-Oct-2020 10:49:54 AM	Ledergerber, Alissa	: Printed By
Print Expiration:	This copy is VIEW ONLY. It is not to be printed or distributed in any manner.		
Notice:	The information contained herein is the property of BioSpectra. The recipient is responsible for its safe-keeping and the prevention of unauthorized appropriation, use, disclosure and copying.		

5. EQUIPMENT:

- 5.1. Lambda 25 UV/VIS Spectrophotometer
- 5.2. Analytical Balance
- 5.3. Perkin Elmer Spectrum Two UATR
- 5.4. Metrohm 907 Auto Titrator
- 5.5. Perkin Elmer NexION 350X ICP-MS
- 5.6. MP50 Melting Range Apparatus
- 5.7. XL200 pH/mV/Conductivity Meter SOP

6. ANALYTICAL PROCEDURES:**IN-PROCESS TESTING****6.1. MOTHER LIQUOR ASSAY MONITOR:**

- 6.1.1. Accurately weigh 0.36g of sample and transfer to a suitable beaker.
- 6.1.2. Add 10 mL of purified water, 10 mL of glacial acetic acid, 10 mL of 0.2% polyvinyl alcohol, and 100 mL methanol to the sample beaker.
- 6.1.3. Titrate with 0.1N AgNO₃ to a potentiometric endpoint utilizing the Metrohm Titrando 907.

$$\% \text{ Guanidine Thiocyanate} = \frac{\text{mL of AgNO}_3 \cdot N \text{ of AgNO}_3 \cdot 11.82}{\text{Sample Weight (g)}}$$

6.1.4. Alternate Manual Titration

- 6.1.4.1. Accurately weigh 0.36g of sample and transfer to a beaker and dissolve in 100mL of purified water.
- 6.1.4.2. Add 5mL of (10g/100mL) Ferric Ammonium Sulfate aqueous, 5mL of USP Dilute Nitric Acid, and 5mL of dibutyl phthalate.
- 6.1.4.3. Titrate to a white/colorless endpoint.

$$\% \text{ Guanidine Thiocyanate} = \frac{\text{mL of AgNO}_3 \cdot N \text{ of AgNO}_3 \cdot 11.82}{\text{Sample Weight (g)}}$$

6.2. ML ABSORBANCE ≤1.00 a.u. @ 280, 300, and 340 nm:

- 6.2.1. Prepare a 1:1 dilution of the mother liquor with purified water in a LOD vial or small beaker.
- 6.2.2. Swirl to homogenize.
- 6.2.3. Refer to the Lambda 25 UV/VIS Operation and Calibration procedure to determine the Absorbance of the sample.

6.3. WC ABSORBANCE 0.03 max. @ 340nm; 0.05max. @ 300; and 0.3max. @ 280 nm:

- 6.3.1. Accurately weigh 5.0 g of sample and accurately transfer the weighed sample to a graduated cylinder.
- 6.3.2. Dilute to 25 mL with purified water.
- 6.3.3. Dissolve completely.
- 6.3.4. Refer to the Lambda 25 UV/VIS Operation and Calibration procedure to determine the absorbance of the sample.

Printed On:	19-Oct-2020 10:49:54 AM	Ledergerber, Alissa	: Printed By
Print Expiration:	This copy is VIEW ONLY. It is not to be printed or distributed in any manner.		
Notice:	The information contained herein is the property of BioSpectra. The recipient is responsible for its safe-keeping and the prevention of unauthorized appropriation, use, disclosure and copying.		

6.4. Dry Crystal LOSS ON DRYING @ 105°C 0.3% max.:

- 6.4.1. Tare an LOD vial that has been previously dried for 30 minutes under the same conditions to be employed in the determination. Cool in desiccator for at least 15 minutes before weighing.
- 6.4.2. Transfer approximately 3 g of the sample to be tested to the LOD vial, and accurately weigh the LOD vial and contents. By gentle, sidewise shaking, distribute the sample as evenly as possible.
- 6.4.3. Place the LOD vial containing the sample into the oven.
- 6.4.4. Dry the sample at 105°C ± 2°C for 3 hours.
- 6.4.5. Allow to cool to room temperature in a desiccator for at least 15 minutes before weighing.
- 6.4.6. Calculate the loss according to the following calculation:

$$\%LOD = \frac{(Initial\ Sample\ Weight\ (g) + Vial(g)) - Final\ Vial\ Weight\ (g)}{Initial\ Sample\ Weight\ (g)} \times 100$$

FINISHED GOOD TESTING**6.1. ABSORBANCE (1.7M) Refer to Summary Sheet:**

- 6.1.1. Accurately weigh 5.0 g of sample.
- 6.1.2. Transfer accurately weighed sample to a graduated cylinder and Q.S. to 25 mL with purified water.
- 6.1.3. Dissolve completely.
- 6.1.4. Refer to Lambda 25 UV/VIS Operation and Calibration procedure to determine the Absorbance of the sample.

6.2. ABSORBANCE (70%) Refer to Summary Sheet:

- 6.2.1. Accurately weigh 17.5 g of sample.
- 6.2.2. Transfer accurately weighed sample to a graduated cylinder and Q.S. to 25 mL with purified water.
- 6.2.3. Dissolve completely.
- 6.2.4. Refer to Lambda 25 UV/VIS Operation and Calibration procedure to determine the Absorbance of the sample.

6.3. ACETONE TEST (20% W/W) Refer to Summary Sheet:

- 6.3.1. Weigh 2 g of sample in a beaker.
- 6.3.2. Add 8 g of acetone to sample.
- 6.3.3. Mix to dissolve completely.
- 6.3.4. Solution should be clear and free of particles to pass the test.

Printed On:	19-Oct-2020 10:49:54 AM	Ledergerber, Alissa	: Printed By
Print Expiration:	This copy is VIEW ONLY. It is not to be printed or distributed in any manner.		
Notice:	The information contained herein is the property of BioSpectra. The recipient is responsible for its safe-keeping and the prevention of unauthorized appropriation, use, disclosure and copying.		

6.4. APPEARANCE AND COLOR**White / Crystals:**

- 6.4.1. Sample Size:
 - 6.4.1.1. For Raw Material: Inspect the entire testing sample for appearance and color.
 - 6.4.1.2. For Finished Goods: Use a suitable amount of sample.
- 6.4.2. Place sample into a clean, dry glass beaker.
- 6.4.3. In an area with sufficient lighting, view the sample from all sides and gently sift through the crystals inspecting for nonconforming matter, color, and structure.
- 6.4.4. The sample should be white in color and characteristic of crystals.

6.5. ASSAY**Refer to Summary Sheet:**

- 6.5.1. Note: Raw Material may be analyzed as-is. Refer to product code analysis required.
- 6.5.2. As-Is: Accurately weigh 0.36g of sample. Transfer to a beaker.
- 6.5.3. Dried Basis: Accurately weigh 0.36g of sample, previously dried at $105 \pm 2^{\circ}\text{C}$ for 3 hours (utilize LOD sample, if available). Transfer to a beaker.
- 6.5.4. Add 10 mL of purified water, 10 mL of glacial acetic acid, 10 mL of 0.2% polyvinyl alcohol, and 100 mL methanol to the sample beaker.
- 6.5.5. Titrate with 0.1N AgNO_3 to a potentiometric endpoint utilizing the Metrohm Titrando 907.

$$\% \text{ Guanidine Thiocyanate} = \frac{(\text{mL of AgNO}_3)(N \text{ of AgNO}_3)(11.82 \text{ mg})}{\text{Sample Weight (g)}}$$

Printed On:	19-Oct-2020 10:49:54 AM	Ledergerber, Alissa	: Printed By
Print Expiration:	This copy is VIEW ONLY. It is not to be printed or distributed in any manner.		
Notice:	The information contained herein is the property of BioSpectra. The recipient is responsible for its safe-keeping and the prevention of unauthorized appropriation, use, disclosure and copying.		

6.5.1. Alternate Manual Titration

- 6.5.1.1. As-Is: Accurately weigh 0.36g of sample. Transfer to a beaker.
- 6.5.1.2. Dried Basis: Accurately weigh 0.36g of sample, previously dried at $105 \pm 2^\circ\text{C}$ for 3 hours (utilize LOD sample, if available). Transfer to a beaker.
- 6.5.1.3. Dissolve with 100mL of purified water.
- 6.5.1.4. Add 5mL of (10g/100mL) Ferric Ammonium Sulfate aqueous, 5mL of USP Dilute Nitric Acid, and 5mL of dibutyl phthalate.
- 6.5.1.5. Titrate to a white endpoint.

$$\% \text{ Guanidine Thiocyanate} = \frac{(\text{mL of AgNO}_3)(N \text{ of AgNO}_3) (11.82 \text{ mg})}{\text{Sample Weight (g)}}$$

6.6. CHLORIDE 5 ppm max.:

- 6.6.1. Accurately weigh 14.3 g of sample, transfer to a Nessler Color Comparison Tube, and Q.S. to 50 mL using purified water.
- 6.6.2. Prepare a standard by pipetting 101 μL of 0.02 N HCl into a Nessler Color Comparison Tube, and Q.S. to 50 mL using purified water.
- 6.6.3. To the sample and standard add 1 mL of concentrated nitric acid and 1 mL of 0.1N silver nitrate.
- 6.6.4. The sample should be less turbid than the standard.

6.7. ENZYME ACTIVITY None Detected:

- 6.7.1. RNase, DNase, and Protease performed as per procedures outlined in section 3. Analysis should be performed in the Analytical Procedure for Gel Assays and Analytical Procedure for Protease Assay Packets.

6.8. IDENTIFICATION (IR) Passes Test:

- 6.8.1. Note: Raw Material may be analyzed as-is. Refer to product code analysis required.
- 6.8.2. Follow Spectrum Two UATR SOP to perform IR analysis.

6.9. LOSS ON DRYING @ 105°C Refer to Summary Sheet:

- 6.9.1. Tare an LOD vial that has been previously dried for 30 minutes under the same conditions to be employed in the determination. Cool in desiccator for at least 15 minutes before weighing.
- 6.9.2. Transfer approximately 3 g of the sample to be tested to the LOD vial, and accurately weigh the LOD vial and contents. By gentle, sidewise shaking, distribute the sample as evenly as possible.
- 6.9.3. Place the LOD vial containing the sample into the oven.
- 6.9.4. Dry the sample at $105^\circ\text{C} \pm 2^\circ\text{C}$ for 3 hours.
- 6.9.5. Allow to cool to room temperature in a desiccator for at least 15 minutes before weighing.
- 6.9.6. Calculate the loss according to the following calculation:

$$\% \text{ LOD} = \frac{(\text{Initial Sample Weight (g)} + \text{Vial(g)}) - \text{Final Vial Weight (g)}}{\text{Initial Sample Weight (g)}} \times 100$$

Printed On:	19-Oct-2020 10:49:54 AM	Ledergerber, Alissa	: Printed By
Print Expiration:	This copy is VIEW ONLY. It is not to be printed or distributed in any manner.		
Notice:	The information contained herein is the property of BioSpectra. The recipient is responsible for its safe-keeping and the prevention of unauthorized appropriation, use, disclosure and copying.		

6.10. MELTING RANGE **Refer to Summary Sheet:**

- 6.10.1. Note: Raw Material may be analyzed as-is.
- 6.10.2. Follow MP50 Melting Range Operation and Calibration procedure to determine melting range.

6.11. pH OF A 5% SOLUTION @25° ±2°C **Refer to Summary Sheet:**

- 6.11.1. Accurately weigh 5.0 g of sample.
- 6.11.2. Dissolve the sample in 100mL of purified water. If necessary, utilize a stir plate and Teflon encapsulated magnetic stirring bar to achieve solubility.
- 6.11.3. Follow the appropriate SOP to calibrate and record the pH measurement of the solution at the 25 ± 2°C.

6.12. SOLUBILITY OF A 6M SOLUTION **Clear:**

- 6.12.1. Weigh 35.2 g of sample and transfer to a graduated cylinder.
- 6.12.2. Dissolve in and dilute to 50 mL with purified water.
- 6.12.3. Observe from all sides, under sufficient lighting.
- 6.12.4. Solution should be clear when compared to purified water.

6.13. SOLUBILITY OF A 35% SOLUTION **Clear:**

- 6.13.1. Weigh 17.5 g of sample and record weight and transfer to a graduated cylinder.
- 6.13.2. Dissolve and dilute to 50mL with purified water.
- 6.13.3. Observe from all sides under sufficient lighting.
- 6.13.4. Solution should be clear when compared to purified water.

6.14. SOLUBILITY (COLOUR) **Colorless:**

- 6.14.1. Prepare a 100mg/mL sample, dissolve and dilute with purified water.
- 6.14.2. Observe from all sides under sufficient lighting.
- 6.14.3. Solution must be colorless when compared to purified water to pass test.

6.15. SOLUBILITY (TURBIDITY) 100mg/mL H₂O **Clear:**

- 6.15.1. Prepare a 100mg/mL sample, dissolve and dilute with purified water.
- 6.15.2. Observe from all sides under sufficient lighting.
- 6.15.3. Solution must be clear when compared to purified water to pass test.

6.16. SULFATE **5 ppm max:**

- 6.16.1. Accurately weigh 2 g of sample.
- 6.16.2. Dissolve in ~5 mL of purified water.
- 6.16.3. Add 1 mL of HCl (1 in 20).
- 6.16.4. Dilute to 10 mL with purified water.
- 6.16.5. Prepare a standard using 1 mL of Sulfate Standard Solution (0.01mg/mL), 1 mL HCl (1 in 20), and dilute to 10 mL with purified water.
- 6.16.6. To the sample and standard add 1 mL of 12% barium chloride and allow to stand for 10 minutes utilizing a calibrated timer.
- 6.16.7. The sample should be less turbid than the standard.

6.17. TRACE METALS **Refer to Summary Sheet:**

- 6.17.1. Refer to NexION 350X ICP-MS SOP for sample preparation and analysis.

Printed On:	19-Oct-2020 10:49:54 AM	Ledergerber, Alissa	: Printed By
Print Expiration:	This copy is VIEW ONLY. It is not to be printed or distributed in any manner.		
Notice:	The information contained herein is the property of BioSpectra. The recipient is responsible for its safe-keeping and the prevention of unauthorized appropriation, use, disclosure and copying.		

6.18. TRACE METALS**Ba, K, Na:**

6.18.1. Refer to NexION 350X ICP-MS SOP for sample preparation and analysis.

6.18.2. Refer to summary sheet for specifications.

6.19. WATER INSOLUBLES**Refer to Summary Sheet:**

6.19.1. Accurately weigh 20 grams of sample utilizing an analytical balance.

6.19.2. Dissolve in 200 mL of purified water.

6.19.3. Heat to boiling and digest, covered, on a hot plate for 1 hour.

6.19.4. Prepare a Gooch filtering crucible and 10-15µm filter by drying at 105±2°C for 1 hour.
Allow to cool in ambient air at least 15 minutes and weigh on analytical balance.

6.19.5. Filter solution through conditioned filtering crucible and 10-15 µm filter. Rinse thoroughly with hot purified water.

6.19.6. Dry the crucible at 105±2°C for 1 hour.

6.19.7. Cool in ambient air for at least 15 minutes and reweigh.

6.19.8. Calculate the water insoluble content using the following calculation:

$$\% \text{ Water Insolubles} = \frac{\text{Initial Crucible Weight}_{(g)} - \text{Final Crucible Weight}_{(g)}}{\text{Sample Weight}_{(g)}} \cdot 100$$

VIEW ONLY

Printed On:	19-Oct-2020 10:49:54 AM	Ledergerber, Alissa	: Printed By
Print Expiration:	This copy is VIEW ONLY. It is not to be printed or distributed in any manner.		
Notice:	The information contained herein is the property of BioSpectra. The recipient is responsible for its safe-keeping and the prevention of unauthorized appropriation, use, disclosure and copying.		