



100 Majestic Way, Bangor, PA 18013 / www.biospectra.us

To Whom It May Concern,

INTRODUCTION: The following analyses are conducted for Sodium Chloride 5M Solution, product code NC3150, in accordance with the Sodium Chloride 5M Testing Methods DCN: 19-003065 v.1.0 and Certificate of Analysis DCN: 19-003007 v.1.1. Specific details for the procedures were also obtained from NexION 350X ICP-MS SOP DCN: 16-001923 v.2.0.

1. APPEARANCE AND COLOR Clear Colorless Solution:

- 1.1. Transfer 50mL of the sample in a clean, dry color Nessler tube.
- 1.2. In order to pass, test solution is a clear, colorless solution. Verify the solution appearance against a clear and colorless reference solution, such as purified water, and view against a color comparison plate with suitable lighting.

2. ASSAY 4.9-5.1M:

- 2.1. Standardize or perform a daily check of 0.1N AgNO₃ as per Standardization of Titrants.
- 2.2. Transfer 10mL of the sample solution into a 100mL volumetric flask and dilute to volume with purified water.
- 2.3. Accurately transfer 3mL of the diluted solution to a suitable titration vessel. Add 7mL of purified water, 10mL of glacial acetic acid, 75mL of methanol, and Eosin Y TS as the indicator.
- 2.4. While stirring, titrate with 0.1N Silver Nitrate TS to a pink endpoint. Each mL of 0.1N Silver Nitrate is equivalent to 5.844mg of Sodium Chloride.

$$M NaCl = \frac{(mL \text{ of titer})(N \text{ of titer})(10)}{3}$$

3. ENDOTOXIN ≤2.5 EU/mL:

- 3.1. Add 0.34mL of sample into a sterile capped conical tube or other suitable sterile vessel.
- 3.2. Dilute to 10mL with LAL reagent water.
 - 3.2.1. The sample preparation may be scaled to desired volume.
- 3.3. Enter specification of 2.5EU/mL.
- 3.4. Ensure all data is accurate before proceeding to load sample.
- 3.5. Load 25µL of sample in to each well and run test.
- 3.6. Ensure all suitability requirements meet specification.
- 3.7. Report result directly from instrument printout.

4. IDENTIFICATION CHLORIDE (USP<191>) Meets Requirements:

- 4.1. Transfer 5mL of sample solution to a beaker.
- 4.2. To the sample solution add 1mL of 0.1N silver nitrate.
- 4.3. Mix thoroughly.

Printed On:	31-Aug-2020 01:12:01 PM	Ledergerber, Alissa	: Printed By
Authored By:	Baun, Cassie	27-Aug-2020 02:54:14 PM	
Approval:	Hamelburg, Crystal	31-Aug-2020 11:49:47 AM	
Notice:	The information contained herein is the property of BioSpectra. The recipient is responsible for its safe-keeping and the prevention of unauthorized appropriation, use, disclosure and copying.		

- 4.4. A curdled, white precipitate is formed.
- 4.5. The precipitate is insoluble when 1 mL of concentrated Nitric acid is added but soluble when 6 mL of 6N Ammonium hydroxide is added.

5. IDENTIFICATION SODIUM (USP<191>) Meets Requirements:

- 5.1. Transfer 0.34mL of sample solution into a suitable beaker.
- 5.2. Add 2mL of purified water.
- 5.3. Add 2mL of 15% potassium carbonate and heat to a boil. No precipitate is formed.
- 5.4. Add 4mL of Potassium Pyroantimonate TS and heat to a boil.
- 5.5. Allow to cool in an ice bath, and if necessary, rub the inside of the test tube with a glass rod. A dense precipitate is formed to pass test.

6. MICROBIAL CONTENT TAMC ≤50 CFU/g; TYMC ≤150 CFU/g:

- 6.1. Package ~20g into a sterile container and send to an approved Outside Testing Facility. The analysis request form should include TAMC and TYMC. In order to pass the Total Aerobic Microbial Count must be ≤50CFU/g, the Total Yeast and Mold Count must be ≤150CFU/g.

7. TRACE METALS As, Cu, Fe, and Pb ≤2 ppm:

- 7.1. Refer to NexION 350X ICP-MS SOP for Instrument Set-Up and Use.
 - 7.1.1. Sample Preparation:
 - 7.1.2. General Notes:
 - 7.1.2.1. Before use, all plasticware that is not rated as “metal-free”, should first be rinsed with purified water, rinsed with 15% Nitric Acid, and then rinsed again with purified water. Plasticware rated as “metal-free” may be used as-is.
 - 7.1.2.2. Glass should be avoided as it has high potential for metal and mineral contaminations.
 - 7.1.2.3. Standard and sample solutions should be prepared in 50mL centrifuge tubes.
 - 7.1.3. 1% Nitric Acid
 - 7.1.3.1. Measure 14.5 mL of Trace Metal Grade Nitric Acid and transfer to a rinsed plastic 1000 mL volumetric flask. QS to 1000 mL with purified water.
 - 7.1.4. 15% Nitric Acid
 - 7.1.4.1. Dilute approximately 110 mL of Trace Metal Grade Nitric Acid to 500 mL with purified water.
 - 7.1.4.2. The solution is only used to rinse glassware and plasticware.
 - 7.1.5. BioSpectra Daily Method:
 - 7.1.5.1. Sample Solutions:
 - 7.1.5.1.1. Weigh 0.10g of sample on an analytical balance. Add 100 µL of Environmental Standard Mix 6 and QS to 50.0 g with 1% Nitric Acid.
 - 7.1.5.2. Standard Curve Preparation:
 - 7.1.5.2.1. 2 ppm Stock
 - 7.1.5.2.1.1. Weigh 1.00 g of Instrument Calibration Standard 2 and QS to 50.0 g with 1% Nitric Acid.

Printed On:	31-Aug-2020 01:12:01 PM	Ledergerber, Alissa	: Printed By
Authored By:	Baun, Cassie	27-Aug-2020 02:54:14 PM	
Approval:	Hamelburg, Crystal	31-Aug-2020 11:49:47 AM	
Notice:	The information contained herein is the property of BioSpectra. The recipient is responsible for its safe-keeping and the prevention of unauthorized appropriation, use, disclosure and copying.		

- 7.1.5.2.2. 100 ppb Stock
 - 7.1.5.2.2.1. Weigh 2.50 g of 2 ppm Stock and QS to 50.0 with 1% Nitric Acid.
- 7.1.5.2.3. Blank
 - 7.1.5.2.3.1. Pipette 100 µL of Environmental Standard Mix 6 into the centrifuge tube. QS to 50.0 g with 1% Nitric Acid.
- 7.1.5.2.4. 1 ppb Standard
 - 7.1.5.2.4.1. Pipette 0.50 mL of the 100 ppb Stock, add 100 µL of Environmental Standard Mix 6 and QS to 50.0 g with 1% Nitric Acid.
- 7.1.5.2.5. 2 ppb Standard (also used as a Continuing Check Verification Sample (CCV))
 - 7.1.5.2.5.1. Pipette 1.00 mL of the 100 ppb Stock, add 100 µL of Environmental Standard Mix 6 and QS to 50.0 g with 1% Nitric Acid.
- 7.1.5.2.6. 4 ppb Standard
 - 7.1.5.2.6.1. Pipette 2.00 mL of the 100 ppb Stock, add 100 µL of Environmental Standard Mix 6 and QS to 50.0 g with 1% Nitric Acid.
- 7.1.5.2.7. 6 ppb Standard
 - 7.1.5.2.7.1. Pipette 3.00 mL of the 100 ppb Stock, add 100 µL of Environmental Standard Mix 6 and QS to 50.0g with 1% Nitric Acid.

If there are any questions or concerns, please feel free to contact ra@biospectra.us.

Sincerely,



Cassie Baun
 Compliance Specialist
BioSpectra, Inc.
 100 Majestic Way
 Bangor, PA 18013

Printed On:	31-Aug-2020 01:12:01 PM	Ledergerber, Alissa	: Printed By
Authored By:	Baun, Cassie	27-Aug-2020 02:54:14 PM	
Approval:	Hamelburg, Crystal	31-Aug-2020 11:49:47 AM	
Notice:	The information contained herein is the property of BioSpectra. The recipient is responsible for its safe-keeping and the prevention of unauthorized appropriation, use, disclosure and copying.		