



100 Majestic Way, Bangor, PA 18013 / www.biospectra.us

SODIUM DECANOATE TESTING METHODS

The information contained herein is the confidential property of BioSpectra. The recipient is responsible for its safe-keeping and the prevention of unauthorized appropriation, use, disclosure and copying.

Page 1 of 11

Controlled Copy Number: 6, Controlled Copy Location: website placement
Printed By: ALI.TISON, on 22 Sep 2022;

TABLE OF CONTENTS

1. PURPOSE:.....	3
2. SCOPE:.....	3
3. RESPONSIBILITIES:	3
4. REFERENCES:.....	3
5. EQUIPMENT:.....	3
6. ANALYTICAL PROCEDURES:.....	4
TABLE 1: OVEN TEMPERATURE GRADIENT.....	8

1. PURPOSE:

- 1.1. To provide the Quality Control (QC) Laboratory personnel with procedures for testing Sodium Decanoate Raw Material, In-Process and Finished Goods at the Bangor, PA facility

2. SCOPE:

- 2.1. Applies to the testing of Sodium Decanoate Raw Material, In-Process and Finished Goods in the QC Laboratory at Biospectra. Methods include testing for all grades of Sodium Decanoate sold by BioSpectra; only the specific tests required for the requested grade must be tested.

3. RESPONSIBILITIES:

- 3.1. The Quality Control Manager is responsible for the control, training, maintenance and implementation of this procedure.
- 3.2. The QC Analysts are responsible for compliance with the terms of this procedure. This includes notifying the appropriate personnel if any analyses fail to meet their respective specifications.
- 3.3. All QC Laboratory personnel are responsible for reviewing the appropriate SDS's prior to handling any chemicals used in this procedure.

4. REFERENCES:

- 4.1. BSI-ATM-0082, Analytical Method of Analysis: Sodium in Sodium Decanoate by ICP-OES
- 4.2. BSI-ATM-0083, Analytical Method of Analysis: Elemental Impurities in Sodium Decanoate via ICP-MS
- 4.3. BSI-RPT-0441, Analytical Method Validation Report: Sodium Decanoate Organic Impurities Quantitative
- 4.4. BSI-SOP-0019, Result Reporting
- 4.5. BSI-SDS-0021, Sodium Decanoate SDS
- 4.6. BSI-SOP-0098, Balance SOP
- 4.7. BSI-SOP-0126, Laboratory Notebooks
- 4.8. BSI-SOP-0133, Blue M Convection Oven Operation and Calibration SOP
- 4.9. BSI-SOP-0143, Metrohm Titrando 907 Auto-Titrator SOP
- 4.10. BSI-SOP-0144, Metrohm 914 pH Conductometer Operation and Calibration
- 4.11. BSI-SOP-0242, Bangor Portable Turbidimeter and Calibration
- 4.12. BSI-SOP-0244, VWR Gravity Convection Oven and Calibration (Model Number 414005-106)
- 4.13. BSI-SOP-0254, Spectrum Two UATR SOP
- 4.14. BSI-SOP-0255, XL200 pH/mV/Conductivity Meter SOP
- 4.15. BSI-SOP-0316, Shimadzu QP2010S GC/MS SOP

5. EQUIPMENT:

- 5.1. Analytical Balance
- 5.2. Calibrated Oven
- 5.3. Metrohm 907 Autotitrator SOP (pH)
- 5.4. Metrohm 914 pH Conductometer Operation and Calibration
- 5.5. Spectrum Two UATR
- 5.6. XL200 pH/mV/Conductivity Meter
- 5.7. Hach Portable Turbidimeter Model 2100Q
- 5.8. Shimadzu QP2010S GC/MS

The information contained herein is the confidential property of BioSpectra. The recipient is responsible for its safe-keeping and the prevention of unauthorized appropriation, use, disclosure and copying.

6. ANALYTICAL PROCEDURES:**6.1. PRE-FILTRATION IN-PROCESS ANALYSIS MONITOR ONLY:**

6.1.1. Perform Visual Inspection and Quantitative NTU result. For NTU result, follow the SOP for Bangor Portable Turbidimeter and Calibration, DCN: BSI-SOP-0242.

6.2. POST-FILTRATION IN-PROCESS SAMPLES CLEAR, TO SLIGHTLY TURBID:

6.2.1. Perform Visual Inspection and Quantitative NTU result. For NTU result, follow the SOP for Bangor Portable Turbidimeter and Calibration, DCN: BSI-SOP-0242.

6.3. DRY CRYSTAL KARL FISCHER 1.5 – 3.0%:

6.3.1. Weigh ~1g of sample into a glass weighing spoon and tare the balance.

6.3.2. Transfer the sample to the Karl Fischer vessel by removing the rubber septum and adding the sample into the titration vessel.

6.3.2.1. Do not leave the rubber septum open for long periods of time as this will allow moisture to enter the titration vessel.

6.3.3. Return the weighing spoon to the balance, making sure not to lose any sample that was left behind. Once the weight stabilizes, Record the sample weight and transfer to instrument.

6.3.4. Check to make sure there is no residual sample stuck to the sides of the titration vessel.

6.3.4.1. If any sample is stuck to the side, stop the stir bead and swirl the Karl Fischer vessel to rinse the sides.

6.3.5. Once the method begins, check to ensure the sample is fully dissolved before the titration begins (i.e. before the stir command completes).

6.3.6. The moisture content will then be determined by the Metrohm Auto Titrando 907.

$$\% \text{ Moisture} = \frac{(mL \text{ of Composite 5}) \left(\frac{mg}{mL} \text{ of Composite 5} \right) (0.1)}{\text{Sample Weight (g)}}$$

6.4. APPEARANCE White to off-white / Powder:

6.4.1. Place a suitable amount of the sample in a clean, dry glass beaker.

6.4.2. In an area with sufficient lighting, view the sample from all sides.

6.4.3. The sample should be white to off-white in color and characteristic of powder.

6.4.4. Report any foreign matter or nonconformity of sample immediately to the appropriate personnel.

6.5. ASSAY (Dried Basis) Refer to Summary Sheet:

6.5.1. **Note: Raw Material Assay analysis requires the Loss on drying analysis and result. Perform LOD on the Raw material sample; there is no raw material specification requirement.**

6.5.2. **Standardization of 0.1N Perchloric Acid:**

6.5.2.1. Note: Before each titration, record the temperature of the 0.1N Perchloric acid.

6.5.2.2. In a well-ventilated area, accurately weigh about 0.7g of previously dried Potassium hydrogen phthalate and transfer into a suitable beaker (150mL size is ideal), record weight immediately after stabilization.

6.5.2.3. Ensure stir bar is in the beaker before the addition of the glacial acetic acid to avoid splashing.

6.5.2.4. Dissolve in 50mL of Glacial Acetic Acid.

6.5.2.5. Add 0.1mL of a 1% crystal violet solution prepared by dissolving 100mg of crystal violet in 10mL of glacial acetic acid.

The information contained herein is the confidential property of BioSpectra. The recipient is responsible for its safe-keeping and the prevention of unauthorized appropriation, use, disclosure and copying.

- 6.5.2.6. Prepare a burette with 0.1N Perchloric Acid in glacial acetic acid and titrate with 0.1N Perchloric Acid to a Blue-Green end point.
- 6.5.2.7. Note: The crystal violet solution will turn blue before blue green, the end point is immediately AFTER the blue end point when the first appearance of a green color is indicated.
- 6.5.2.8. Calculate normality of the titrant utilizing the following equation:

$$N \text{ Perchloric Acid} = \frac{g \text{ KHP} \left(\frac{\text{Reference Material Assay}}{100} \right)}{0.20423 \times (\text{Perchloric Acid Volume (mL)} - \text{Blank (mL)})}$$

- 6.5.2.9. Perform standardization in triplicate with a single blank.
- 6.5.3. Sample preparation:
 - 6.5.3.1. Note: Due to the nature of non-aqueous titrations minimize sample and titrant exposure to the environment (i.e. prepare samples one at a time and only immediately before titration).
 - 6.5.3.2. Note: Before each titration, record the temperature of the 0.1N Perchloric acid.
 - 6.5.3.3. Accurately weigh 0.5g of Sodium Decanoate and transfer to a suitable beaker.
 - 6.5.3.4. Dissolve in ~50mL of Glacial Acetic Acid.
 - 6.5.3.5. Add 0.1mL of a 1% crystal violet solution prepared by dissolving 100mg of crystal violet in 10mL of glacial acetic acid.
 - 6.5.3.6. Prepare a burette with 0.1N Perchloric Acid in glacial acetic acid and titrate with 0.1N Perchloric Acid to a Blue-Green end point.
 - 6.5.3.7. Note: The crystal violet solution will turn blue before blue green; the end point is immediately AFTER the blue end point when the first appearance of a green color is indicated.
 - 6.5.3.8. To calculate the result: Reference the LOD result from section 6.9, calculate the temperature correction refer to section 6.5.4. and apply the calculation in section 6.5.5.
- 6.5.4. Temperature Correction:
 - 6.5.4.1. 0.1N Perchloric acid in glacial acetic acid has a high coefficient of expansion. Due to the high coefficient of expansion of the titrant and inherent variability of temperature, care should be taken to correct for the difference in temperature from when the titration is carried out (T_i) and from when the titrant was standardized (T_s).
 - 6.5.4.1.1. To Calculate the Temperature Correction Factor Coefficient (C_f) use the following equation:

$$C_f = [1 + a(T_s - T_i)]$$
 - 6.5.4.2. Where:
 - 6.5.4.2.1. α = Thermal expansion coefficient of titrant (Acetic Acid= 1.07×10^{-3})
 - 6.5.4.2.1.1. T_i = temperature at sample titration
 - 6.5.4.2.1.2. T_s = temperature at titrant standardization.
- 6.5.5. Assay Result Calculation:
 - 6.5.5.1. Calculate the result using the following equation:

$$\text{Result (\%)} = \frac{(C_f)(EP1 - EPb)(N)(19.425)}{W} \times \left(\frac{100}{100 - LOD} \right)$$

- 6.5.5.2. Where:
- 6.5.5.2.1. C_f = Temperature Correction Factor Coefficient (Refer to 6.4.4.2.1.)
 - 6.5.5.2.2. EP1 = Sample Titration End Point (mL)
 - 6.5.5.2.3. EPb = Blank Titration End Point (mL)
 - 6.5.5.2.4. N = Normality of Titrant as Standardized
 - 6.5.5.2.5. W = Sample Weight (g)
 - 6.5.5.2.6. LOD = Loss on Drying (%)

6.6. ELEMENTAL IMPURITIES IRON, AND SODIUM QUANTIFICATION Refer to Summary Sheet:

- 6.6.1. Primary Method: Refer to Analytical Method of Analysis: Sodium in Sodium Decanoate by ICP-OES, BSI-ATM-0082, and Analytical Method of Analysis: Elemental Impurities in Sodium Decanoate via ICP-MS, BSI-ATM-0083.
- 6.6.2. Alternate Method: Package and send no less than 10 grams of sample to an approved outside testing facility for Elemental Impurities, Iron, and Sodium Quantification.

6.7. IDENTIFICATION (IR) Refer to Summary Sheet:

- 6.7.1. For UATR analysis, follow Spectrum Two UATR SOP.

6.8. IDENTIFICATION (GC) Refer to Summary Sheet:

- 6.8.1. Refer to Single Impurities and Total Impurities Analysis.

6.9. LOSS ON DRYING Refer to Summary Sheet:

- 6.9.1. Tare an LOD vial that has been previously dried for 30 minutes under the same conditions to be employed in the determination.
- 6.9.2. Allow to come to room temperature in a desiccator for at least 15 minutes before weighing.
- 6.9.3. Transfer approximately 2 g of the sample to be tested to the LOD vial, and accurately weigh the vial and contents. By gentle, sidewise shaking, distribute the sample as evenly as possible in the weighing bottle.
- 6.9.4. Place the LOD vial containing the sample into the oven.
- 6.9.5. Dry the sample at 105°C +/- 2°C for 3 hours.
- 6.9.6. Allow to come to room temperature in a desiccator for at least 15 minutes before weighing.
- 6.9.7. Calculate result using the equation below:

$$\%LOD = \frac{Initial\ Sample\ Weight(g) - Final\ Sample\ Weight(g)}{Initial\ Sample\ Weight(g)} \times 100$$

6.10. pH 10% **9.0-11.0:**

- 6.10.1. Accurately weigh 2.5 g of sample.
- 6.10.2. Dissolve the sample in 25 mL of freshly boiled and cooled water.
- 6.10.3. Measure the pH of the solution @ 25°C±2°C.
- 6.10.4. Follow the appropriate SOP for pH calibration and measurement.

6.11. RESIDUAL SOLVENTS **Refer to Summary Sheet:**

- 6.11.1. Package and send no less than 30 grams of sample to an approved outside testing facility for Ethanol.

6.12. SINGLE IMPURITIES AND TOTAL IMPURITIES **Refer to Summary Sheet:**

- 6.12.1. Dilute Sulfuric acid: Add 1mL of concentrated sulfuric acid and sufficient water to make 35mL (1N equivalent).
- 6.12.2. Sample solution (10mg/mL Decanoic Acid): Dissolve 113 mg of Sodium Decanoate in 5 mL of water, add 1 mL of dilute sulfuric acid and extract with 10 mL of ethyl acetate. Separate and collect the organic layer, dry organic layer over anhydrous sodium sulfate and then filter before use.
- 6.12.3. Working Standard: Prepare a 0.75mg/mL solution of decanoic acid in ethyl acetate, record precise sample weight and dilution factor. (e.g. 100mL volumetric flask is equal to 100 dilution factor, 25mL volumetric flask is equal to a 25dilution factor).
- 6.12.4. Related compound solution (0.1mg/mL): Prepare a related compound standard solution containing 0.1mg/mL of octanoic acid, 0.1mg/mL nonanoic acid, and 0.2mg/mL of decanoic acid in ethyl acetate.
 - 6.12.4.1. Accurately weigh 100mg of each related substance and 200mg of decanoic acid and transfer each quantitatively to a 25mL volumetric flask. Dissolve and dilute to volume with ethyl acetate.
 - 6.12.4.2. Pipette 0.25mL of the solution into a 10mL volumetric flask, mix and QS to volume.
- 6.12.5. Blank: Pipette 5mL of water, add 1mL of dilute sulfuric acid (1N) and extract with 10mL of ethyl acetate. Separate and collect the organic layer, dry organic layer over anhydrous sodium sulfate and then filter before use.
- 6.12.6. The GC-MS is used with the following parameters:
 - 6.12.6.1. Method File: Sodium Decanoate Single Impurity Quantitative.qgm
 - 6.12.6.2. Column: 0.25-mm x 30-m fused silica; coated with a 0.25-mm layer of phase G25
 - 6.12.6.3. Injection Port: 250°C
 - 6.12.6.4. Detector: 250°C
 - 6.12.6.5. Flow rate: 1.5 mL/min
 - 6.12.6.6. Carrier gas: Helium
 - 6.12.6.7. Injection volume: 1 µL
 - 6.12.6.8. Autosampler: AOC-20i
 - 6.12.6.9. Injection type: Split
 - 6.12.6.10. Split ratio: 5
 - 6.12.6.11. Oven Temperature: See Table 1

TABLE 1: OVEN TEMPERATURE GRADIENT

Initial Temperature (°C)	Temperature Ramp (°C/min)	Final Temperature (°C)	Hold Time at Final Temperature (min)
100	-	100	1
100	5	220	10

- 6.12.7. Determine the peak responses for system suitability.
- 6.12.8. System suitability: The system suitability is determined using *related compound solution (0.1mg/mL)* and the signal-to-noise ratio shall be not less than 14.5 for each related substance.
- 6.12.8.1. Open Post Run Software
- 6.12.8.2. Open associated suitability file.
- 6.12.8.3. Click Qualitative
- 6.12.8.3.1. Select Peak *Integrate for TIC (All Group)*
- 6.12.8.3.2. *Set Min/Area to 100,000 and width to 5 seconds.*
- 6.12.8.3.3. *Click Program* and Set Program to:
- 6.12.8.3.3.1. T=0 Integration off
- 6.12.8.3.3.2. T=10 Integration on
- 6.12.8.3.3.3. T=22 Integration off
- 6.12.8.3.3.4. (These times may change if a new column is used and retention times change, ensure the suitability peak is bracketed by integration)
- 6.12.8.3.4. Click OK
- 6.12.8.3.5. Click OK
- 6.12.8.3.6. Click *OK*
- 6.12.8.4. Click Compound Table
- 6.12.8.4.1. Select Wizard (New...)
- 6.12.8.4.2. Click Next (1)
- 6.12.8.4.3. Click Next (2)
- 6.12.8.4.4. Click Next (3)
- 6.12.8.4.5. Click Next (4)
- 6.12.8.4.6. Click Next (5)
- 6.12.8.4.7. Click Next (6)
- 6.12.8.4.8. Click Finish
- 6.12.8.5. Select “View” from the right hand side of the screen.
- 6.12.8.6. Assign peak ID based upon area count and retention time. The largest peak will be assigned as decanoic acid. Assign nonanoic acid as the first adjacent impurity with RRT <1 to the decanoic acid peak. Assign octanoic acid as the first adjacent impurity peak to nonanoic acid with a RRT of <1.
- 6.12.9. Click Quantitative
- 6.12.9.1. Select “Peak Integrate for all IDs”
- 6.12.10. Click “Report” from the left hand side bar.
- 6.12.11. Click “Open”
- 6.12.12. Open “System Suitability Report” from data files sodium decanoate OI folder.
- 6.12.13. Ensure the S/N ratio is displayed for octanoic acid and nonanoic acid.
- 6.12.13.1. Print report.
- 6.12.13.2. Sample analysis:
- 6.12.13.2.1. Prepare a sequence including:
- 6.12.13.2.2. Blank Sample Single Injection (1)
- 6.12.13.2.3. System Suitability Single Injection (1)
- 6.12.13.2.4. Inject Working Standard Solution A

The information contained herein is the confidential property of BioSpectra. The recipient is responsible for its safe-keeping and the prevention of unauthorized appropriation, use, disclosure and copying.

- 6.12.13.2.5. Inject Working Standard Solution B
- 6.12.13.2.6. Inject Working Standard Solution A
- 6.12.13.2.7. Inject Working Standard Solution B
- 6.12.13.2.8. Inject Sample(s)
- 6.12.13.2.9. Inject Working Standard Solution A
- 6.12.13.2.10. Inject Working Standard Solution B
- 6.12.13.2.11. Blank Sample Single Injection
- 6.12.13.3. % RSD of the working standards should not exceed 12.0%
- 6.12.14. Calculate the percentage of each impurity using the equation below.

$$Result \% = \frac{(A_u)(D_u)}{(RF_{Ave})(W_{Sample})} \times 100\%$$

6.12.14.1. Where

- 6.12.14.1.1. A_u = Peak response of the individual impurity
- 6.12.14.1.2. D_u = Dilution Factor of Sample Solution (10mL)
- 6.12.14.1.3. RF_{Ave} = Average response factor
- 6.12.14.1.4. RF = Response Factor = $((A_s)(D_s))/(W_{std})(P)$
- 6.12.14.1.5. A_s = Peak area of decanoic acid in working standard solution
- 6.12.14.1.6. D_s = Dilution Factor of Working Standard (mL)
- 6.12.14.1.7. W_{std} = Weight of Decanoic Acid Reference Material (mg)
- 6.12.14.1.8. P = Purity of Decanoic Acid Reference Material Expressed as a Percentage
- 6.12.14.1.9. W_{sample} = Weight of Sodium Decanoate Sample (mg)
- 6.12.15. To perform the integration with the GC-MS software follow the steps below:
 - 6.12.15.1. Open Post Run Software
 - 6.12.15.2. Open associated sample file.
 - 6.12.15.3. Click Qualitative
 - 6.12.15.3.1. Select Peak Integrate for TIC (All Group)
 - 6.12.15.3.2. Set Min/Area to: $(A_s * 0.0013)$ and width to 10 seconds.
 - 6.12.15.3.3. 0.01% Reporting Threshold: $0.75\text{mg/mL} * 0.0013 = <0.01\%$ of 10mg/mL
 - 6.12.15.3.4. Calculation for Reporting Threshold: $(0.75 * 0.0013) = 0.000975$
 $(0.000975/10) * 100 = 0.00975$ or 0.01%
 - 6.12.15.3.5. Click Program and set program to:
 - 6.12.15.3.5.1. T=0 Integration off
 - 6.12.15.3.5.2. T=3 Integration on (These times may change if a new column is used and retention times change, ensure the suitability peak is bracketed by integration)
 - 6.12.15.3.6. Click OK
 - 6.12.15.3.7. Click OK
 - 6.12.15.3.8. Click OK
 - 6.12.15.4. Click Compound Table
 - 6.12.15.4.1. Select Wizard (New...)
 - 6.12.15.4.2. Click Next (1)
 - 6.12.15.4.3. Click Next (2)
 - 6.12.15.4.4. Click Next (3)
 - 6.12.15.4.5. Click Next (4)
 - 6.12.15.4.6. Click Next (5)
 - 6.12.15.4.7. Click Next (6)
 - 6.12.15.4.8. Click Finish

The information contained herein is the confidential property of BioSpectra. The recipient is responsible for its safe-keeping and the prevention of unauthorized appropriation, use, disclosure and copying.

- 6.12.16. Click “Report” from the left hand side bar.
- 6.12.17. Click “Open”
- 6.12.18. Open “Sodium Decanoate Report” from data files sodium decanoate OI folder.
- 6.12.19. Ensure all impurities are visible in the report chart; refer to each by retention time.
- 6.12.20. Print report.
- 6.12.21. Calculate each integrated impurity as %w/w using the calculation in 6.11.14 and report total impurity as the sum of all calculated impurity.
- 6.12.22. In the case where an impurity is unknown and requires identification, run the sample solution under “Sodium Decanoate Single Impurity.qgm” and notify a supervisor to identify the impurity.

6.13. SODIUM

Passes Test:

- 6.13.1. Dissolve 0.1 g of sample in 2 mL of water.
- 6.13.2. Add 2 mL of 15% potassium carbonate and heat to boiling. No precipitation is formed.
- 6.13.3. Add 4 mL of potassium pyroantimonate TS and heat to boiling.
- 6.13.4. Allow to cool in ice water and, if necessary, rub the inside of the glass vessel with a glass rod. A passing sample will produce a dense precipitation.

6.14. SOLUBILITY IN WATER

Passes Test:

- 6.14.1. Weigh at least 2g of sample directly in to a clean and dry beaker. Dissolve in 40mL of purified water by using a magnetic stir bar at a low to medium stir speed at room temperature, do not heat.
- 6.14.2. Inspect the sample solution for complete solubility and any foreign matter.
 - 6.14.2.1. Note: the following steps are for Monitoring only:
 - 6.14.2.1.1. The sample is expected to be slightly turbid and not completely colorless.
 - 6.14.2.1.2. Measure NTU result and document. Follow the SOP for Bangor Portable Turbidimeter and Calibration DCN: 16-001317. Monitor specification is <50NTU.
 - 6.14.2.1.3. Filter the solution through a white round filter paper, inspect for insoluble matter.
 - 6.14.2.1.4. After this solution is analyzed for Solubility, analyze for particulates.
 - 6.14.2.1.5. Inspect the magnetic stir bar and verify that no particulates are present on the stir bar. If present notify the appropriate personnel immediately.

6.15. WATER (KARL FISCHER)**Refer to Summary Sheet:**

- 6.15.1. Weigh ~1g of sample into a glass weighing spoon and tare the balance.
- 6.15.2. Transfer the sample to the Karl Fischer vessel by removing the rubber septum and adding the sample into the titration vessel.
 - 6.15.2.1. Do not leave the rubber septum open for long periods of time as this will allow moisture to enter the titration vessel.
- 6.15.3. Return the weighing spoon to the balance, making sure not to lose any sample that was left behind. Once the weight stabilizes, Record the sample weight and transfer to instrument.
- 6.15.4. Check to make sure there is no residual sample stuck to the sides of the titration vessel.
 - 6.15.4.1. If any sample is stuck to the side, stop the stir bead and swirl the Karl Fischer vessel to rinse the sides.
- 6.15.5. Once the method begins, check to ensure the sample is fully dissolved before the titration begins (i.e. before the stir command completes).
- 6.15.6. The moisture content will then be determined by the Metrohm Auto Titrando 907.

$$\% \text{ Moisture} = \frac{(\text{mL of Composite 5})\left(\frac{\text{mg}}{\text{mL}} \text{ of Composite 5}\right)(0.1)}{\text{Sample Weight (g)}}$$

Signature Manifest

Document Number: BSI-ATM-0029

Revision: 6.2

Title: Sodium Decanoate Testing Methods

Effective Date: 19 Sep 2022

All dates and times are in US/Eastern.

Sodium Decanoate Testing Methods

Change Request

Name/Signature	Title	Date	Meaning/Reason
Amy Yencho (AMY.YENCHO)	Vice President, Laboratory Services	08 Sep 2022, 11:09:21 AM	Approved
Virginia Pena (VIRGINIA.PENA)	Document Control Technician I	08 Sep 2022, 11:11:44 AM	Approved

Departmental Approval

Name/Signature	Title	Date	Meaning/Reason
Amy Yencho (AMY.YENCHO)	Vice President, Laboratory Services	16 Sep 2022, 12:06:43 PM	Approved

Author Approval

Name/Signature	Title	Date	Meaning/Reason
Amy Yencho (AMY.YENCHO)	Vice President, Laboratory Services	16 Sep 2022, 12:04:15 PM	Approved

Quality Approval

Name/Signature	Title	Date	Meaning/Reason
Carissa Albert (CARISSA.ALBERT)	Quality Assurance Manager	19 Sep 2022, 08:27:06 AM	Approved

Training Checkpoint

Name/Signature	Title	Date	Meaning/Reason
Kimberly Hartzell (KIMBERLY.HARTZELL)	Training Manager	19 Sep 2022, 09:09:00 AM	Approved

Set Date

Name/Signature	Title	Date	Meaning/Reason
Virginia Pena (VIRGINIA.PENA)	Document Control Technician II	19 Sep 2022, 09:24:50 AM	Approved