




TROMETHAMINE



BIO ACTIVE GRADE REGULATORY PACKET

Signature/Date:

The information contained herein is the confidential property of BioSpectra. The recipient is responsible for its safe-keeping and the prevention of unauthorized appropriation, use, disclosure and copying.



TABLE OF CONTENTS

1. TROMETHAMINE (TRIS) BIO ACTIVE GRADE: 3

1.1. GENERAL PRODUCT INFORMATION: 3

1.2. MANUFACTURING, PACKAGING RELEASE SITE AND SUPPLIER INFORMATION: 3

1.3. PHYSICO-CHEMICAL INFORMATION: 3

1.4. REGULATORY INFORMATION: 6

1.5. MISCELLANEOUS PRODUCT INFORMATION: 8

1.6. CONTACT INFORMATION: 9

The information contained herein is the confidential property of BioSpectra. The recipient is responsible for its safe-keeping and the prevention of unauthorized appropriation, use, disclosure and copying.



1. TROMETHAMINE (TRIS) BIO ACTIVE GRADE:

1.1. General Product Information:

- 1.1.1. Product Name:
 - 1.1.1.1. Tromethamine
- 1.1.2. Product Code:
 - 1.1.2.1. Current Code: TRIS-2220, TRIS-2250, TRIS-2255, TRIS-2256
- 1.1.3. Scope:
 - 1.1.3.1. This regulatory packet will provide the quality and regulatory information regarding the manufacturing, testing, packaging, storage, release, shipping, and handling of Bio Active Grade Tromethamine manufactured by and at the BioSpectra, Bangor, PA facility.
- 1.1.4. Molecular Formula:
 - 1.1.4.1. $\text{NH}_2\text{C}(\text{CH}_2\text{OH})_3$
- 1.1.5. Molecular Weight:
 - 1.1.5.1. 121.14 g/mol

1.2. Manufacturing, Packaging Release Site and Supplier Information:

- 1.2.1. General Information:
 - 1.2.1.1. BioSpectra manufactures Tromethamine in its Bangor, PA facility. Tromethamine is manufactured, packaged, stored, tested and released at BioSpectra's Bangor, PA facility.
- 1.2.2. Manufacturing:
 - 1.2.2.1. The manufacturing of Tromethamine is performed at BioSpectra's Bangor, PA facility utilizing multiuse equipment. Equipment used in the manufacturing of Tromethamine is cleaned in accordance with BioSpectra's Cleaning Worksheet Procedure.
- 1.2.3. Packaging:
 - 1.2.3.1. The packaging of Tromethamine occurs in the following BioSpectra site:
BioSpectra Bangor, PA Facility: 100 Majestic Way, Bangor, PA 18013.
- 1.2.4. Testing for Release:
 - 1.2.4.1. Testing and release of Tromethamine is performed at the BioSpectra Bangor, PA Facility: 100 Majestic Way, Bangor, PA 18013.
- 1.2.5. GMP Compliance Statement:
 - 1.2.5.1. Bio Active Grade Tromethamine is suitable for use as a non-Sterile Active Pharmaceutical Ingredient manufactured in accordance with the ICH Q7 Good Manufacturing Practice Guide. This grade of Tromethamine is not suitable to be used as a Sterile or Injectable Active Pharmaceutical Ingredient, Drug Product or Household Item.


1.3. Physico-Chemical Information:

- 1.3.1. CAS Number:
 - 1.3.1.1. CAS # 77-86-1



- 1.3.2. Origin:
 - 1.3.2.1. The origin of Tromethamine is through synthetic chemical manufacturing using approved raw materials, which are further purified in accordance with ICH Q7 guidelines. Only raw materials of synthetic origin are used in the synthesis and purification of Tromethamine.
- 1.3.3. Synonyms:
 - 1.3.3.1. Trometamol
 - 1.3.3.2. THAM
 - 1.3.3.3. Tris (Hydroxymethyl) Aminomethane
 - 1.3.3.4. Trisaminol
 - 1.3.3.5. Trismethylolaminomethane
 - 1.3.3.6. 2-Amino-2-(Hydroxymethyl) propane-1,3-diol
 - 1.3.3.7. Tromethamine
- 1.3.4. Morphological Form:
 - 1.3.4.1. White / Crystals
- 1.3.5. Manufacturing Process:
 - 1.3.5.1. The Tromethamine, Bio Active Grade manufacturing process is performed by the following:


Approved Supplier Raw Material Synthesis



The information contained herein is the confidential property of BioSpectra. The recipient is responsible for its safe-keeping and the prevention of unauthorized appropriation, use, disclosure and copying.



BioSpectra Tris Bio Active Manufacturing Process



1.3.6. Specifications:

1.3.6.1. Available upon request.



1.4. Regulatory Information:

- 1.4.1. Compendial Compliance:
 - 1.4.1.1. USP Monograph
- 1.4.2. Master File:
 - 1.4.2.1. Drug Master File (DMF) is available for this product.
 - 1.4.2.2. EDQM Certificate of Suitability is currently not available for this product.
- 1.4.3. REACH:
 - 1.4.3.1. Refer to the Tromethamine Safety Data Sheet for the REACH Number or contact your Commercial Team Representative for further information.
- 1.4.4. BSE/TSE Statement:
 - 1.4.4.1. Tromethamine, Bio Active Grade is a synthetic chemical and has been evaluated for the source of the raw materials used in its production through the Supplier Qualification Program. BioSpectra can state that BSE/TSE is not a concern based on this evaluation. Tromethamine, Bio Active Grade and its raw materials are not of animal origin. Tromethamine, Bio Active Grade complies with the note for guidance in minimizing the risk of transmitting animal spongiform encephalopathy agents via human and veterinary medicinal products (EMA/410/01).
- 1.4.5. Allergens Statement:
 - 1.4.5.1. Tromethamine, Bio Active Grade manufactured by BioSpectra and its raw materials are not manufactured using with or using any of the following substances: Cereals containing gluten (i.e., wheat, rye, barley, oats, spelt, kamut or their hybridized strains) and products thereof, Crustaceans and products thereof, Eggs and products thereof, Fish and products thereof, Peanuts and products thereof, Soybeans and products thereof, Milk and products thereof (including lactose), Celery and products thereof, Mustard and products thereof, Sesame seeds and products thereof, Lupin and products thereof, Molluscs and products thereof, Sulphur dioxide and sulphites at concentrations of more than 10 mg/kg or 10 mg/litre expressed as SO₂, Nuts i.e., Almond (*Amygdalus communis* L.), Hazelnut (*Corylus avellana*), Walnut (*Juglans regia*), Cashew (*Anacardium occidentale*), Pecan nut (*Carya illinoensis* (Wangenh.) K. Koch), Brazil nut (*Bertholletia excelsa*), Pistachio nut (*Pistacia vera*), Macadamia nut and Queensland nut (*Macadamia ternifolia*) and products thereof, Beef, Chicken, Pork, Azo Dyes, Benzoic Acid, Tartrazine, Vanillin, Cocoa, Cinnamon, Coriander, Yeast, Glutamate, Legumes, and Corn. BioSpectra has evaluated the Raw Material Supply through the Supplier Qualification Program. Tromethamine, Bio Active Grade does not contain any of the substances causing allergies or intolerances as listed under Annex II of Regulation (EU) No. 1169/2011.
- 1.4.6. Genetically Modified Organisms (GMO) Statement:
 - 1.4.6.1. Tromethamine, Bio Active Grade is a synthetic chemical and has been evaluated for the source of the raw materials used in its production through the Supplier Qualification Program. BioSpectra can state that genetic modification is not a concern based on this evaluation.

The information contained herein is the confidential property of BioSpectra. The recipient is responsible for its safe-keeping and the prevention of unauthorized appropriation, use, disclosure and copying.



1.4.7. Residual Solvents Statement:

1.4.7.1. BioSpectra can state based on the manufacturing process and the controlled handling, storage, and analysis of this product, that the Tromethamine, Bio Active Grade manufactured by BioSpectra complies with the requirements and specifications listed in the current USP method <467> Tables 1, 2, 3, or 4. BioSpectra has evaluated the raw material supply through the Supplier Qualification Program and can state that the Tris raw material supplied to BioSpectra is manufactured using the Class 2 residual solvents Methanol and Nitromethane, with expected concentrations meeting the requirements of ICH Q3C for residual solvents. BioSpectra analyzes Tromethamine, Bio Active Grade and its raw materials for Methanol and Nitromethane during Degradation and Impurity Profiling as part of process validation, with results meeting the requirements. Tromethamine is additionally analyzed annually for Methanol and Nitromethane, with specifications NMT 3000ppm and NMT 50ppm respectively.

1.4.8. Metal Catalyst and Metal Reagent Residues Statement:

1.4.8.1. BioSpectra does not intentionally add or use any metal catalysts in the manufacturing process of Tromethamine, Bio Active Grade. BioSpectra has evaluated the raw material supply through the Supplier Qualification Program and can state that the Tris raw material is made from an ingredient which uses a very small amount of a nickel compound in the manufacture of that ingredient. BioSpectra additionally analyzes Tromethamine for elemental impurities during Degradation and Impurity Profiling as part of process validation.

1.4.9. Pallet Statement:

1.4.9.1. BioSpectra can state that the wooden pallets used in the packaging and shipping of Tromethamine, Bio Active Grade manufactured at BioSpectra are ISPM 15 compliant.

1.4.10. Elemental Impurities Statement:

1.4.10.1. BioSpectra's Tromethamine, Bio Active Grade material has been profiled for elemental impurities via ICP utilizing USP <232> and <233> in accordance with ICH Q3D. The results are reported in the associated Elemental Impurity Profile and are available upon request.

1.4.11. Melamine Statement:

1.4.11.1. Tromethamine, Bio Active Grade has been evaluated for the source of the raw materials used in its production through the Supplier Qualification Program. BioSpectra can state that melamine is not a risk based on this evaluation. BioSpectra additionally analyzes Tromethamine annually for melamine, with a specification of 2.5 mg/kg max. Tromethamine, Bio Active Grade is not considered at-risk per the US FDA Guidance for Industry, Pharmaceutical Components at Risk for Melamine Contamination.

The information contained herein is the confidential property of BioSpectra. The recipient is responsible for its safe-keeping and the prevention of unauthorized appropriation, use, disclosure and copying.



1.5. Miscellaneous Product Information:

1.5.1. Description of Batch:

1.5.1.1. The Tromethamine manufacturing process is a batch or continuous process where expected batch yields are established during validation in accordance with the Manufacturing Process Validation Master Plan. Individual batch yield is additionally determined for each manufactured batch and documented in the respective batch record.

1.5.2. Lot/Batch Numbering System:

1.5.2.1. The lot numbering system at BioSpectra employs the following format per BSI-DGM-0009 BioSpectra Lot Number Identification:

1.5.2.2. A sample lot number would appear as:

1.5.2.2.1. QS6: TRIS-0124-00001

1.5.2.2.1.1. The first four digits are alpha digits which indicate the material manufactured, where TRIS represents Tromethamine. The fifth and sixth digits are numeric digits which indicate the site of final packaging, where 01 represents the Bangor, PA facility. The seventh and eighth digits are numeric digits which indicate the year the batch record was issued, where 24 represents 2024. The final five digits are numeric digits which indicate the sequential batch number, where 00001 represents the first Tromethamine batch of 2024 that is automatically generated by the ERP system. The sequential batch number automatically resets on the first of the new calendar year.

1.5.2.2.2. QS7: TRIS-E04-0124-0001

1.5.2.2.2.1. The first four digits are alpha digits which indicate the material, where TRIS represents Tromethamine. The fifth, sixth, and seventh digits are alphanumeric digits which indicate the location of manufacturing. The eighth and ninth digits are numerical digits which indicate the month of work order issuance, where 01 represents January. The tenth and eleventh digits are numerical digits which indicate the year of work order issuance, where 24 represents 2024. The final four digits are numerical digits which indicate the sequential batch number, where 0001 represents the first Tromethamine batch of 2024 that is automatically generated by the ERP system. The sequential batch number automatically resets on the first day of each calendar year.



1.5.3. Expiration Date and/or Recommended Re-evaluation Interval:

1.5.3.1. The current recommended retest or expiration date for Tromethamine, Bio Active Grade is available in the BioSpectra Product Retest and Expiration Date List, DCN: BSI-LST-0239, and is based on current available stability data in accordance with BioSpectra's Stability Testing Program. Additionally, the recommended Retest or Expiration Date will be available on the Product Specific Certificate of Analysis, as applicable.

1.5.4. Storage and Shipping Conditions:

1.5.4.1. Store in a cool, dry place. Store in a cool, dry, well-ventilated area away from incompatible substances. Keep away from acids. Keep containers tightly closed.

1.5.5. Packaging:

1.5.5.1. Packaging information is available through the following:
<https://biospectra.us/packaging>

1.6. Contact Information:

1.6.1. <https://www.biospectra.us/about-us/commercial-marketing-team>