

May 7th, 2021 Revision 2

ELEMENTAL IMPURITY STATEMENT

Sodium Decanoate GMP

BioSpectra's Sodium Decanoate Bio Excipient Grade material was profiled for elemental impurity via ICP-MS utilizing USP <232> and <233> in accordance with ICHQ3D with reference to parenteral application and associated limits. Osmium is not intentionally added and therefore is not required for parenteral use and was excluded from these analyses. In addition to Osmium being excluded from parenteral limits as per USP <232>, USP <233> and ICH Q3D, the standard preparation of Osmium stock solution could potentially produce Osmium Tetroxide and for safety considerations, should not be attempted to be analyzed utilizing a nitric acid based ICP-MS system. The results comply with the limits as per USP <232>, USP <233> and ICH Q3D.

Current Product Number	Historic Product Number	
NDEC-3201	ND3201	
NDEC-3220	ND3220	
NDEC-3320	ND3320	



Table 1. Parenteral limits as per USP <232> and <233>

Element	Class	Impurity Limit (μg/g)	Result:
Cadmium (Cd)	1	0.2	<0.2 μg/g
Lead (Pb)	1	0.5	<0.5 μg/g
Arsenic (As)	1	1.5	<1.5 μg/g
Mercury (Hg)	1	0.3	<0.3 μg/g
Cobalt (Co)	2A	0.5	<0.5 μg/g
Vanadium (V)	2A	1	<1 μg/g
Nickel (Ni)	2A	2	<2 μg/g
Thallium (TI)	2B	0.8	<0.8 μg/g
Gold (Au)	2B	10	<10 μg/g
Palladium (Pd)	2B	1	<1 μg/g
Iridium (Ir)	2B	1	<1 μg/g
Rhodium (Rh)	2B	1	<1 μg/g
Ruthenium (Ru)	2B	1	<1 μg/g
Selenium (Se)	2B	8	<8 μg/g
Silver (Ag)	2B	1	<1 μg/g
Platinum (Pt)	2B	1	<1 μg/g
Lithium (Li)	3	25	<25 μg/g
Antimony (Sb)	3	9	<9 μg/g
Barium (Ba)	3	70	<70 μg/g
Molybdenum (Mo)	3	150	<150 μg/g
Copper (Cu)	3	30	<30 μg/g
Tin (Sn)	3	60	<60 μg/g
Chromium (Cr)	3	110	<110 μg/g

For further information, please contact info@biospectra.us

Cassie Baun

Compliance Specialist