


HEPES, Sodium Salt

INTENDED FOR USE IN COMMERCIAL MANUFACTURING, DIAGNOSTIC AND LABORATORY APPLICATIONS

HEPES Sodium Salt is a zwitterionic salt which is used in the preparation of a HEPES buffer solution. Its pH range makes an ideal buffer for pH maintenance. HEPES Sodium Salt interferes with the Folin protein assay. The tendency of this buffer to form radicals makes it unsuitable for redox studies.



CAS #: 75277-39-3

Molecular Formula:

$C_8H_{17}N_2NaO_4S$

Solubility in Water (g/L): 990

F.W.: 260.29 g/mol

pH @ 20°C: 4.5 - 6.0

Useful pH: 7.8 - 8.2

pK_a @ 20°C: 7.5

BIO ULTRA GRADE | Product Code: HEPN-5251 | Previously: HN5251

$C_8H_{17}N_2NaO_4S$ · F.W. 260.29 g/mol. · CAS# 75277-39-3

These are general specifications. BioSpectra will customize our products to meet your quality based requirements.

ANALYSIS	SPECIFICATIONS
Appearance (Color)	White
Appearance (Form)	Powder
Infrared Spectrum	Conforms to Reference
Water (by Karl Fischer)	≤ 3.0%
Titration with H ₂ SO ₄	≥ 96.0 %
RNase Detection	Pass
pH 1.0% aqueous	9.5 – 10.5
Endotoxin	≤ 50 EU/g

Key Compliance Attributes of Bio Ultra Grade

Suitable for Research and Diagnostic Application

C of A Available for Each Lot

BioSpectra Supply Chain Audit Trail

Product Origin Statement

Customer Quality Audits

Customized Additional Specifications

Enzyme Tested

General Product Description:

- HEPES Sodium Salt is a White Crystalline product.
- Molecular Formula: $C_8H_{17}N_2NaO_4S$
- Molecular Weight: 260.29 g/mol.
- CAS Number: 75277-39-3.
- Synonyms: 4-(2-Hydroxyethyl)piperazine-1-ethanesulfonic acid sodium salt.

Compliance:

Material represented by this grade is suitable to be used only as the following: Quality System Manufactured Process Chemical for use in further manufacturing or as a reagent chemical for commercial manufacturing, diagnostic and laboratory research. The material represented by this grade is not suitable to be used as an Excipient, Active Pharmaceutical Ingredient, Drug, Drug Product or household item.

Retest Date:

Two years from date of manufacture, testing and packaging

Storage and Shipping Conditions:

Ship and Store in ambient temperature.

Package Sizes:

10kg, 25kg and 50kg pails.