

RECEIVED  
JUN 10 2020

## Product Data


Life Latch® Drum  
12.0 Gallon

UN Certification: 1H2/Y53/S & 1H2/Y40/S  
No Gasket      Gasket



### Cover Specification | Part No. 600104


Origin of Materials	100% US Manufactured
Material	High Density Polyethylene (HDPE)
Mfg. Method	Injection Molding
Weight	1.85 lbs / 840g
Height	2.51"
Min Wall Thickness	.100"
Bottom OD	16.31"
Gasket	Neoprene / Available w/o Gasket
Standard Color	White <i>Please inquire about custom colors</i>
Case Quantity	180 per pallet



Available with  
Buna Vent

### Pail Specification | Part No. 500034

Origin of Materials	100% US Manufactured
Material	High Density Polyethylene (HDPE)
Mfg. Method	Injection Molding
Weight	5.83 lbs / 2646 g
Min. Wall Thickness	.140"
Exterior Height	21.38" With Cover 21.5"
Top Inside Diameter	14.04"
Bottom Outside Diameter	13"
Shipping Quantity	xx Pails and Covers Per Case 54 Units per Pallet 2,160 Units per Truckload
Standard Color	White <i>Please inquire about custom colors</i>
Top Load Rating	450 lbs



Bottom OD 13"

