



Attn: Sarah Demaio
Biospectra
1474 Rockdale Lane
Stroudsburg, PA 18360

RECEIVED
APR 28 2020

SPECIFICATION SHEET

ITEM DESCRIPTION: Sierra 3100 Linear Low Density Hexene Blend Non-Slip

DIMENSION: 34 x 60 (34" wide x 60" Length)

TYPE: Clear Bags on Rolls Non-Slip

GAUGE: .003 ± 10%

COLOR: Clear

PACKAGING: 200 Bags/Roll

FDA STATUS: .21 CFR 177, 1520 Olefin Polymers Par © 3.2