



July 11th, 2023
Revision 1

LLDPE NON-SLIP LINER EU COMPLIANCE

Dextran Sulfate Sodium Salt, 8000MW GMP

BioSpectra has evaluated the packaging component supply through the Supplier Qualification Program, and can state that the current packaging utilized for Dextran Sulfate Sodium Salt, 8000MW GMP, Bio Pharma Grade manufactured by BioSpectra does not contain slip agents based on this evaluation. The Linear Low-Density Polyethylene (LLDPE) liners utilized by BioSpectra for material packaging are manufactured from virgin polyethylene FDA approved material, conforming to 21 CFR 1520. The polyethylene resins adhere to Regulation EC 1935/2004 and Regulation EC 2023/2006, as listed under EU No. 10/2011.

Product Number	Liner Size
COMP-0014	24"x36"
COMP-0012	30"x48"

For further information, please contact info@biospectra.us

Cassie Baun
Senior Compliance Specialist