

STABILITY STATEMENT

HEPES GMP

BioSpectra can state that HEPES, Bio Excipient Grade, and Bio Pharma Grade manufactured by BioSpectra undergoes a 36-month stability testing program. A 24-month retest date is currently assigned based on the Long-Term Stability Data with initial regulated conditions for temperature at 10-40°C and humidity (monitor), and temperature at 23-27°C and humidity at 55-65% for the duration of the study. The packaging configurations included in this study are Poly/Poly, 2 Poly/Fiber, and Labline.

Current Product Number	Historic Product Number
HEPE-3220	HE3220
HEPE-3221	HE3221
HEPE-3250	HE3250
HEPE-3251	HE3251
HEPE-4220	HE4220

For further information, please contact info@biospectra.us

Cassie Baun

Cassie Baun
Compliance Specialist