

UOM	Product Code Suffix	Package Size	Description
kg	01	1 kg	1500ml White HDPE bottle, blue lid
kg	05	5kg	3.5-gallon White HDPE round pail, blue lid, LDPE liner
kg	10	10 kg	5-gallon White HDPE round pail, blue lid, LDPE liner
kg	12	12 kg	5-gallon White HDPE round pail, blue lid, LDPE liner
kg	25	25 kg	12-gallon White HDPE round pail, blue lid, LDPE liner
kg	50	50 kg	30-gallon Dark blue HDPE drum, LDPE liner
L	69	200 L	200 Liter Sterile Single Use Bio-Processing Container (BPC Shell/Bag)
L	70	1L	1L Non-fluorinated bottle, black lid, unfluorinated, single shipper
L	71	6 x 1L	Case containing 6 x 1L Nonfluorinated bottle, black lid
L	72	200 L	55-Gallon Black HDPE drum
L	73	4 x 3.78L	Case containing 4 x 1-gallon Non-fluorinated bottle, black lid
L	75	Up to 1005 L	275-gallon Tote
L	76	Up to 1025 L	1000 Liter Sterile Single Use Bio-Processing Container (BPC Shell/Bag)
L	78	10 L	2.5-gallon carboy
L	79	Up to 1135 L	330-gallon tote
L	82	Up to 100L	100L Sterile Single Use Bio- Processing Container (BPC Shell/Bag)
g	89	200 gr	1000 ml White HDPE bottle, blue lid- samples
g	91	50 g	250 ml White HDPE bottle, blue lid- samples
g	92	100 g	250 ml White HDPE bottle, blue lid- samples
g	93	500 g	1000 ml White HDPE bottle, blue lid- samples
g	94	250gr	1000 ml White HDPE bottle, blue lid- samples
g	95	10gr	250 ml White HDPE bottle, blue lid- samples
L	96	Up to 14,000L (TBD)	Bulk Tank Wagon
kg	97	180 kg	55 gallons Black HDPE drum (sold in kg, not liters)

Special Consideration- product (TRED, NDEC) needs larger pack size due to morphology- may want to double check

UOM	Product Code Suffix	Package Size	Description
kg	03	3 kg	3.5-gallon White HDPE round pail, blue lid, LDPE liner
kg	05	5kg	5-gallon White HDPE round pail, blue lid, LDPE liner
kg	10	10 kg	12-gallon White HDPE round pail, blue lid, LDPE liner