

### Packaging Dimensions, Stack and FCL Quantities

Package Type	Dimensions of Package (inches)	Number of Package Types per Pallet	Dimensions of a Full Pallet
2.5 Gallon/ 10 Liter Carboy	7.75" W x 9.13" L x 12.25" H	3 Stack HT US Pallet: 90 3 Stack Euro Pallet: 80	HT US Pallet: 48" L x 42" W x 42.75" H (+/- 0.5")
1 Liter HDPE Bottles, Non-Fluorinated	1.48" – 3.25" W x 10.455" – 10.615" H	3 Stack HT US Pallet: 72 3 Stack Euro Pallet: 64	HT US Pallet: 48" L x 42" W x 40.5" H (+/- 0.5")
4 X 4 Liter Bottles	13.75" ± 0.5" H x 5.947" – 6.09" W	3 Stack HT US Pallet: 27 3 Stack Euro Pallet: 18	HT US Pallet: 48" L x 42" W x 51.75" H (+/- 0.5")
5 Gallon/10kg White Round Poly Pail	10.25" W x 16" H	2 Stack HT US Pallet: 24 2 Stack Euro Pallet: 24	2 Stack HT US Pallet: 48" W x 42" L x 36.75" H 2 Stack Euro Pallet: 47" W x 31.5" L x 37.75" H
		3 Stack HT US Pallet: 36 3 Stack Euro Pallet: 36	3 Stack HT US Pallet: 48" W x 42" L x 52.75" H 3 Stack Euro Pallet: 47" W x 31.5" L x 53.75" H
12 Gallon/25kg White Round Poly Pail	13" W x 21.5" H	2 Stack HT US Pallet: 16 2 Stack Euro Pallet: 12	2 Stack HT US Pallet: 48" W x 42" L x 47.5" H 2 Stack Euro Pallet: 47" W x 31.5" L x 48.75" H
30 Gallon/50kg Light and Dark Blue Poly Drum	20.75"W x 30.1"H	HT US Pallet: 4 Euro Pallet: 2	HT US Pallet: 48" W x 42" L x 47.5" H Euro Pallet: 47" W x 31.5" L x 48.75" H
55 Gallon/200 Liters Black Poly Drum	23" W x 36" H	HT US Pallet: 4 Euro Pallet 2	HT US Pallet: 48" W x 42" L x 40.75" H Euro Pallet: 47" W x 31.5" L x 41.75" H
55 Gallon/200 Liters Natural Poly Drum	23.3" W x 34.8" H	HT US Pallet: 4 Euro Pallet: 2	HT US Pallet: 48" W x 42" L x 39.55" H Euro Pallet: 47" W x 31.5" L x 40.55" H
275 Gallon/ 940 Liter Tote	39"W x 47"L x 46"H	<p style="text-align: center;">Index 1 INCH = 2.54 CM HT = HEAT TREATED FCL = FULL CONTAINER LOAD</p>	
330 Gallon/ 1135 Liter Tote	40"W x 48"L x 54"H		
HT US Pallet	48" W x 42" L x 4.75" H		
Euro Pallet	47" W x 31.5" L x 5.75" H		
Wooden CRATE(Height Includes HT US Pallet)	48" L x 48" W x 48" H		
Wooden Crate for Euro Pallet	48" L X 30" W X 30" H		

Package Type	Number of Package Types Per FCL- 20 Foot Container	Number of Package Types Per FCL- 40 Foot Container
HT US Pallet	8 Pallets	18 Pallets
Euro Pallet	8 Pallets	18 Pallets
275 Gallon/940 Liter Tote	8 Totes	20 Totes
330 Gallon/ 1135 Liter Tote	8 Totes	18 Totes
5 Gallon/10kg White Round Poly Pail	2 Stack HT US Pallet: 192 2 Stack Euro Pallet: 192	2 Stack HT US Pallet: 432 2 Stack Euro Pallet:432
	3 Stack HT US Pallet: 288 3 Stack Euro Pallet: 288	3 Stack HT US Pallet: 648 3 Stack Euro Pallet: 648
12 Gallon/25kg White Round Poly Pail	2 Stack HT US Pallet: 128 2 Stack Euro Pallet: 96	2 Stack HT US Pallet: 288 2 Stack Euro Pallet: 216
30 Gallon/50kg Light and Dark Blue Poly Drum	HT US Pallet: 32 Euro Pallet: 16	HT US Pallet: 72 Euro Pallet: 36
55 Gallon/200 Liter Black Poly Drum	HT US Pallet: 32 Euro Pallet: 16	HT US Pallet: 72 Euro Pallet: 36
55 Gallon/200 Liter Natural Poly Drum	HT US Pallet: 32 Euro Pallet: 16	HT US Pallet: 72 Euro Pallet: 36

Package Type	Case Dimensions	Number of Packages Per CRATE
2.5 Gallon/ 10 Liter Carboy	N/A	3 Stack 48x48 Crate: 90 2 Stack EURO Crate: 40
1 Liter HDPE Bottles, Non-Fluorinated	11.5" H x 7" W x 10.25" L	3 Stack 48x48 Crate: 72 4 Stack EURO Crate: 32
4 X 4 Liter Bottles	15.25" H x 12.75" W x 12.75" L	2 Stack 48x48 Crate: 18 3 Stack EURO Crate: 12
3.5 Gallon / 5kg White Round Poly Pail	N/A	3 Stack 48x48 Crate: 27 2 Stack EURO Crate: 12
5 Gallon/10kg White Round Poly Pail	N/A	3 Stack 48x48 Crate: 27 1 Stack EURO Crate: 7
12 Gallon/25kg White Round Poly Pail	N/A	2 Stack 48x48 Crate: 16 1 Stack EURO Crate: 4
30 Gallon/50kg Light and Dark Blue Poly Drum	N/A	QTY on 48x48 Crate: 4 Euro Crate: N/A