

June 21<sup>st</sup>, 2022 Revision 1

## LABLINE BOTTLE EU COMPLIANCE

BioSpectra has evaluated the packaging component supply through the Supplier Approval Program. The Labline bottles utilized by BioSpectra for material packaging are manufactured from polyethylene (HDPE) compliant with US FDA 21 CFR 177.1520(1)(3)(i)(c1) concerning specification 2.1 and 2.2. These bottles are intended to come into contact with all types of food groups in terms of the Regulation (EU) No. 10/2011.

These bottles additionally comply with the following Regulations, and are UN rated as indicated in the table below.

- 94/62/EC
- RoHS 2011/65/EU

Product Number	Bottle Size	UN Rating
COMP-0025	100g / 250mL / 8.5oz.	3H1/X0.8/S
COMP-0033	500g / 1000mL / 34oz.	3H1/X2.5/S
COMP-0028	1kg / 1.5L / 50oz.	3H1/X4/S

For further information, please contact info@biospectra.us

**Cassie Baun** 

Senior Compliance Specialist