

June 21st, 2022 Revision 2

2.5 GALLON CARBOY EU COMPLIANCE

BioSpectra has evaluated the packaging component supply through the Supplier Approval Program. The 2.5 Gallon (10L) HDPE Jerry Can utilized by BioSpectra for material packaging is made from FDA approved high density polyethylene plastic. These carboys are additionally UN certified for the shipping and handling of hazardous materials as indicated in the table below.

Product Number	Size	UN Rating
COMP-0040	2.5 Gallon	3H1/Y1.2/100
COMP-0041	Cap for 2.5 Gallon	

For further information, please contact info@biospectra.us

Cassie Baun

Senior Compliance Specialist