

June 21st, 2022 Revision 1

HDPE PAIL EU COMPLIANCE

BioSpectra has evaluated the packaging component supply through the Supplier Approval Program. The HDPE pails utilized by BioSpectra for material packaging are manufactured from FDA compliant polyethylene resins in full compliance with the relevant sections of 21 CFR for direct food contact.

Part Number	Packaging Type	UN Rating
COMP-0114	12 Gallon White Pail	1H2/Y53/S (no gasket)
COMP-0103	12 Gallon Blue Lid	1H2/Y40/S (with gasket)
COMP-0113	5 Gallon White Pail	
COMP-0112	5 Gallon White Pall	1H2/Y30/S
COMP-0107	5 Gallon Blue Lid	1H2/X11.5/S
COMP-0179	5 Gallon White Lid	
COMP-0115	3.5 Gallon White Pail	1112/1/10/5
COMP-0105	3.5 Gallon Blue Lid	1H2/Y19/S
COMP-0116	2.5 Gallon White Pail	1112/012/5
COMP-0117	2.5 Gallon Blue Lid	1H2/Y13/S
COMP-0206	1.25 Gallon Pail and Blue Lid	1H2/Y6/S
COMP-0111	0.6 Gallon White Pail	1H2/Y4/S
COMP-0102	0.6 Gallon Blue Lid	

These pails are additionally UN rated as detailed in the table below.

For further information, please contact info@biospectra.us

Cassie Baun Senior Compliance Specialist

The information contained herein is the property of BioSpectra. The recipient is responsible for its safe-keeping and the prevention of unauthorized appropriation, use, disclosure, copying or alteration. Decision based on the use of this information is the sole responsibility and liability of the recipient. If you would like a controlled version of this document, please contact info@biospectra.us