

BIOSPECTRA

100 Majestic Way, Bangor, PA 18013 / www.biospectra.us

MES MONOHYDRATE SPECIFICATION SHEET

BIO EXCIPIENT GRADE / ME3220

$C_6H_{13}NO_4S \cdot H_2O$ * F.W. 213.25 g/mol. * CAS# 145224-94-8

ANALYSIS		SPECIFICATIONS
Absorbance	260 nm	0.1000 a.u. max.
	280 nm	0.1000 a.u. max.
Appearance and Color		White / Crystals
Assay		99.0% min.
Chloride		0.005% max.
Color (1M, Alkaline)		Colorless
Enzymes	DNase	None Detected
	RNase	None Detected
	Protease	None Detected
Heavy Metals		2 ppm max.
Identification (IR)		Passes Test
pH (5% solution)		3.1 – 3.5
pH (0.5M)		2.5 – 4.5
pK _a		5.9 – 6.3
Residue on Ignition		0.05% max.
Solubility		Passes Test
Sulfate		0.005% max.
Trace Elements	Arsenic (As)	5 ppm max.
	Copper (Cu)	5 ppm max.
	Iron (Fe)	5 ppm max.
	Lead (Pb)	5 ppm max.