

STABILITY STATEMENT

Guanidine Hydrochloride & 6M Solution Bio Excipient Bio Pharma

GH3101 GH3201 GH3202 GH3220 GH3250 GH4101 GH4120 GH4121 GH4201 GH4220 GH4221 GH4222 GH4250

BioSpectra can state that Guanidine Hydrochloride and Guanidine Hydrochloride 6M Solution, Bio Excipient Grade, and Bio Pharma Grade product codes GH3101, GH3201, GH3202, GH3220, GH3250, GH4101, GH4120, GH4121, GH4201, GH4220, GH4221, GH4222, GH4250 manufactured at BioSpectra undergoes a 36-month stability testing program. A 24-month retest date is currently assigned based on the Long-Term Stability Data with regulated conditions for temperature at 15-30°C and humidity (monitor). The packaging configurations included in this study are Poly, Labline, and Poly/Fiber.

For further information, please contact info@biospectra.us

Cassie Baun

Cassie Baun
Compliance Specialist