

## QUALITY CONTROL TESTING AND VALIDATION STATEMENT

### L-Arginine Hydrochloride

L-Arginine Hydrochloride, Bio Pharma Grade for BioBuffer Solutions Product line is analyzed by BioSpectra's Quality Control laboratory. Analytical methods utilized for the analysis of L-Arginine Hydrochloride, Bio Pharma Grade for BioBuffer Solutions Product line are validated as applicable in accordance with BioSpectra's Analytical Methods Validation Master Plan. The L-Arginine Hydrochloride testing methods listed in the table below are validated in-house in accordance with Compendial Monograph or General Chapters.

Current Product Number
LARH-4220

Analysis Name
Ammonium and Ninhydrin-Positive Substances (USP/EP)
Arsenic, Heavy Metals, Iron
Assay (dried basis/dried substance)

For further information, please contact [info@biospectra.us](mailto:info@biospectra.us)

*Cassie Baum*

**Cassie Baum**

Senior Compliance Specialist