

STABILITY STATEMENT

Trehalose, Dihydrate Bio Excipient

TE3250 TE3251 TE3252 TE3253

BioSpectra certifies that Trehalose, Dihydrate, Bio Excipient Grade product codes TE3250, TE3251, TE3252, TE3253 manufactured at BioSpectra undergoes a 36-month stability testing program. A 24-month retest date is currently assigned based on the Long-Term Stability Data with initial regulated conditions for temperature at 10-40°C and humidity (monitor), and temperature at 23-27°C and humidity at 55-65% for the duration of the study. The packaging configurations included in this study are Poly/Poly and Poly Labline.

For further information, please contact info@biospectra.us

Cassie Baun

Cassie Baun
Compliance Specialist