

## STABILITY STATEMENT

Tris Bio Active

TR2220 TR2221 TR2250

BioSpectra can state that Tris, Bio Active Grade, product codes TR2220, TR2221, TR2250 manufactured at BioSpectra undergoes a 36-month stability testing program. A 24-month retest date is currently assigned based on the Long-Term Stability Data with regulated conditions for temperature at 15-30°C and humidity (monitor). The packaging configurations included in this study are Poly/Poly and Poly Labline.

For further information, please contact [info@biospectra.us](mailto:info@biospectra.us)



**Cassie Baun**  
Compliance Specialist