

August 2nd, 2021 Revision 3

SUPPLY CHAIN STATEMENT

Tris API

The cGMP manufacturing of Tris, Bio Active Grade is performed at BioSpectra's Bangor, PA facility. The packaging, storage, and shipping of Tris, Bio Active Grade occurs in BioSpectra's Bangor, PA Facility: 100 Majestic Way, Bangor, PA 18013.

Current Product Number	Historic Product Number
TRIS-2220	TR2220
TRIS-2250	TR2250
TRIS-2255	TR2255
TRIS-2256	TR2256

For further information, please contact info@biospectra.us

Cassie Baun

Compliance Specialist