

SUPPLY CHAIN STATEMENT

L-Arginine Hydrochloride

L-Arginine Hydrochloride, Bio Qualified Grade is manufactured and packaged by BioSpectra's approved supplier in accordance with the approved supplier's ISO 9001:2015 certified management system. Additional blending and packaging of L-Arginine Hydrochloride, Bio Qualified Grade is performed at BioSpectra's Bangor, PA Facility: 100 Majestic Way, Bangor, PA 18013. The primary storage and shipment location for L-Arginine Hydrochloride, Bio Qualified Grade is BioSpectra's Supply Chain Center: 51 North 3rd Street, Stroudsburg, PA 18360. BioSpectra may additionally utilize the Bangor, PA Facility as a qualified secondary storage and shipment location, as needed.

Current Product Number
LARH-6350

For further information, please contact info@biospectra.us



Cassie Baun

Senior Compliance Specialist