

RESIDUAL SOLVENTS STATEMENT

L-Arginine Hydrochloride

L-Arginine Hydrochloride, Bio Qualified Grade is manufactured using the solvent Ethyl Acetate, and complies with the ICH Q3C Guideline for Residual Solvents required concentration limit of 5000ppm for a Class 3 residual solvent. BioSpectra has evaluated the material supply through the Supplier Qualification Program.

Current Product Number
LARH-6350

For further information, please contact info@biospectra.us



Cassie Baun

Senior Compliance Specialist