

DEXTRAN SULFATE SODIUM SALT, MW 8,000

LOW BIOBURDEN, GMP

CAS #: 9011-18-1

Formula: $(C_6H_7Na_3O_{14}S_3)_n$

F.W.: 8,000 g/mol

DXSE-4252

BIO PHARMA GRADE

ANALYSIS		SPECIFICATIONS
Appearance and Color		Off White to Light Yellow Powder
Chloride		≤ 1000 ppm
Absorbance (20%)	360 nm	≤ 0.9 OD Unit
Enzymes	DNase – NICKase	None Detected
	DNase – Exonuclease	None Detected
	Protease	None Detected
Free Inorganic Sulfate		$\leq 0.2\%$
Glucose		35 – 48%
Identification, Colorimetric		Passes Test
Insoluble Iron		$\leq 2.0\%$
Loss on Drying		$\leq 10.0\%$
Microbial Content	TAMC	≤ 10 CFU/g
	TYMC	≤ 10 CFU/g
pH (1%)		5.0 – 7.5
Residue on Ignition		35 – 50%
Pyridine		$\leq 2.0\%$
Specific Rotation, 20°C		+75° to +105°
Sulfur		17.0 – 20.0%
Specific Viscosity (In 1.0M NaCl, 25°C)		0.018 – 0.032

Country of Origin: USA

Standard Shelf Life Policy

Each Certificate of Analysis (CofA) will contain a 2-year retest/recertification date supported by a 3-year ICH Q1 Stability Study (if one is completed).

Package Sizes

1kg, 5kg, 10kg, 25kg, 50kg

Intended for Use in Biopharmaceutical & Biotechnological Applications and Products

This is a high purity, reagent grade product, purified under full cGMP conditions. Dextran Sulfate Sodium Salt, MW 8,000, is a polyanionic derivative of dextran, produced by the esterification of the glucose polymer (glucan polysaccharide) with chlorosulphonic acid. Dextran fractions are characterized by their average MW and MW distribution, with predominantly 1-6 glycosidic bonds and 5% or less of 1-3 branching. Dextran Sulfate Sodium Salt, MW 8,000 is used in the solubilization and purification process of protein molecules intended for use in a final drug product. Dextran is neutral in pH, soluble in water, easily filtered and biodegradable.

General Product Description

- Appears as an off-white to light yellow powder
- Manufactured under an IPEC Quality Managed cGMP System
- Manufactured in an enzyme free, hormone free and animal free environment
- Has no major food allergens (as defined by FDA and WHO)
- The final product nor its raw materials are not derived nor come into contact with animal parts, animal products, and/or animal byproducts or derivatives.
- Is not subject to genetic modification
- Synonyms: Dextran Hydrogen Sulfate, Sodium Salt
- Visit the product page on our website (www.biospectra.us) for additional information, supporting regulatory documents, and CofAs.

Storage and Shipping Conditions

Refer to SDS.

Standard Lead Time

1-2 weeks

This is not considered a controlled document. We are not responsible for any errors or omissions.