

12/19/2025

Version 1.1

## 50KG FIBER DRUM SPECIFICATIONS

BSI COMPONENT NUMBER: COMP-0069



26 Gallon / 50kg Fiber Drum	
Test	Specifications
Appearance and Color	Brown fiberboard drum, lid, and bottom with no water stains. Metal Chime and Lock Ring with no corrosion.
Outside Height with Lid (Nominal)	33"
Inside Top Diameter (Nominal)	15.5"
Stamp Markings	Minimum Test = 600 psi Net Wt. Limit = 150 pounds
Capacity	26-gallons

For further information, please contact [marketing@biospectra.us](mailto:marketing@biospectra.us).

**Connie Augustine | Director**  
**Business Integration & Retention**  
**BioSpectra Inc.**

*This is not considered a controlled document. We are not responsible for any errors or omissions, and the user is responsible for any decisions based on the information herein.*