



100 Majestic Way, Bangor, PA 18013 / www.biospectra.us

|                      |  |                    |                       |
|----------------------|--|--------------------|-----------------------|
| Effective Date:      | 05-Feb-2025                            | 05-Feb-2028        | : Date of Next Review |
| Prepared By:         | Taylor Yurick                          | BSI-COA-0225 v.3.1 | : Supersedes          |
| QA/QC Approval:      | Jessica DeMaio                         | Hannah Kuchmas     | : Management Approval |
| Reason for Revision: | See Revision History in MasterControl. |                    |                       |

**CERTIFICATE OF ANALYSIS**  
**TRIS HYDROCHLORIDE**  
**BIO EXCIPIENT GRADE / THCL-3259**  
**LOT: THCL-S03-1125-0164**

$\text{NH}_2\text{C}(\text{CH}_2\text{OH})_3 \cdot \text{HCl}$   $\wedge$  F.W. 157.60 g/mol.  $\wedge$  CAS# 1185-53-1

Manufacturing Date: 11/19/25      Expiration Date: 11/30/28

Manufacturing Site: 1474 Rockdale Lane, Stroudsburg, PA 18360

Packaging Site: 100 Majestic Way, Bangor PA, 18013

| ANALYSIS                 | SPECIFICATION           | TEST RESULT   |
|--------------------------|-------------------------|---------------|
| Absorbance (1M)          | 260 nm $\leq 0.06$ a.u. | $< 0.06$ a.u. |
|                          | 280 nm $\leq 0.06$ a.u. | $< 0.06$ a.u. |
|                          | 400 nm $\leq 0.01$ a.u. | $< 0.01$ a.u. |
| Appearance and Color     | White / Crystals        | Passes Test   |
| Assay, Dried             | 99.5 – 101.0%           | 99.8%         |
| Bioburden                | $\leq 100$ CFU/g        | $< 10$ CFU/g  |
| Endotoxin                | $\leq 2.5$ EU/g         | $< 1.9$ EU/g  |
| Enzymes                  | DNase                   | None Detected |
|                          | RNase                   | None Detected |
|                          | Protease                | None Detected |
| Heavy Metals             | 2 ppm max.              | $< 2$ ppm     |
| Identification           | (IR)                    | Passes Test   |
|                          | (Chloride)              | Passes Test   |
| Insoluble Matter         | 0.001% max.             | $< 0.001\%$   |
| Loss on Drying @ 105°C   | $\leq 0.5\%$            | 0.1%          |
| Melting Range            | 150 – 152 °C            | 151-152°C     |
| pH (1% Aqueous Solution) | 4.0 – 5.0               | 4.7           |
| pH (0.5M) @ 25°C         | 3.5 – 5.0               | 4.3           |
| pK <sub>a</sub>          | 8.0 – 8.4               | 8.2           |
| Residue on Ignition      | 0.1% max.               | $< 0.1\%$     |

The information contained herein is the confidential property of BioSpectra. The recipient is responsible for its safe-keeping and the prevention of unauthorized appropriation, use, disclosure and copying.

| ANALYSIS             | SPECIFICATION  | TEST RESULT  |
|----------------------|----------------|--------------|
| Solubility 35%       | Passes Test    | Passes Test  |
| Sulfated Ash (EP)    | $\leq 300$ ppm | $< 100$ ppm  |
| Arsenic (As)         | 1 ppm max.     | $< 0.45$ ppm |
| Cadmium (Cd)         | 1 ppm max.     | $< 0.30$ ppm |
| Calcium (Ca)         | 1 ppm max.     | $< 0.60$ ppm |
| Trace Metals         | Copper (Cu)    | $< 0.15$ ppm |
|                      | Iron (Fe)      | $< 0.30$ ppm |
|                      | Lead (Pb)      | $< 0.30$ ppm |
|                      | Magnesium (Mg) | $< 0.60$ ppm |
| Water (Karl Fischer) | 0.5% max.      | 0.3 %        |

COUNTRY OF ORIGIN: U.S.A.

TEST METHOD REFERENCE: DCN: BSI-ATM-0002

INTENDED USE: Material represented by this Certificate of Analysis is suitable for use as an excipient. It is manufactured in accordance with the ICH Q7 Good Manufacturing Practice Guide. The material represented by this Certificate of Analysis is not suitable to be used as an Active Pharmaceutical Ingredient, Drug Product or Household Item.

RESIDUAL SOLVENTS: Based on the manufacturing process and the controlled handling, storage and analysis of this product, this product complies with the requirements and specifications listed in the current USP method <467> Tables 1, 2, 3, or 4.

Prepared by: Shail McCall Date: 12/17/25 Job Title: GATEch III

Reviewed by: Joan Bingham Date: 12/17/25 Job Title: QA Supervisor