

Product		Description 1	Product Status	Product size	Product UOM	Package size	Package UOM	Package type	Liner size	Liner type	Desiccant	Interior cable ties	Exterior locks/seals
BTRI-3250-92	Bis Tris, 100 G	Active	100.00 G		250.00 ML			White HDPE bottle, blue lid					
BTRI-3250-93	Bis Tris, 500 G	Active	500.00 G		1000.00 ML			White HDPE bottle, blue lid					
BTRI-4250-92	Bis-Tris Biotech, 100 G	Active	100.00 G		250.00 ML			White HDPE bottle, blue lid					
BTRI-4250-93	Bis-Tris Biotech, 500 G	Active	500.00 G		1000.00 ML			White HDPE bottle, blue lid					
CSMH-3250-91	Cysteamine HCl (2-MEA), 50 G	In development	50.00 G		0.60 GAL			White HDPE round pail, blue lid	(2x) 9 x 15	LDPE	4 Unit	(2x) Large plastic	Blue numbered BSI metal cable lock
CSMH-3250-92	Cysteamine HCl (2-MEA), 100 G	In development	100.00 G		0.60 GAL			White HDPE round pail, blue lid	(2x) 9 x 15	LDPE	4 Unit	(2x) Large plastic	Blue numbered BSI metal cable lock
CSMH-8250-92	Cysteamine HCl (2-MEA), 100 G	In development	100.00 G		0.60 GAL			White HDPE round pail, blue lid	9 x 15	LDPE	4 Unit	(2x) Large plastic	Blue numbered BSI metal cable lock
DXSE-4250-92	Dextran Sulfate, 8000, GMP, 100g	Active	100.00 G		0.60 GAL			White HDPE round pail, blue lid	(2x) 9 x 15	LDPE	4 Unit	Large plastic	Blue numbered BSI metal cable lock
DXSE-4250-93	Dextran Sulfate, 8000, GMP, 20g	Active	20.00 G		0.60 GAL			White HDPE round pail, blue lid	(2x) 9 x 15	LDPE	4 Unit	(2x) Large plastic	Blue numbered BSI metal cable lock
DXSE-4250-94	Dextran Sulfate, 8000, GMP, 500g	Active	500.00 G		0.60 GAL			White HDPE round pail, blue lid	(2x) 9 x 15	LDPE	4 Unit	(2x) Large plastic	Blue numbered BSI metal cable lock
DXSE-5201-92	Dextran Sulfate, 8000, 100g	Active	100.00 G		0.60 GAL			White HDPE round pail, blue lid	(2x) 9 x 15	LDPE	4 Unit	(2x) Large plastic	Blue numbered BSI metal cable lock
GHCL-3202-93	Guanidine HCl, 500 G	Active	500.00 G		1000.00 ML			White HDPE bottle, blue lid					
GHCL-3220-92	Guanidine HCl, 100 G	Active	100.00 G		250.00 ML			White HDPE bottle, blue lid					
GHCL-3220-93	Guanidine HCl, 500 G	Active	500.00 G		1000.00 ML			White HDPE bottle, blue lid					
GHCL-3224-93	Guanidine HCl, 500 G	Active	500.00 G		1000.00 ML			White HDPE bottle, blue lid					
GHCL-3225-92	Guanidine HCl, 100G	Active	100.00 G		250.00 ML			White HDPE bottle, blue lid					
GHCL-3253-93	Guanidine HCl, 500 G	Active	500.00 G		1000.00 ML			White HDPE bottle, blue lid					
GHCL-4223-92	Guanidine Hydrochloride Biotech, 100 G	Active	100.00 G		250.00 ML			White HDPE bottle, blue lid					
GHCL-4223-93	Guanidine Hydrochloride Biotech, 500 G	Active	500.00 G		1000.00 ML			White HDPE bottle, blue lid					
GHCL-7201-93	Guanidine HCl, 500 G	In development	500.00 G		1000.00 ML			White HDPE bottle, blue lid					
GHCL-7202-93	Guanidine HCl, 500 G	Active	500.00 G		1000.00 ML			White HDPE bottle, blue lid					
GHCL-7203-93	Guanidine HCl, 300 G	Active	300.00 G		1000.00 ML			White HDPE bottle, blue lid					
GTHI-3220-92	Guanidine Thiocyanate, 100 G	Active	100.00 G		250.00 ML								

GTHI-3220-95	Guanidine Thiocyanate, 10 G	In development	10.00 G	60.00 ML	
GTHI-4221-92	Guanidine Thiocyanate Biotech, 100 G	Active	100.00 G	250.00 ML	White HDPE bottle, blue lid
GTHI-4221-93	Guanidine Thiocyanate Biotech, 500 G	Active	500.00 G	1000.00 ML	White HDPE bottle, blue lid
GTHI-5221-92	Guanidine Thiocyanate, 100 G	Active	100.00 G	250.00 ML	
HEPE-3220-91	HEPES, 250 G	Active	250.00 G	1000.00 ML	White HDPE bottle, blue lid
HEPE-3220-92	HEPES, 100 G	Active	100.00 G	250.00 ML	White HDPE bottle, blue lid
HEPE-3220-93	HEPES, 500 G	Active	500.00 G	1000.00 ML	White HDPE bottle, blue lid
HEPE-3251-92	HEPES, 100 G	Active	100.00 G	250.00 ML	White HDPE bottle, blue lid
HEPE-3320-92	HEPES, 100 G	In development	100.00 G	250.00 ML	White HDPE bottle, blue lid
HEPE-3351-92	HEPES, 100 G	In development	100.00 G	250.00 ML	White HDPE bottle, blue lid
HEPE-4250-92	HEPES Biotech, 100 G	Active	100.00 G	250.00 ML	White HDPE bottle, blue lid
HEPE-4250-93	HEPES Biotech, 500 G	Active	500.00 G	1000.00 ML	White HDPE bottle, blue lid
LARH-6250-92	Arginine (L) Hydrochloride Biotech, 100 G	In development	100.00 G	250.00 ML	White HDPE bottle, blue lid
LARH-6250-93	Arginine (L) Hydrochloride Biotech, 500 G	In development	500.00 G	1000.00 ML	White HDPE bottle, blue lid
LARH-6350-92	Arginine (L) Hydrochloride Biotech, 100 G	Active	100.00 G	250.00 ML	White HDPE bottle, blue lid
LARH-6350-93	Arginine (L) Hydrochloride Biotech, 500 G	Active	500.00 G	1000.00 ML	White HDPE bottle, blue lid
LCYS-6250-92	Cystine (L) Hydrochloride Biotech, 100 G	Active	100.00 G	250.00 ML	White HDPE bottle, blue lid
LCYS-6250-93	Cystine (L) Hydrochloride Biotech, 500 G	Active	500.00 G	1000.00 ML	White HDPE bottle, blue lid
LCYS-6350-92	Cystine (L) Hydrochloride Biotech, 100 G	Active	100.00 G	250.00 ML	White HDPE bottle, blue lid
LCYS-6350-93	Cystine (L) Hydrochloride Biotech, 500 G	Active	500.00 G	1000.00 ML	White HDPE bottle, blue lid
LGLM-6201-92	Glutamine (L) Biotech, 100 G	Active	100.00 G	250.00 ML	White HDPE bottle, blue lid
LGLM-6201-93	Glutamine (L) Biotech, 500 G	Active	500.00 G	1000.00 ML	White HDPE bottle, blue lid
LGLM-6301-92	Glutamine (L) Biotech, 100 G	In development	100.00 G	250.00 ML	White HDPE bottle, blue lid
LGLM-6301-93	Glutamine (L) Biotech, 500 G	In development	500.00 G	1000.00 ML	White HDPE bottle, blue lid
LHMM-6250-92	Histidine (L) MonoHCl Monohydrate Biotech, 100 G	Active	100.00 G	250.00 ML	White HDPE bottle, blue lid
LHMM-6250-93	Histidine (L) MonoHCl Monohydrate Biotech, 500 G	Active	500.00 G	1000.00 ML	White HDPE bottle, blue lid
LHMM-6350-92	Histidine (L) MonoHCl Monohydrate Biotech, 100 G	In development	100.00 G	250.00 ML	White HDPE bottle, blue lid
LHMM-6350-93	Histidine (L) MonoHCl Monohydrate Biotech, 500 G	In development	500.00 G	1000.00 ML	White HDPE bottle, blue lid

MESH-3201-92	MES Hydrate, 100 G	Active	100.00 G	250.00 ML	White HDPE bottle, blue lid
MESH-4250-92	MES Hydrate Biotech, 100 G	Active	100.00 G	250.00 ML	White HDPE bottle, blue lid
MESH-4250-93	MES Hydrate Biotech, 500 G	Active	500.00 G	1000.00 ML	White HDPE bottle, blue lid
MESH-9250-91	MES Hydrate, 50 G	In development	50.00 G	250.00 ML	White HDPE bottle, blue lid
MESM-3220-91	MES Monohydrate, 50G	Active	50.00 G	250.00 ML	White HDPE bottle, blue lid
MESM-3220-92	MES, Monohydrate, 100G	Active	100.00 G	250.00 ML	White HDPE bottle, blue lid
MESM-3220-93	MES, Monohydrate, 500G	Active	500.00 G	1000.00 ML	White HDPE bottle, blue lid
MESM-3221-93	MES, Monohydrate, 500G	Active	500.00 G	1000.00 ML	White HDPE bottle, blue lid
MESM-3222-91	MES Monohydrate, 50G	Active	50.00 G	250.00 ML	White HDPE bottle, blue lid
MESM-3222-92	MES, Monohydrate, 100G	Active	100.00 G	250.00 ML	White HDPE bottle, blue lid
MESM-3222-93	MES, Monohydrate, 500G	In development	500.00 G	1000.00 ML	White HDPE bottle, blue lid
MESM-3250-92	MES, Monohydrate, 100G	Active	100.00 G	250.00 ML	White HDPE bottle, blue lid
MESM-3250-93	MES, Monohydrate, 500G	Active	500.00 G	1000.00 ML	White HDPE bottle, blue lid
MESM-3250-95	MES, Monohydrate, 10G	In development	10.00 G	250.00 ML	White HDPE bottle, blue lid
MESM-3251-92	MES, Monohydrate, 100G	Active	100.00 G	250.00 ML	White HDPE bottle, blue lid
MESM-3320-92	MES, Monohydrate, 100G	In development	100.00 G	250.00 ML	White HDPE bottle, blue lid
MESM-3321-93	MES, Monohydrate, 500G	In development	500.00 G	1000.00 ML	White HDPE bottle, blue lid
MESM-3350-92	MES, Monohydrate, 100G	In development	100.00 G	250.00 ML	White HDPE bottle, blue lid
MESM-4221-92	MES, Monohydrate Biotech, 100 G	Active	100.00 G	250.00 ML	White HDPE bottle, blue lid
MESM-4221-93	MES, Monohydrate Biotech, 500 G	Active	500.00 G	1000.00 ML	White HDPE bottle, blue lid
MESM-5220-93	MES, Monohydrate, 500 G	Active	500.00 G	1000.00 ML	White HDPE bottle, blue lid
MESN-5220-93	MES Sodium, 300 G	Active	300.00 G	1000.00 ML	White HDPE bottle, blue lid
MOPS-3220-92	MOPS, 100 G	Active	100.00 G	250.00 ML	White HDPE bottle, blue lid
MOPS-3220-93	MOPS, 500 G	Active	500.00 G	1000.00 ML	White HDPE bottle, blue lid
MOPS-4221-92	MOPS Biotech, 100 G	Active	100.00 G	250.00 ML	White HDPE bottle, blue lid
MOPS-4221-93	MOPS Biotech, 500 G	Active	500.00 G	1000.00 ML	White HDPE bottle, blue lid
NDEC-3201-92	Sodium Decanoate, 100 gr	Active	100.00 G	1000.00 ML	White HDPE bottle, blue lid

NDEC-3220-91	Sodium Decanoate, 500 G	Active	500.00 G	0.60 GAL	White HDPE round pail, blue lid	9 x 15	LDPE	Large plastic	Blue numbered BSI metal cable lock
NDEC-3220-92	Sodium Decanoate, 100 G	Active	100.00 G	1000.00 ML	White HDPE bottle, blue lid				
NDEC-4201-92	Sodium Decanoate Biotech, 100 G	Active	100.00 G	1000.00 ML	White HDPE bottle, blue lid				
NDEC-4201-93	Sodium Decanoate Biotech, 500 G	Active	500.00 G	0.60 GAL	White HDPE round pail, blue lid	9 x 15	LDPE	Large plastic	Blue numbered BSI metal cable lock
THCL-3203-92	Tris Hydrochloride, 100 G	Active	100.00 G	250.00 ML	White HDPE bottle, blue lid				
THCL-3203-93	Tris Hydrochloride, 500 G	Active	500.00 G	1000.00 ML	White HDPE bottle, blue lid				
THCL-3220-92	Tris Hydrochloride, 100 G	Active	100.00 G	250.00 ML	White HDPE bottle, blue lid				
THCL-3220-93	Tris Hydrochloride, 500 G	Active	500.00 G	1000.00 ML	White HDPE bottle, blue lid				
THCL-3221-92	Tris Hydrochloride, 100 G	Active	100.00 G	250.00 ML	White HDPE bottle, blue lid				
THCL-3250-91	Tris Hydrochloride, 50 G	Active	50.00 G	250.00 ML	White HDPE bottle, blue lid				
THCL-3250-92	Tris Hydrochloride, 100 G	Active	100.00 G	250.00 ML	White HDPE bottle, blue lid				
THCL-3250-93	Tris Hydrochloride, 500 G	Active	500.00 G	1000.00 ML	White HDPE bottle, blue lid				
THCL-3253-92	Tris Hydrochloride, 100 G	Active	100.00 G	250.00 ML	White HDPE bottle, blue lid				
THCL-3253-93	Tris Hydrochloride, 500 G	Active	500.00 G	1000.00 ML	White HDPE bottle, blue lid				
THCL-3254-92	Tris Hydrochloride, 100 G	Active	100.00 G	250.00 ML	White HDPE bottle, blue lid				
THCL-3254-93	Tris Hydrochloride, 500 G	Active	500.00 G	1000.00 ML	White HDPE bottle, blue lid				
THCL-3257-93	Tris Hydrochloride, 500 G	Active	500.00 G	1000.00 ML	White HDPE bottle, blue lid				
THCL-3259-91	Tris Hydrochloride, 250 G	Active	250.00 G	1000.00 ML	White HDPE bottle, blue lid				
THCL-3259-92	Tris Hydrochloride, 100 G	Active	100.00 G	250.00 ML	White HDPE bottle, blue lid				
THCL-3259-93	Tris Hydrochloride, 500 G	Active	500.00 G	1000.00 ML	White HDPE bottle, blue lid				
THCL-3260-92	Tris Hydrochloride, 100 G	Active	100.00 G	250.00 ML	White HDPE bottle, blue lid				
THCL-4251-92	Tris Hydrochloride Biotech, 100 G	Active	100.00 G	250.00 ML	White HDPE bottle, blue lid				
THCL-4251-93	Tris Hydrochloride Biotech, 500 G	Active	500.00 G	1000.00 ML	White HDPE bottle, blue lid				
THCL-7201-93	Tris Hydrochloride, 500G	Active	500.00 G	1000.00 ML	White HDPE bottle, blue lid				
THCL-7205-93	Tris Hydrochloride, 500G	Active	500.00 G	1000.00 ML	White HDPE bottle, blue lid				
TRIS-1201-91	Tris/Tromethamine, 250 G	Active	250.00 G	1000.00 ML	White HDPE bottle, blue lid				

TRIS-1201-92	Tris/Tromethamine, 100 G	Active	100.00 G	250.00 ML	White HDPE bottle, blue lid			
TRIS-1201-93	Tris/Tromethamine, 500 G	Active	500.00 G	1000.00 ML	White HDPE bottle, blue lid			
TRIS-1220-91	Tris/Tromethamine, 250 G	Active	250.00 G	1000.00 ML	White HDPE bottle, blue lid			
TRIS-1220-92	Tris/Tromethamine, 100 G	Active	100.00 G	250.00 ML	White HDPE bottle, blue lid			
TRIS-1220-93	Tris/Tromethamine, 500 G	Active	500.00 G	1000.00 ML	White HDPE bottle, blue lid			
TRIS-2255-91	Tris/Tromethamine, 250 G	Active	250.00 G	1000.00 ML	White HDPE bottle, blue lid			
TRIS-2255-92	Tris/Tromethamine, 100 G	Active	100.00 G	250.00 ML	White HDPE bottle, blue lid			
TRIS-2255-93	Tris/Tromethamine, 500 G	Active	500.00 G	1000.00 ML	White HDPE bottle, blue lid			
TRIS-2257-91	Tris/Tromethamine, 250 G	Active	250.00 G	1000.00 ML	White HDPE bottle, blue lid			
TRIS-2257-92	Tris/Tromethamine, 100 G	Active	100.00 G	250.00 ML	White HDPE bottle, blue lid			
TRIS-2257-93	Tris/Tromethamine, 500 G	Active	500.00 G	1000.00 ML	White HDPE bottle, blue lid			
TRIS-3201-92	Tris/Tromethamine, 100 G	Active	100.00 G	250.00 ML	White HDPE bottle, blue lid			
TRIS-3201-93	Tris/Tromethamine, 500 G	Active	500.00 G	1000.00 ML	White HDPE bottle, blue lid			
TRIS-3220-92	Tris/Tromethamine, 100 G	Active	100.00 G	250.00 ML	White HDPE bottle, blue lid			
TRIS-3220-93	Tris/Tromethamine, 500 G	Active	500.00 G	1000.00 ML	White HDPE bottle, blue lid			
TRIS-3251-92	Tris/Tromethamine, 100 G	Active	100.00 G	250.00 ML	White HDPE bottle, blue lid			
TRIS-3251-95	Tris/Tromethamine, 80 G	Active	80.00 G	1000.00 ML	White HDPE bottle, blue lid	9 x 15	LDPE	Large plastic
TRIS-3254-92	Tris/Tromethamine, 100 G	Active	100.00 G	250.00 ML	White HDPE bottle, blue lid			
TRIS-3254-93	Tris/Tromethamine, 500 G	Active	500.00 G	1000.00 ML	White HDPE bottle, blue lid			
TRIS-3255-91	Tris/Tromethamine, 50 G	Active	50.00 G	250.00 ML	White HDPE bottle, blue lid			
TRIS-3255-92	Tris/Tromethamine, 100 G	Active	100.00 G	250.00 ML	White HDPE bottle, blue lid			
TRIS-3255-93	Tris/Tromethamine, 500 G	Active	500.00 G	1000.00 ML	White HDPE bottle, blue lid			
TRIS-3255-94	Tris/Tromethamine, 250 G	Active	250.00 G	1000.00 ML	White HDPE bottle, blue lid			
TRIS-4223-92	Tris (Tromethamine) Biotech, 100 G	Active	100.00 G	250.00 ML	White HDPE bottle, blue lid			
TRIS-4223-93	Tris (Tromethamine) Biotech, 500 G	Active	500.00 G	1000.00 ML	White HDPE bottle, blue lid			
TRIS-5223-89	Tris/Tromethamine, 200 G	Active	200.00 G	1000.00 ML	White HDPE bottle, blue lid			
TRIS-7201-93	Tris/Tromethamine, 500 G	Active	500.00 G	1000.00 ML	White HDPE bottle, blue lid			

URAC-4202-91	Uracil, 50 G	Active	50.00 G	250.00 ML	White HDPE bottle, blue lid
URAC-4202-92	Uracil, 100 G	Active	100.00 G	250.00 ML	White HDPE bottle, blue lid
URAC-4202-95	Uracil, 10g	Active	10.00 G	60.00 ML	White polypropylene bottle
URAC-4250-95	Uracil, 10g	Active	10.00 G	60.00 ML	White polypropylene bottle
UREA-3220-89	Urea, 200g	Active	200.00 G	1000.00 ML	White HDPE bottle, blue lid
UREA-3220-92	Urea, 100g	Active	100.00 G	250.00 ML	White HDPE bottle, blue lid
UREA-3220-93	Urea, 500 G	Active	500.00 G	1000.00 ML	White HDPE bottle, blue lid
UREA-3250-92	Urea, 100g	Active	100.00 G	250.00 ML	White HDPE bottle, blue lid
UREA-3250-93	Urea, 500 G	Active	500.00 G	1000.00 ML	White HDPE bottle, blue lid
UREA-4221-92	Urea Biotech, 100 G	Active	100.00 G	250.00 ML	White HDPE bottle, blue lid
UREA-4221-93	Urea Biotech, 500 G	Active	500.00 G	1000.00 ML	White HDPE bottle, blue lid
UREA-7201-93	Urea, 500 G	Active	500.00 G	1000.00 ML	White HDPE bottle, blue lid
URID-3250-91	Uridine, 250 G	Active	250.00 G	1000.00 ML	White HDPE bottle, blue lid
URID-3250-92	Uridine, 100 G	Active	100.00 G	250.00 ML	White HDPE bottle, blue lid
URID-3250-93	Uridine, 500 G	Active	500.00 G	1000.00 ML	White HDPE bottle, blue lid
URID-4250-92	Uridine Biotech, 100 G	Active	100.00 G	250.00 ML	White HDPE bottle, blue lid
URID-4250-93	Uridine Biotech, 500 G	Active	500.00 G	1000.00 ML	White HDPE bottle, blue lid
URID-4250-95	Uridine Biotech, 50 G	Active	100.00 G	250.00 ML	White HDPE bottle, blue lid