

TRIS HYDROCHLORIDE

LBLE

CAS #: 1185-53-1

Formula: $\text{NH}_2\text{C}(\text{CH}_2\text{OH})_3 \cdot \text{HCl}$

F.W.: 157.60 g/mol

THCL-6251
BIO QUALIFIED GRADE

ANALYSIS		SPECIFICATIONS
Absorbance (1.0M)	260 nm	< = 0.05 a.u.
	280 nm	< = 0.05 a.u.
	400 nm	< = 0.01 a.u.
Appearance		White Crystalline
Assay, Dried Basis		> = 99.5%
Assay, by Titration (0.1N NaOH)		99.5 – 101.0%
Chloride		Passes Test
Cytotoxicity		Passes Test
Endotoxin		< = 2.5 EU/g
Enzymes	RNase	Passes Test
	DNase	Passes Test
	Exonuclease	Passes Test
	Protease	Passes Test
Heavy Metal		< = 10 ppm
Identification, IR		Conforms to Reference Standard
Insoluble Matter		< = 10 ppm
Loss on Drying		< = 0.05%
Melting Point		149 – 153C
Microbial Content	TAMC	< 50 CFU/g
	TYMC	< 50 CFU/g
pH (1%)		4.2 – 4.9
Residue on Ignition		< = 0.1%
Residual Solvents	Methanol	< = 300 ppm
	Nitromethane	< = 50 ppm
	Ethanol	< = 5000 ppm
	Ethylmethylketon	< = 5000 ppm
Solubility (1.0M)		Clear, Colorless Solution

ANALYSIS		SPECIFICATIONS
Trace Metals	Aluminum (Al)	< = 5 ppm
	Arsenic (As)	< = 1 ppm
	Barium (Ba)	< = 5 ppm
	Bismuth (Bi)	< = 5 ppm
	Cadmium (Cd)	< = 5 ppm
	Calcium (Ca)	< = 5 ppm
	Chromium (Cr)	< = 5 ppm
	Cobalt (Co)	< = 5 ppm
	Copper (Cu)	< = 1 ppm
	Iron (Fe)	< = 1 ppm
	Lead (Pb)	< = 5 ppm
	Lithium (Li)	< = 5 ppm
	Magnesium (Mg)	< = 5 ppm
	Manganese (Mn)	< = 1 ppm
	Mercury (Hg)	< = 5 ppm
	Molybdenum (Mo)	< = 5 ppm
	Nickel (Ni)	< = 5 ppm
	Potassium (K)	< = 50 ppm
	Silver (Ag)	< = 5 ppm
	Sodium (Na)	< = 50 ppm
	Strontium (Sr)	< = 5 ppm
	Tin (Sn)	< = 5 ppm
	Zinc (Zn)	< = 5 ppm
Water, KF		< = 0.5%

Refer to page 2 for additional product details.

This is not considered a controlled document. We are not responsible for any errors or omissions, and the user is responsible for any decisions based on the information herein.

General Product Overview

Tris Hydrochloride is a stabilizing buffer in biological applications such as electrochromatography, UV analysis and HPLC. It is used to adjust and stabilize the pH ranges for gels used in electrophoresis applications. Tris Hydrochloride is extensively used as a biological buffer or a component of buffer solutions.

Industry Application

Suitable for use as a cGMP chemical in pharmaceutical manufacturing processes and biological/biotech process applications from R&D through scale cGMP production.

Key Product Features

- Appears as a white crystalline product
- Manufactured in accordance with the ISO Certified Quality Management System
- Manufactured in an enzyme free, hormone free and animal free environment
- Contains no known major food allergens (as defined by FDA and WHO)
- The final product and its raw materials are not derived from nor come into contact with animal parts, animal products, and/or animal byproducts or derivatives.
- Is not subject to genetic modification
- Synonyms: 2-Amino-2-(Hydroxymethyl)-1,3-Propanediol Hydrochloride; Tri(Hydroxymethyl) Aminomethane Hydrochloride.

Storage and Shipping Conditions

Refer to SDS.

Standard Shelf-Life Policy

Please inquire for information regarding shelf-life.

Package Sizes

1kg, 5kg, 10kg, 25kg, 50kg

[Click here to view SDS, CoAs and other supporting regulatory documents on our website.](#)

This is not considered a controlled document. We are not responsible for any errors or omissions, and the user is responsible for any decisions based on the information herein.