

BIOSPECTRA

100 Majestic Way, Bangor, PA 18013 / www.biospectra.us

Effective Date:	31-Jan-2025	31-Jan-2028	: Date of Next Review
Prepared By:	Taylor Yurick	BSI-COA-0197 v.4.3	: Supersedes
QA/QC Approval:	Jaron Hughes	Carissa Albert	: Management Approval
Reason for Revision:	See Revision History in MasterControl.		

CERTIFICATE OF ANALYSIS

TREHALOSE, DIHYDRATE

BIO EXCIPIENT GRADE / TRED-3252

LOT: TRED-N02-0925-0026

C₁₂H₂₂O₁₁ · 2H₂O ★ F.W. 378.33 g/mol. ★ CAS# 6138-23-4

Manufacturing Date: 09/09/25 Retest Date: 09/30/27

Manufacturing Site: 100 Majestic Way, Bangor PA, 18013

Packaging Site: 100 Majestic Way, Bangor PA, 18013

Meets or Exceeds NF, EP, and JP Specifications

NF COMPENDIA

ANALYSIS	SPECIFICATION	TEST RESULT
Assay ¹	98.0 - 101.0% ³	99.6%
Chloride and Sulfate, <i>Chloride</i>	≤ 0.0125%	≤ 0.0125 %
Color and Clarity of Solution	A720 ≤ 0.050	< 0.003
	A420 – A720 ≤ 0.100	0.011
Endotoxins ²	≤ 0.3 EU/g ³	≤ 0.2 EU/g
Identification A ²	Conforms to Standard	Conforms to standard
Identification B ²	Passes Test	Passes Test
Identification C ²	Passes Test	Passes Test
	<i>Escherichia coli</i>	Absent/g
	<i>Salmonella species</i>	Absent/10g
Microbial Content ²	TAMC ≤ 50 CFU/g	< 10 CFU/g
	TYMC ≤ 20 CFU/g	< 10 CFU/g
Nitrogen Determination ²	≤ 0.005%	≤ 0.005%
Optical Rotation, Specific Rotation @ 20°C ²	+197° to +201°	+200°
pH @ 25°C ²	4.5 - 6.5	5.8
Related Substances ¹	Total Impurities with RRT <1.0 ≤ 0.5%	0.21%
	Total Impurities with RRT >1.0 ≤ 0.5%	< 0.01%
Residue on Ignition ²	≤ 0.1%	< 0.1%

The information contained herein is the confidential property of BioSpectra. The recipient is responsible for its safe-keeping, and the prevention of unauthorized appropriation, use, disclosure and copying.

ANALYSIS	SPECIFICATION	TEST RESULT
Soluble Starch ²	Passes Test	Passes Test
Chloride and Sulfate, <i>Sulfate</i>	≤ 0.0200%	< 0.0200%
Water Determination ²	9.0% to 11.0%	9.9%

EP COMPENDIA

ANALYSIS	SPECIFICATION	TEST RESULT
Assay ¹	98.0 – 101.0% ³	99.6%
Appearance of Solution	Clear, colorless	Clear, colorless
Chlorides	≤ 0.0125%	< 0.0125%
Endotoxins ²	≤ 0.3 EU/g ³	< 0.2 EU/g
Identification A ²	Conforms to Standard	Conforms to standard
Identification B ²	Passes Test	Passes Test
Identification C ²	Passes Test	Passes Test
Related Substances ¹	Impurity A	< 0.10%
	Impurity B	< 0.10%
	Unspecified Impurities	0.2%
	Total Impurities	0.2%
Microbial Content ²	<i>Escherichia coli</i>	Absent/g
	<i>Salmonella species</i>	Absent/10g
	TAMC	< 10 CFU/g
	TYMC	< 10 CFU/g
pH @ 25°C ²	4.5 – 6.5	5.8
Soluble Starch ²	Passes Test	Passes Test
Specific, Optical Rotation @ 20°C ²	+197° to +201°	+200°
Sulfated Ash	≤ 0.1%	< 0.1%
Sulfate	≤ 0.0200%	< 0.0200%
Water ²	9.0% to 11.0%	9.9%

JP COMPENDIA

ANALYSIS	SPECIFICATION	TEST RESULT
Assay ¹	98.0 – 101.0%	99.6%
Chloride	≤ 0.018%	< 0.018%
Dextrin, Soluble Starch, Sulfite ²	Passes Test	Passes Test

ANALYSIS	SPECIFICATION	TEST RESULT
Heavy Metals (as Pb)	≤ 5 ppm	< 5 ppm
Identification 1 ²	Passes Test	Passes Test
Identification 2 ²	Passes Test	Passes Test
Identification 3 ²	Conforms to Standard	Conforms to standard
Nitrogen ²	≤ 0.005%	< 0.005%
Optical Rotation @ 20°C ²	+197° to +201°	+200°
pH @ 25°C ²	4.5 - 6.5	5.8
Residue on Ignition ²	≤ 0.1%	< 0.1%
Related Substances ¹	Total Impurities with RRT <1.0	0.21%
	Total Impurities with RRT >1.0	< 0.01%
Sulfate	≤ 0.024%	< 0.024%
Water ²	9.0% to 11.0%	9.9%

NON-COMPENDIAL ANALYSES

ANALYSIS	SPECIFICATION	TEST RESULT
Appearance and Color	White to Off White Crystalline Powder	White to Off White Crystalline Powder
Microbial Content	<i>Staphylococcus aureus</i>	Absent/g
	<i>Pseudomonas aeruginosa</i>	Absent/g
Residual Ethanol ¹	≤ 200 ppm	< 95ppm
Residual Isopropyl Alcohol ¹	≤ 250 ppm	< 130ppm
Residual Methanol ¹	≤ 50 ppm	< 25ppm
Trace Metals	Cadmium (Cd)	≤50 ppb
	Arsenic (As)	≤50 ppb
	Mercury (Hg)	≤50 ppb
	Lead (Pb)	≤50 ppb
	Nickel (Ni)	≤100 ppb
	Molybdenum (Mo)	≤100 ppb
	Copper (Cu)	≤100 ppb
	Chromium (Cr)	≤100 ppb
	Iron (Fe)	≤100 ppb
	Aluminum (Al)	≤100 ppb
	Zinc (Zn)	≤100 ppb

¹Alternate Validated Method²Analyses are Harmonized³Specifications is more stringent than Compendia Monograph

COUNTRY OF ORIGIN: U.S.A.

TEST METHOD REFERENCE: DCN: BSI-ATM-0027

INTENDED USE: Material represented by this Certificate of Analysis is suitable for use as an excipient. It is manufactured in accordance with the ICH Q7 Good Manufacturing Practice Guide. The material represented by this Certificate of Analysis is not suitable to be used as an Active Pharmaceutical Ingredient, Drug Product or Household Item.

Prepared by: David McCall Date: 10/15/26 Job Title: QA Tech III

Reviewed by: Car Albert Date: 10/19/26 Job Title: Senior Quality Manager