

SUPPLY CHAIN STATEMENT

Tris

Tris, Bio Qualified Grade is manufactured and packaged in accordance with an ISO 9001:2015 Certified Quality Management System. Additional packaging of Tris, Bio Qualified Grade is performed at BioSpectra's Majestic Facility. The primary storage and shipment location for Tris, Bio Qualified Grade is BioSpectra's McConnell Facility: 51 North 3rd Street, Stroudsburg, PA 18360. BioSpectra may additionally utilize the Majestic Facility: 100 Majestic Way, Bangor, PA 18013 as a qualified secondary storage and shipment location as needed.

Current Product Number
TRIS-6201
TRIS-6202
TRIS-6250
TRIS-6251
TRIS-6252

For further information, please contact Customer.Service@BioSpectra.us



Hailey Glissmann

Compliance Support Technician