


## PROCESS FLOW DIAGRAM

### L-Arginine Hydrochloride

L-Arginine Hydrochloride, Bio Qualified Grade is manufactured by BioSpectra's approved Supplier, and blended by BioSpectra through the following process.



For further information, please contact [info@biospectra.us](mailto:info@biospectra.us)

**Cassie Baun**

Senior Compliance Specialist

The information contained herein is the property of BioSpectra. The recipient is responsible for its safe-keeping and the prevention of unauthorized appropriation, use, disclosure, copying or alteration. Decision based on the use of this information is the sole responsibility and liability of the recipient. If you would like a controlled version of this document, please contact [info@biospectra.us](mailto:info@biospectra.us)