

GUANIDINE THIOCYANATE

BIOTECH, GMP

CAS #: 593-84-0

Formula: $C_2H_6N_4S$

F.W.: 118.16 g/mol

GTHI-4221

BIO PHARMA GRADE

ANALYSIS		SPECIFICATIONS
Absorbance (1.7M)	280 nm	≤ 0.300 a.u.
	300 nm	≤ 0.050 a.u.
	340 nm	≤ 0.030 a.u.
Appearance and Color		White / Crystals
Assay		99.0 – 101.0%
Enzymes	DNase	None Detected
	Protease	None Detected
	RNase	None Detected
Identification, IR		Conforms to Reference Standard
Loss on Drying		$\leq 0.5\%$
Melting Range		115 – 121°C
pH (5% solution)		5.0 – 7.0
Sulfated Ash		$\leq 0.5\%$
Sulfate		≤ 50 ppm
Trace Metals	Arsenic (As)	≤ 5 ppm
	Aluminum (Al)	≤ 5 ppm
	Barium (Ba)	≤ 5 ppm
	Bismuth (Bi)	≤ 5 ppm
	Calcium (Ca)	≤ 10 ppm
	Cadmium (Cd)	≤ 5 ppm
	Chromium (Cr)	≤ 5 ppm
	Cobalt (Co)	≤ 5 ppm
	Copper (Cu)	≤ 5 ppm
	Iron (Fe)	≤ 5 ppm
	Potassium (K)	≤ 50 ppm
	Lead (Pb)	≤ 5 ppm
	Lithium (Li)	≤ 5 ppm
	Magnesium (Mg)	≤ 5 ppm
	Manganese (Mn)	≤ 5 ppm
	Molybdenum (Mo)	≤ 5 ppm
	Sodium (Na)	≤ 5000 ppm
	Nickel (Ni)	≤ 5 ppm
	Strontium (Sr)	≤ 5 ppm
	Zinc (Zn)	≤ 5 ppm

Country of Origin: USA

Intended for Use in Biopharmaceutical & Biotechnological Applications and Products

This is a high purity, GMP, Reagent Grade product, purified under full cGMP conditions. Guanidine Thiocyanate is a chaotropic agent. As one of the strongest protein denaturants, it inactivates nucleases more than twice as fast as Guanidine Hydrochloride. Guanidine Thiocyanate is preferred for the purification of RNA because it dissociates the RNA into its nucleic acids and protein forms. Higher molarity solutions of Guanidine Thiocyanate irreversibly inactivate RNase. Guanidine Thiocyanate causes proteins to dissolve readily, thus disintegrating cellular structures.

General Product Description

- Appears as a white crystalline product
- Manufactured under an ICH-Q7 Quality Managed cGMP System
- Manufactured in an enzyme free, hormone free and animal free environment
- No known major food allergens (as defined by FDA and WHO)
- The final product and its raw materials are not derived from nor come into contact with animal parts, animal products, and/or animal byproducts/derivatives.
- Is not subject to genetic modification
- Synonyms: Guanidine Hydro Thiocyanate, Guanidine Isothiocyanate, Guanidine Mono Thiocyanate
- Visit the product page on our website (www.biospectra.us) for additional information, supporting regulatory documents, and CofAs.

Storage and Shipping Conditions

Refer to SDS.

Standard Shelf Life Policy

Each Certificate of Analysis will contain a 2-year retest/recertification date supported by a 3-year ICH Q1 Stability Study (if one is completed).

Package Sizes

100g, 500g, 1kg, 5kg, 10kg, 25kg, 50kg

Standard Lead Time

1-2 weeks

This is not considered a controlled document. We are not responsible for any errors or omissions.