

STORAGE AND SHIPPING CONDITIONS STATEMENT

Urea GMP

Urea, Bio Excipient Grade manufactured by BioSpectra should be stored as indicated in the Urea Safety Data Sheet, BSI-SDS-0003 v.5.4: "Preserve in well-closed containers. Store at 25°C, excursions permitted between 15°C and 30°C. Store in a tightly closed container. Store in a cool, dry, well-ventilated area away from incompatible substances."

Current Product Number
UREA-3220
UREA-3221
UREA-3222
UREA-3223
UREA-3224
UREA-3250

For further information, please contact Customer.Service@BioSpectra.us

Cassie Baun

Cassie Baun

Senior Compliance Specialist