

RESIDUAL SOLVENTS STATEMENT

Tris GMP

BioSpectra can state based on the manufacturing process and the controlled handling, storage, and analysis of this product that the Tris, Bio Excipient Grade, and Bio Pharma Grade manufactured by BioSpectra complies with the requirements and specifications listed in the current USP method <467> Tables 1, 2, 3, or 4. Tris is analyzed annually for residual solvents, with specifications for Methanol $\leq$ 3000ppm, and Nitromethane $\leq$ 50 ppm.

Current Product Number
TRIS-3220
TRIS-3251
TRIS-3252
TRIS-3254
TRIS-3255
TRIS-3257
TRIS-4224
TRIS-4252

For further information, please contact Customer.Service@BioSpectra.us



Cassie Baum

Senior Compliance Specialist