

2/10/2026

Version 1.2

12.0 GAL / 25KG PAIL SPECIFICATIONS

BSI COMPONENT NUMBER: COMP-0114



12 Gallon / 25kg HDPE Round Screw Top Pharma Pail and Lid	
Test	Specifications
Appearance and Color	White HDPE round pail with tab insert, blue tamper evident lid with gasket, report slot location (if applicable)
Inside Height	21.0" \pm 0.5"
Inside Top Diameter	14.0" \pm 0.5"
Outside Height, with Lid	21.5" \pm 0.5"
Outside Bottom Diameter	13.0" \pm 0.5"
UN Rating	1H2/Y40/S/YR with gasket 1H2/Y53/S/YR without gasket

For further information, please contact marketing@biospectra.us.



**Connie Augustine | Director
Business Integration & Retention
BioSpectra Inc.**

This is not considered a controlled document. We are not responsible for any errors or omissions, and the user is responsible for any decisions based on the information herein.