

# BIOBUFFER SOLUTIONS

100 Majestic Way, Bangor, PA 18013 / [www.biospectra.us](http://www.biospectra.us)

May 4<sup>th</sup>, 2026

Revision 1

## SUPPLY CHAIN STATEMENT

### L-Arginine Hydrochloride

L-Arginine Hydrochloride, Bio Ultra Grade is manufactured and packaged in accordance with an ISO 9001:2015 certified management system. Additional blending and packaging of L-Arginine Hydrochloride, Bio Ultra Grade is performed at BioSpectra's Majestic Facility: 100 Majestic Way, Bangor, PA 18013. The primary storage and shipment location for L-Arginine Hydrochloride, Bio Ultra Grade is BioSpectra's McConnell Facility: 51 North 3<sup>rd</sup> Street, Stroudsburg, PA 18360. BioSpectra may additionally utilize the Majestic Facility as a qualified secondary storage and shipment location, as needed.

Current Product Number
LARH-5350

For further information, please contact [Customer.Service@BioSpectra.us](mailto:Customer.Service@BioSpectra.us)



**Cassie Baun**

Senior Compliance Specialist

The information contained herein is the property of BioSpectra. The recipient is responsible for its safe-keeping and the prevention of unauthorized appropriation, use, disclosure, copying or alteration. Decision based on the use of this information is the sole responsibility and liability of the recipient. If you would like a controlled version of this document, please contact [Customer.Service@BioSpectra.us](mailto:Customer.Service@BioSpectra.us)

100 Majestic Way | Bangor, PA 18013 | [www.BioSpectra.us](http://www.BioSpectra.us) | 1.610.599.3400