

12/19/2025

Version 1.1

## 5.0 GAL / 10/12KG PAIL SPECIFICATIONS

BSI COMPONENT NUMBER: COMP-0112



5 Gallon / 10/12kg HDPE Round Screw Top Pharma Pail and Lid	
Test	Specifications
Appearance and Color	White HDPE round pail with hole(s) for tamper evident seal, blue tamper evident lid with gasket, red trigger, and slot on top of lid
Inside Height	15.0" $\pm$ 0.5"
Inside Top Diameter	11.25" $\pm$ 0.5"
Outside Height, with Lid	16.0" $\pm$ 0.5"
Outside Bottom Diameter	10.5" $\pm$ 0.5"
UN Rating	1H2/Y30/S/YR 1H2/X11.5/S/YR

For further information, please contact [marketing@biospectra.us](mailto:marketing@biospectra.us).

**Connie Augustine | Director  
Business Integration & Retention  
BioSpectra Inc.**

*This is not considered a controlled document. We are not responsible for any errors or omissions, and the user is responsible for any decisions based on the information herein.*