

RESIDUAL SOLVENTS STATEMENT

Tris

Tris, Bio Qualified Grade is manufactured using the Class 2 residual solvents Methanol and Nitromethane. Tris, Bio Qualified Grade complies with the specifications of ≤ 300 ppm Methanol, and ≤ 50 ppm Nitromethane.

Current Product Number
TRIS-6201
TRIS-6202
TRIS-6250
TRIS-6251
TRIS-6252

For further information, please contact Customer.Service@BioSpectra.us



Hailey Glissmann
Compliance Support Technician