

Product	Product Description	Product Status	Product size	Product UOM	Package size	Package UOM	Package type	Liner size	Liner type	Desiccant
BTRI-3250-01	Bis Tris, 1 KG	Active	1.00 KG		1500.00 ML		White HDPE bottle, blue lid			
BTRI-3250-05	Bis Tris, 5 KG	Active	5.00 KG		3.50 GAL		White HDPE round pail, blue lid	24 x 36	LDPE	
BTRI-3250-10	Bis Tris, 10 KG	Active	10.00 KG		5.00 GAL		White HDPE round pail, blue lid	24 x 36	LDPE	
BTRI-3250-25	Bis Tris, 25 KG	Active	25.00 KG		12.00 GAL		White HDPE round pail, blue lid	30 x 48	LDPE	
BTRI-3250-92	Bis Tris, 100 G	Active	100.00 G		250.00 ML		White HDPE bottle, blue lid			
BTRI-3250-93	Bis Tris, 500 G	Active	500.00 G		1000.00 ML		White HDPE bottle, blue lid			
BTRI-4250-01	Bis-Tris Biotech, 1 KG	Active	1.00 KG		1500.00 ML		White HDPE bottle, blue lid			
BTRI-4250-05	Bis-Tris Biotech, 5 KG	Active	5.00 KG		3.50 GAL		White HDPE round pail, blue lid	24 x 36	LDPE	
BTRI-4250-10	Bis-Tris Biotech, 10 KG	In development	10.00 KG		5.00 GAL		White HDPE round pail, blue lid	24 x 36	LDPE	
BTRI-4250-25	Bis-Tris Biotech, 25 KG	Active	25.00 KG		12.00 GAL		White HDPE round pail, blue lid	30 x 48	LDPE	
BTRI-4250-50	Bis-Tris Biotech, 50 KG	Active	50.00 KG		30.00 GAL		Dark blue HDPE drum	34 x 60	LDPE	
BTRI-4250-92	Bis-Tris Biotech, 100 G	Active	100.00 G		250.00 ML		White HDPE bottle, blue lid			
BTRI-4250-93	Bis-Tris Biotech, 500 G	Active	500.00 G		1000.00 ML		White HDPE bottle, blue lid			
CSMH-3250-01	Cysteamine HCl (2-MEA), 1 KG	In development	1.00 KG		1.25 GAL		White HDPE round pail, blue lid	(2x) 24 x 36	LDPE	4 Unit
CSMH-3250-05	Cysteamine HCl (2-MEA), 5 KG	Active	5.00 KG		3.50 GAL		White HDPE round pail, blue lid	(2x) 24 x 36	LDPE	4 Unit
CSMH-3250-10	Cysteamine HCl (2-MEA), 10 KG	Active	10.00 KG		5.00 GAL		White HDPE round pail, blue lid	(2x) 24 x 36	LDPE	4 Unit
CSMH-3250-25	Cysteamine HCl (2-MEA), 25 KG	Active	25.00 KG		12.00 GAL		White HDPE round pail, blue lid	(2x) 30 x 48	LDPE	
CSMH-3250-90	Cysteamine HCl (2-MEA), 100 G	Active	100.00 G		3.50 GAL		White HDPE round pail, blue lid	(2x) 9 x 15	LDPE	4 Unit
CSMH-3250-91	Cysteamine HCl (2-MEA), 50 G	In development	50.00 G		0.60 GAL		White HDPE round pail, blue lid	(2x) 9 x 15	LDPE	4 Unit
CSMH-3250-92	Cysteamine HCl (2-MEA), 100 G	Active	100.00 G		0.60 GAL		White HDPE round pail, blue lid	(2x) 9 x 15	LDPE	4 Unit
CSMH-3250-94	Cysteamine HCl (2-MEA), 1 KG	In development	1.00 KG		0.60 GAL		White HDPE round pail, blue lid	(2x) 9 x 15	LDPE	4 Unit
CSMH-3250-96	Cysteamine HCl (2-MEA), 1 KG	Active	1.00 KG		3.50 GAL		White HDPE round pail, blue lid	(2x) 24 x 36	LDPE	4 Unit
CSMH-3350-10	Cysteamine HCl (2-MEA), 10 KG	In development	10.00 KG		5.00 GAL		White HDPE round pail, blue lid	(2x) 24 x 36	LDPE	4 Unit
CSMH-3350-25	Cysteamine HCl (2-MEA), 25 KG	In development	25.00 KG		12.00 GAL		White HDPE round pail, blue lid	(2x) 30 x 48	LDPE	4 Unit
D010-3201-01	Dextran Powder MW 10,000, 1 KG	Active	1.00 KG		0.60 GAL		White HDPE round pail, blue lid	(2x) 9 x 15	LDPE	4 Unit
D010-3201-30	Dextran Powder MW 10,000, 30 KG	Active	30.00 KG		30.00 GAL		Dark blue HDPE drum	(2x) 34 x 60	LDPE	4 Unit
DXS4-3250-05	Dextran Sulfate Sodium Salt MW 40,000	In development	5.00 KG		3.50 GAL		White HDPE round pail, blue lid	(2x) 24 x 36	LDPE	4 Unit
DXSE-4250-01	Dextran Sulfate, 8000, GMP, 1kg	Active	1.00 KG		0.60 GAL		White HDPE round pail, blue lid	(2x) 9 x 15	LDPE	4 Unit
DXSE-4250-05	Dextran Sulfate, 8000, GMP, 5kg	Active	5.00 KG		3.50 GAL		White HDPE round pail, blue lid	(2x) 24 x 36	LDPE	4 Unit
DXSE-4250-10	Dextran Sulfate, 8000, GMP, 10kg	Active	10.00 KG		5.00 GAL		White HDPE round pail, blue lid	(2x) 24 x 36	LDPE	4 Unit
DXSE-4250-15	Dextran Sulfate, 8000, GMP, 15kg	Active	15.00 KG		12.00 GAL		White HDPE round pail, blue lid	(2x) 30 x 48	LDPE	4 Unit
DXSE-4250-25	Dextran Sulfate, 8000, GMP, 25kg	Active	25.00 KG		12.00 GAL		White HDPE round pail, blue lid	(2x) 30 x 48	LDPE	4 Unit
DXSE-4250-92	Dextran Sulfate, 8000, GMP, 100g	Active	100.00 G		0.60 GAL		White HDPE round pail, blue lid	(2x) 9 x 15	LDPE	4 Unit
DXSE-4250-93	Dextran Sulfate, 8000, GMP, 20g	Active	20.00 G		0.60 GAL		White HDPE round pail, blue lid	(2x) 9 x 15	LDPE	4 Unit
DXSE-4250-94	Dextran Sulfate, 8000, GMP, 500g	Active	500.00 G		0.60 GAL		White HDPE round pail, blue lid	(2x) 9 x 15	LDPE	4 Unit
DXSE-4252-01	Dextran Sulfate Sodium Salt MW 8000 Biotech, 1 KG	Active	1.00 KG		0.60 GAL		White HDPE round pail, blue lid	(2x) 9 x 15	LDPE	4 Unit
DXSE-4252-05	Dextran Sulfate Sodium Salt MW 8000 Biotech, 5 KG	Active	5.00 KG		3.50 GAL		White HDPE round pail, blue lid	(2x) 24 x 36	LDPE	4 Unit
DXSE-4252-10	Dextran Sulfate Sodium Salt MW 8000 Biotech, 10 KG	Active	10.00 KG		5.00 GAL		White HDPE round pail, blue lid	(2x) 24 x 36	LDPE	4 Unit
DXSE-4252-25	Dextran Sulfate Sodium Salt MW 8000 Biotech, 25 KG	Active	25.00 KG		12.00 GAL		White HDPE round pail, blue lid	(2x) 30 x 48	LDPE	4 Unit
DXSE-4252-50	Dextran Sulfate Sodium Salt MW 8000 Biotech, 50 KG	Active	50.00 KG		30.00 GAL		Dark blue HDPE drum	(2x) 34 x 60	LDPE	4 Unit
DXSE-4252-92	Dextran Sulfate Sodium Salt MW 8000 Biotech, 100 G	Active	100.00 G		0.60 GAL		White HDPE round pail, blue lid	(2x) 9 x 15	LDPE	4 Unit
DXSE-4252-93	Dextran Sulfate Sodium Salt MW 8000 Biotech, 500 G	Active	500.00 G		0.60 GAL		White HDPE round pail, blue lid	(2x) 9 x 15	LDPE	4 Unit
DXSE-5201-92	Dextran Sulfate, 8000, 100g	Active	100.00 G		0.60 GAL		White HDPE round pail, blue lid	(2x) 9 x 15	LDPE	4 Unit
GALP-3224-06	D-Galactose, Plant Derived, 1KG	Active	1.00 KG		0.60 GAL		White HDPE round pail, blue lid	(2x) 9 x 15	LDPE	
GALP-3224-25	D-Galactose, Plant Derived, 20kg	In development	20.00 KG		12.00 GAL		White HDPE round pail, blue lid	30 x 48	LDPE	

Product	Product Description	Product Status	Product size	Product UOM	Package size	Package UOM	Package type	Liner size	Liner type	Desiccant
GALP-3250-05	D-Galactose, Plant Derived, 5kg	In development	5.00	KG	3.50	GAL	White HDPE round pail, blue lid	24 x 36	LDPE	
GALP-3250-06	D-Galactose, Plant Derived, 1KG	Active	1.00	KG	0.60	GAL	White HDPE round pail, blue lid	(2x) 9 x 15	LDPE	
GALP-3250-10	D-Galactose, Plant Derived, 10kg	Active	10.00	KG	5.00	GAL	White HDPE round pail, blue lid	24 x 36	LDPE	
GALP-3250-25	D-Galactose, Plant Derived, 20kg	Active	20.00	KG	12.00	GAL	White HDPE round pail, blue lid	30 x 48	LDPE	
GALP-3250-27	D-Galactose, Plant Derived, 25kg	Active	25.00	KG	30.00	GAL	Dark blue HDPE drum	34 x 60	LDPE	
GALP-3250-92	D-Galactose, Plant Derived, 100 G	Active	100.00	G	250.00	ML	White HDPE bottle, blue lid			
GALP-3251-06	D-Galactose, Plant Derived, 1KG	Active	1.00	KG	0.60	GAL	White HDPE round pail, blue lid	(2x) 9 x 15	LDPE	
GALP-3251-10	D-Galactose, Plant Derived, 10kg	Active	10.00	KG	5.00	GAL	White HDPE round pail, blue lid	24 x 36	LDPE	
GALP-3251-21	D-Galactose, Plant Derived, 20kg	Active	20.00	KG	12.00	GAL	White HDPE round pail, blue lid	30 x 48	LDPE	
GALP-3251-27	D-Galactose, Plant Derived, 25kg	Active	25.00	KG	30.00	GAL	Dark blue HDPE drum	34 x 60	LDPE	
GALP-3251-50	D-Galactose, Plant Derived, 50kg	Active	50.00	KG	30.00	GAL	Dark blue HDPE drum	34 x 60	LDPE	
GALP-3252-06	D-Galactose, Plant Derived, 1KG	Active	1.00	KG	0.60	GAL	White HDPE round pail, blue lid	(2x) 9 x 15	LDPE	
GALP-3252-27	D-Galactose, Plant Derived, 25kg	Active	25.00	KG	30.00	GAL	Dark blue HDPE drum	34 x 60	LDPE	
GALP-3252-92	D-Galactose, Plant Derived, 100 G	Active	100.00	G	250.00	ML	White HDPE bottle, blue lid			
GALP-3350-06	D-Galactose, Plant Derived, 1KG	Active	1.00	KG	0.60	GAL	White HDPE round pail, blue lid	(2x) 9 x 15	LDPE	
GALP-3350-25	D-Galactose, Plant Derived, 20kg	Active	20.00	KG	12.00	GAL	White HDPE round pail, blue lid	30 x 48	LDPE	
GALP-3350-27	D-Galactose, Plant Derived, 25kg	Active	25.00	KG	30.00	GAL	Dark blue HDPE drum	34 x 60	LDPE	
GALP-3351-05	D-Galactose, Plant Derived, 5KG	In development	5.00	KG	3.50	GAL	White HDPE round pail, blue lid	24 x 36	LDPE	
GALP-3351-06	D-Galactose, Plant Derived, 1KG	Active	1.00	KG	0.60	GAL	White HDPE round pail, blue lid	(2x) 9 x 15	LDPE	
GALP-3351-27	D-Galactose, Plant Derived, 25kg	Active	25.00	KG	30.00	GAL	Dark blue HDPE drum	34 x 60	LDPE	
GALP-3351-50	D-Galactose, Plant Derived, 50kg	Active	50.00	KG	30.00	GAL	Dark blue HDPE drum	34 x 60	LDPE	
GALP-3351-92	D-Galactose, Plant Derived, 100 G	In development	100.00	G	250.00	ML	White HDPE bottle, blue lid			
GALP-4250-05	Galactose (D) Biotech, 5 KG	Active	5.00	KG	3.50	GAL	White HDPE round pail, blue lid	24 x 36	LDPE	
GALP-4250-06	Galactose (D) Biotech, 1 KG	Active	1.00	KG	0.60	GAL	White HDPE round pail, blue lid	(2x) 9 x 15	LDPE	
GALP-4250-10	Galactose (D) Biotech, 10 KG	Active	10.00	KG	5.00	GAL	White HDPE round pail, blue lid	24 x 36	LDPE	
GALP-4250-21	Galactose (D) Biotech, 20 KG	Active	20.00	KG	12.00	GAL	White HDPE round pail, blue lid	30 x 48	LDPE	
GALP-4250-27	Galactose (D) Biotech, 25 KG	Active	25.00	KG	30.00	GAL	Dark blue HDPE drum	34 x 60	LDPE	
GALP-4250-50	Galactose (D) Biotech, 50 KG	Active	50.00	KG	30.00	GAL	Dark blue HDPE drum	34 x 60	LDPE	
GALP-4250-92	Galactose (D) Biotech, 100 G	Active	100.00	G	250.00	ML	White HDPE bottle, blue lid			
GHCL-3101-02	Guanidine HCl 6M Solution, 1L	Active	1.00	L	1.00	L	Natural HDPE bottle			
GHCL-3101-72	Guanidine HCl 6M Solution, 200 L	Active	200.00	L	55.00	GAL	Black HDPE drum			
GHCL-3101-78	Guanidine HCl 6M Solution, 1135 L	Active	1135.00	L	330.00	GAL	Tote			
GHCL-3101-80	Guanidine HCl 6M Solution, 200 L	Active	200.00	L	55.00	GAL	Black HDPE drum			
GHCL-3201-59	Guanidine HCl, 50 KG	Active	50.00	KG	26.00	GAL	Fiber drum	(2x) 34 x 60	LDPE	4 Unit
GHCL-3202-01	Guanidine HCl, 1 KG	Active	1.00	KG	1500.00	ML	White HDPE bottle, blue lid			
GHCL-3202-50	Guanidine HCl, 50 KG	Active	50.00	KG	30.00	GAL	Dark blue HDPE drum	34 x 60	LDPE	4 Unit
GHCL-3202-93	Guanidine HCl, 500 G	Active	500.00	G	1000.00	ML	White HDPE bottle, blue lid			
GHCL-3220-01	Guanidine HCl, 1 KG	Active	1.00	KG	1500.00	ML	White HDPE bottle, blue lid			
GHCL-3220-10	Guanidine HCl, 10 KG	Active	10.00	KG	5.00	GAL	White HDPE round pail, blue lid	24 x 36	LDPE	
GHCL-3220-25	Guanidine HCl, 25 KG	Active	25.00	KG	12.00	GAL	White HDPE round pail, blue lid	30 x 48	LDPE	
GHCL-3220-50	Guanidine HCl, 50 KG	In development	50.00	KG	30.00	GAL	Dark blue HDPE drum	34 x 60	LDPE	
GHCL-3220-51	Guanidine HCl, 50 KG	Active	50.00	KG	30.00	GAL	Dark blue HDPE drum	34 x 60	LDPE	4 Unit
GHCL-3220-92	Guanidine HCl, 100 G	Active	100.00	G	250.00	ML	White HDPE bottle, blue lid			
GHCL-3220-93	Guanidine HCl, 500 G	Active	500.00	G	1000.00	ML	White HDPE bottle, blue lid			
GHCL-3220-94	Guanidine HCl, 2.5 KG	Active	2.50	KG	3.50	GAL	White HDPE round pail, blue lid	24 x 36	LDPE	

Product	Product Description	Product Status	Product size	Product UOM	Package size	Package UOM	Package type	Liner size	Liner type	Desiccant
GHCL-3221-25	Guanidine HCl, 25 KG	Active	25.00 KG		12.00 GAL		White HDPE round pail, blue lid	30 x 48	LDPE	
GHCL-3224-25	Guanidine HCl, 25 KG	Active	25.00 KG		12.00 GAL		White HDPE round pail, blue lid	30 x 48	LDPE	
GHCL-3224-93	Guanidine HCl, 500 G	Active	500.00 G		1000.00 ML		White HDPE bottle, blue lid			
GHCL-3225-01	Guanidine HCl, 1KG	Active	1.00 KG		1500.00 ML		White HDPE bottle, blue lid			
GHCL-3225-25	Guanidine HCl, 25KG	Active	25.00 KG		12.00 GAL		White HDPE round pail, blue lid	30 x 48	LDPE	
GHCL-3225-92	Guanidine HCl, 100G	Active	100.00 G		250.00 ML		White HDPE bottle, blue lid			
GHCL-3250-01	Guanidine HCl, 1 KG	Active	1.00 KG		1500.00 ML		White HDPE bottle, blue lid			
GHCL-3251-59	Guanidine HCl, 50 KG	Active	50.00 KG		26.00 GAL		Fiber drum	(2x) 34 x 60	LDPE	4 Unit
GHCL-3251-93	Guanidine HCl, 500 G	Active	500.00 G		1000.00 ML		White HDPE bottle, blue lid			
GHCL-3253-25	Guanidine HCl, 25KG	Active	25.00 KG		12.00 GAL		White HDPE round pail, blue lid	30 x 48	LDPE	
GHCL-3253-93	Guanidine HCl, 500 G	Active	500.00 G		1000.00 ML		White HDPE bottle, blue lid			
GHCL-3255-25	Guanidine HCl, 25KG	In development	25.00 KG		12.00 GAL		White HDPE round pail, blue lid	30 x 48	LDPE	
GHCL-4223-01	Guanidine Hydrochloride Biotech, 1 KG	Active	1.00 KG		1500.00 ML		White HDPE bottle, blue lid			
GHCL-4223-05	Guanidine Hydrochloride Biotech, 5 KG	Active	5.00 KG		3.50 GAL		White HDPE round pail, blue lid	24 x 36	LDPE	
GHCL-4223-10	Guanidine Hydrochloride Biotech, 10 KG	Active	10.00 KG		5.00 GAL		White HDPE round pail, blue lid	24 x 36	LDPE	
GHCL-4223-25	Guanidine Hydrochloride Biotech, 25 KG	Active	25.00 KG		12.00 GAL		White HDPE round pail, blue lid	30 x 48	LDPE	
GHCL-4223-50	Guanidine Hydrochloride Biotech, 50 KG	Active	50.00 KG		30.00 GAL		Dark blue HDPE drum	34 x 60	LDPE	
GHCL-4223-92	Guanidine Hydrochloride Biotech, 100 G	Active	100.00 G		250.00 ML		White HDPE bottle, blue lid			
GHCL-4223-93	Guanidine Hydrochloride Biotech, 500 G	Active	500.00 G		1000.00 ML		White HDPE bottle, blue lid			
GHCL-5220-01	Guanidine HCl, 1 KG	Active	1.00 KG		1500.00 ML		White HDPE bottle, blue lid			
GHCL-5220-25	Guanidine HCl, 25 KG	Active	25.00 KG		12.00 GAL		White HDPE round pail, blue lid	30 x 48	LDPE	
GTHI-3220-10	Guanidine Thiocyanate, 10 KG	Active	10.00 KG		5.00 GAL		White HDPE round pail, blue lid	24 x 36	LDPE	
GTHI-3220-14-76	Guanidine Thiocyanate, 14.76 KG	Active	14.76 KG		30.00 GAL		Dark blue HDPE drum	34 x 60	LDPE/ 25L EZI-Flow	
GTHI-3220-15	Guanidine Thiocyanate, 14.76 KG	In development	14.76 KG		12.00 GAL		White HDPE round pail, blue lid	30 x 48	LDPE	
GTHI-3220-25	Guanidine Thiocyanate, 25 KG	Active	25.00 KG		12.00 GAL		White HDPE round pail, blue lid	30 x 48	LDPE	
GTHI-3220-92	Guanidine Thiocyanate, 100 G	Active	100.00 G		250.00 ML					
GTHI-3220-95	Guanidine Thiocyanate, 10 G	In development	10.00 G		60.00 ML					
GTHI-4220-05	Guanidine Thiocyanate, 5 KG	Active	5.00 KG		3.50 GAL		White HDPE round pail, blue lid	24 x 36	LDPE	
GTHI-4220-10	Guanidine Thiocyanate, 10 KG	Active	10.00 KG		5.00 GAL		White HDPE round pail, blue lid	24 x 36	LDPE	
GTHI-4220-25	Guanidine Thiocyanate, 25 KG	Active	25.00 KG		12.00 GAL		White HDPE round pail, blue lid	30 x 48	LDPE	
GTHI-4221-01	Guanidine Thiocyanate Biotech, 1 KG	Active	1.00 KG		1500.00 ML		White HDPE bottle, blue lid			
GTHI-4221-05	Guanidine Thiocyanate Biotech, 5 KG	Active	5.00 KG		3.50 GAL		White HDPE round pail, blue lid	24 x 36	LDPE	
GTHI-4221-10	Guanidine Thiocyanate Biotech, 10 KG	Active	10.00 KG		5.00 GAL		White HDPE round pail, blue lid	24 x 36	LDPE	
GTHI-4221-25	Guanidine Thiocyanate Biotech, 25 KG	Active	25.00 KG		12.00 GAL		White HDPE round pail, blue lid	30 x 48	LDPE	
GTHI-4221-50	Guanidine Thiocyanate Biotech, 50 KG	Active	50.00 KG		30.00 GAL		Dark blue HDPE drum	34 x 60	LDPE	
GTHI-4221-92	Guanidine Thiocyanate Biotech, 100 G	Active	100.00 G		250.00 ML		White HDPE bottle, blue lid			
GTHI-4221-93	Guanidine Thiocyanate Biotech, 500 G	Active	500.00 G		1000.00 ML		White HDPE bottle, blue lid			
GTHI-5203-10	Guanidine Thiocyanate, 10 KG	Active	10.00 KG		5.00 GAL		White HDPE round pail, blue lid	24 x 36	LDPE	
GTHI-5220-10	Guanidine Thiocyanate, 10 KG	Active	10.00 KG		5.00 GAL		White HDPE round pail, blue lid	24 x 36	LDPE	
GTHI-5221-05	Guanidine Thiocyanate, 5 KG	Active	5.00 KG		3.50 GAL		White HDPE round pail, blue lid	24 x 36	LDPE	
GTHI-5221-10	Guanidine Thiocyanate, 10 KG	Active	10.00 KG		5.00 GAL		White HDPE round pail, blue lid	24 x 36	LDPE	
GTHI-5221-25	Guanidine Thiocyanate, 25 KG	Active	25.00 KG		12.00 GAL		White HDPE round pail, blue lid	30 x 48	LDPE	
GTHI-5221-92	Guanidine Thiocyanate, 100 G	Active	100.00 G		250.00 ML					
HEHH-5201-25	HEPES Hemi Sodium, 25 KG	Active	25.00 KG		12.00 GAL		White HDPE round pail, blue lid	30 x 48	LDPE	
HEPE-3220-01	HEPES, 1 KG	Active	1.00 KG		1500.00 ML		White HDPE bottle, blue lid			

Product	Product Description	Product Status	Product size	Product UOM	Package size	Package UOM	Package type	Liner size	Liner type	Desiccant
HEPE-3220-06	HEPES, 1 KG	Active	1.00 KG		0.60 GAL		White HDPE round pail, blue lid	(2x) 9 x 15	LDPE	
HEPE-3220-10	HEPES, 10 KG	Active	10.00 KG		5.00 GAL		White HDPE round pail, blue lid	24 x 36	LDPE	
HEPE-3220-12	HEPES, 12 KG	Active	12.00 KG		5.00 GAL		White HDPE round pail, blue lid	24 x 36	LDPE	
HEPE-3220-20	HEPES, 20 KG	Active	20.00 KG		12.00 GAL		White HDPE round pail, blue lid	(2x) 30 x 48	LDPE	
HEPE-3220-22	HEPES, 2 KG	Active	2.00 KG		1.25 GAL		White HDPE round pail, blue lid	24 x 36	LDPE	
HEPE-3220-25	HEPES, 25 KG	Active	25.00 KG		12.00 GAL		White HDPE round pail, blue lid	30 x 48	LDPE	
HEPE-3220-50	HEPES, 50 KG	Active	50.00 KG		30.00 GAL		Dark blue HDPE drum	34 x 60	LDPE	
HEPE-3220-91	HEPES, 250 G	Active	250.00 G		1000.00 ML		White HDPE bottle, blue lid			
HEPE-3220-92	HEPES, 100 G	Active	100.00 G		250.00 ML		White HDPE bottle, blue lid			
HEPE-3220-93	HEPES, 500 G	Active	500.00 G		1000.00 ML		White HDPE bottle, blue lid			
HEPE-3221-10	HEPES, 10 KG	Active	10.00 KG		5.00 GAL		White HDPE round pail, blue lid	24 x 36	LDPE	
HEPE-3250-05	HEPES, 5 KG	Active	5.00 KG		3.50 GAL		White HDPE round pail, blue lid	24 x 36	LDPE	
HEPE-3250-10	HEPES, 10 KG	Active	10.00 KG		5.00 GAL		White HDPE round pail, blue lid	24 x 36	LDPE	
HEPE-3251-01	HEPES, 1kg	Active	1.00 KG		1500.00 ML		White HDPE bottle, blue lid			
HEPE-3251-25	HEPES, 25 KG	Active	25.00 KG		12.00 GAL		White HDPE round pail, blue lid	30 x 48	LDPE	
HEPE-3251-50	HEPES, 50 KG	Active	50.00 KG		30.00 GAL		Dark blue HDPE drum	34 x 60	LDPE	
HEPE-3251-92	HEPES, 100 G	Active	100.00 G		250.00 ML		White HDPE bottle, blue lid			
HEPE-3254-25	HEPES, 25 KG	In development	25.00 KG		12.00 GAL		White HDPE round pail, blue lid	30 x 48	LDPE	
HEPE-3254-50	HEPES, 50 KG	In development	50.00 KG		30.00 GAL		Dark blue HDPE drum	34 x 60	LDPE	
HEPE-3320-25	HEPES, 25 KG	In development	25.00 KG		12.00 GAL		White HDPE round pail, blue lid	30 x 48	LDPE	
HEPE-3320-92	HEPES, 100 G	In development	100.00 G		250.00 ML		White HDPE bottle, blue lid			
HEPE-3351-25	HEPES, 25 KG	In development	25.00 KG		12.00 GAL		White HDPE round pail, blue lid	30 x 48	LDPE	
HEPE-3351-92	HEPES, 100 G	In development	100.00 G		250.00 ML		White HDPE bottle, blue lid			
HEPE-3354-25	HEPES, 25 KG	In development	25.00 KG		12.00 GAL		White HDPE round pail, blue lid	30 x 48	LDPE	
HEPE-3354-50	HEPES, 50 KG	Active	50.00 KG		30.00 GAL		Dark blue HDPE drum	34 x 60	LDPE	
HEPE-4250-01	HEPES Biotech, 1 KG	Active	1.00 KG		1500.00 ML		White HDPE bottle, blue lid			
HEPE-4250-05	HEPES Biotech, 5 KG	Active	5.00 KG		3.50 GAL		White HDPE round pail, blue lid	24 x 36	LDPE	
HEPE-4250-10	HEPES Biotech, 10 KG	Active	10.00 KG		5.00 GAL		White HDPE round pail, blue lid	24 x 36	LDPE	
HEPE-4250-25	HEPES Biotech, 25 KG	Active	25.00 KG		12.00 GAL		White HDPE round pail, blue lid	30 x 48	LDPE	
HEPE-4250-50	HEPES Biotech, 50 KG	Active	50.00 KG		30.00 GAL		Dark blue HDPE drum	34 x 60	LDPE	
HEPE-4250-92	HEPES Biotech, 100 G	Active	100.00 G		250.00 ML		White HDPE bottle, blue lid			
HEPE-4250-93	HEPES Biotech, 500 G	Active	500.00 G		1000.00 ML		White HDPE bottle, blue lid			
HEPE-5201-26	HEPES, 25 KG	Active	25.00 KG		26.00 GAL		Fiber drum	(2x) 34 x 60	LDPE	
HEPE-5203-40	HEPES, 40 KG	Active	40.00 KG		30.00 GAL		Dark blue HDPE drum	34 x 60	LDPE	
HEPE-5220-01	HEPES, 1 KG	Active	1.00 KG		1500.00 ML		White HDPE bottle, blue lid			
HEPE-5220-10	HEPES, 10 KG	Active	10.00 KG		5.00 GAL		White HDPE round pail, blue lid	24 x 36	LDPE	
HEPN-5201-26	HEPES Sodium, 25 KG	Active	25.00 KG		26.00 GAL		Fiber drum	(2x) 34 x 60	LDPE	
HEPN-5203-50	HEPES Sodium, 50 KG	Active	50.00 KG		30.00 GAL		Dark blue HDPE drum	34 x 60	LDPE	
HEPN-5220-25	HEPES Sodium, 25 KG	Active	25.00 KG		12.00 GAL		White HDPE round pail, blue lid	30 x 48	LDPE	
HEPN-5251-25	HEPES Sodium, 25 KG	Active	25.00 KG		12.00 GAL		White HDPE round pail, blue lid	30 x 48	LDPE	
IHCL-4101-02	6N HCl in IPA Solution, 500ML	Active	500.00 ML		17.00 OZ		Black HDPE UN leakproof bottle			
IHCL-4101-97	6N HCl in IPA Solution, 180 KG	Active	180.00 KG		55.00 GAL		Black HDPE drum			
KBRO-2220-01	Potassium Bromide, 1 KG	Active	1.00 KG		1500.00 ML		White HDPE bottle, blue lid			
KBRO-2220-05	Potassium Bromide, 5 KG	Active	5.00 KG		3.50 GAL		White HDPE round pail, blue lid	(2x) 24 x 36	LDPE	4 Unit
KBRO-2220-10	Potassium Bromide, 10 KG	Active	10.00 KG		5.00 GAL		White HDPE round pail, blue lid	(2x) 24 x 36	LDPE	4 Unit

Product	Product Description	Product Status	Product size	Product UOM	Package size	Package UOM	Package type	Liner size	Liner type	Desiccant
KBRO-2220-25	Potassium Bromide, 25 KG	Active	25.00 KG		12.00 GAL		White HDPE round pail, blue lid	(2x) 30 x 48	LDPE	4 Unit
KBRO-2220-92	Potassium Bromide, 100 G	Active	100.00 G		0.60 GAL		White HDPE round pail, blue lid	9 x 15	LDPE	
KBRO-2301-01	Potassium Bromide, 1 KG	Active	1.00 KG		1500.00 ML		White HDPE bottle, blue lid			
KBRO-2301-25	Potassium Bromide, 25 KG	Active	25.00 KG		12.00 GAL		White HDPE round pail, blue lid	(2x) 30 x 48	LDPE	
KBRO-2302-01	Potassium Bromide, 1 KG	Active	1.00 KG		1500.00 ML		White HDPE bottle, blue lid			
KBRO-2302-25	Potassium Bromide, 25 KG	Active	25.00 KG		12.00 GAL		White HDPE round pail, blue lid	(2x) 30 x 48	LDPE	
LARH-5350-92	Arginine (L) Hydrochloride Biotech, 100 G	Active	100.00 G		250.00 ML		White HDPE bottle, blue lid			
LARH-6250-01	Arginine (L) Hydrochloride Biotech, 1 KG	In development	1.00 KG		1500.00 ML		White HDPE bottle, blue lid			
LARH-6250-05	Arginine (L) Hydrochloride Biotech, 5 KG	In development	5.00 KG		3.50 GAL		White HDPE round pail, blue lid	24 x 36	LDPE	
LARH-6250-10	Arginine (L) Hydrochloride Biotech, 10 KG	In development	10.00 KG		5.00 GAL		White HDPE round pail, blue lid	24 x 36	LDPE	
LARH-6250-25	Arginine (L) Hydrochloride Biotech, 25 KG	In development	25.00 KG		12.00 GAL		White HDPE round pail, blue lid	30 x 48	LDPE	
LARH-6250-50	Arginine (L) Hydrochloride Biotech, 50 KG	In development	50.00 KG		30.00 GAL		Dark blue HDPE drum	34 x 60	LDPE	
LARH-6250-92	Arginine (L) Hydrochloride Biotech, 100 G	In development	100.00 G		250.00 ML		White HDPE bottle, blue lid			
LARH-6250-93	Arginine (L) Hydrochloride Biotech, 500 G	In development	500.00 G		1000.00 ML		White HDPE bottle, blue lid			
LARH-6350-01	Arginine (L) Hydrochloride Biotech, 1 KG	Active	1.00 KG		1500.00 ML		White HDPE bottle, blue lid			
LARH-6350-05	Arginine (L) Hydrochloride Biotech, 5 KG	Active	5.00 KG		3.50 GAL		White HDPE round pail, blue lid	24 x 36	LDPE	
LARH-6350-10	Arginine (L) Hydrochloride Biotech, 10 KG	Active	10.00 KG		5.00 GAL		White HDPE round pail, blue lid	24 x 36	LDPE	
LARH-6350-25	Arginine (L) Hydrochloride Biotech, 25 KG	Active	25.00 KG		12.00 GAL		White HDPE round pail, blue lid	30 x 48	LDPE	
LARH-6350-50	Arginine (L) Hydrochloride Biotech, 50 KG	Active	50.00 KG		30.00 GAL		Dark blue HDPE drum	34 x 60	LDPE	
LARH-6350-92	Arginine (L) Hydrochloride Biotech, 100 G	Active	100.00 G		250.00 ML		White HDPE bottle, blue lid			
LARH-6350-93	Arginine (L) Hydrochloride Biotech, 500 G	Active	500.00 G		1000.00 ML		White HDPE bottle, blue lid			
LCYS-6250-01	Cystine (L) Hydrochloride Biotech, 1 KG	Active	1.00 KG		1500.00 ML		White HDPE bottle, blue lid			
LCYS-6250-05	Cystine (L) Hydrochloride Biotech, 5 KG	Active	5.00 KG		3.50 GAL		White HDPE round pail, blue lid	24 x 36	LDPE	
LCYS-6250-10	Cystine (L) Hydrochloride Biotech, 10 KG	Active	10.00 KG		5.00 GAL		White HDPE round pail, blue lid	24 x 36	LDPE	
LCYS-6250-25	Cystine (L) Hydrochloride Biotech, 25 KG	Active	25.00 KG		12.00 GAL		White HDPE round pail, blue lid	30 x 48	LDPE	
LCYS-6250-50	Cystine (L) Hydrochloride Biotech, 50 KG	Active	50.00 KG		30.00 GAL		Dark blue HDPE drum	34 x 60	LDPE	
LCYS-6250-92	Cystine (L) Hydrochloride Biotech, 100 G	Active	100.00 G		250.00 ML		White HDPE bottle, blue lid			
LCYS-6250-93	Cystine (L) Hydrochloride Biotech, 500 G	Active	500.00 G		1000.00 ML		White HDPE bottle, blue lid			
LCYS-6350-01	Cystine (L) Hydrochloride Biotech, 1 KG	Active	1.00 KG		1500.00 ML		White HDPE bottle, blue lid			
LCYS-6350-05	Cystine (L) Hydrochloride Biotech, 5 KG	Active	5.00 KG		3.50 GAL		White HDPE round pail, blue lid	24 x 36	LDPE	
LCYS-6350-10	Cystine (L) Hydrochloride Biotech, 10 KG	Active	10.00 KG		5.00 GAL		White HDPE round pail, blue lid	24 x 36	LDPE	
LCYS-6350-25	Cystine (L) Hydrochloride Biotech, 25 KG	Active	25.00 KG		12.00 GAL		White HDPE round pail, blue lid	30 x 48	LDPE	
LCYS-6350-50	Cystine (L) Hydrochloride Biotech, 50 KG	Active	50.00 KG		30.00 GAL		Dark blue HDPE drum	34 x 60	LDPE	
LCYS-6350-92	Cystine (L) Hydrochloride Biotech, 100 G	Active	100.00 G		250.00 ML		White HDPE bottle, blue lid			
LCYS-6350-93	Cystine (L) Hydrochloride Biotech, 500 G	Active	500.00 G		1000.00 ML		White HDPE bottle, blue lid			
LGLM-6201-01	Glutamine (L) Biotech, 1 KG	Active	1.00 KG		1500.00 ML		White HDPE bottle, blue lid			
LGLM-6201-05	Glutamine (L) Biotech, 5 KG	Active	5.00 KG		3.50 GAL		White HDPE round pail, blue lid	24 x 36	LDPE	
LGLM-6201-10	Glutamine (L) Biotech, 10 KG	Active	10.00 KG		5.00 GAL		White HDPE round pail, blue lid	24 x 36	LDPE	
LGLM-6201-25	Glutamine (L) Biotech, 25 KG	Active	25.00 KG		12.00 GAL		White HDPE round pail, blue lid	30 x 48	LDPE	
LGLM-6201-50	Glutamine (L) Biotech, 50 KG	Active	50.00 KG		30.00 GAL		Dark blue HDPE drum	34 x 60	LDPE	
LGLM-6201-92	Glutamine (L) Biotech, 100 G	Active	100.00 G		250.00 ML		White HDPE bottle, blue lid			
LGLM-6201-93	Glutamine (L) Biotech, 500 G	Active	500.00 G		1000.00 ML		White HDPE bottle, blue lid			
LGLM-6301-01	Glutamine (L) Biotech, 1 KG	In development	1.00 KG		1500.00 ML		White HDPE bottle, blue lid			
LGLM-6301-05	Glutamine (L) Biotech, 5 KG	In development	5.00 KG		3.50 GAL		White HDPE round pail, blue lid	24 x 36	LDPE	
LGLM-6301-10	Glutamine (L) Biotech, 10 KG	In development	10.00 KG		5.00 GAL		White HDPE round pail, blue lid	24 x 36	LDPE	

Product	Product Description	Product Status	Product size	Product UOM	Package size	Package UOM	Package type	Liner size	Liner type	Desiccant
LGLM-6301-25	Glutamine (L) Biotech, 25 KG	In development	25.00	KG	12.00	GAL	White HDPE round pail, blue lid	30 x 48	LDPE	
LGLM-6301-50	Glutamine (L) Biotech, 50 KG	In development	50.00	KG	30.00	GAL	Dark blue HDPE drum	34 x 60	LDPE	
LGLM-6301-92	Glutamine (L) Biotech, 100 G	In development	100.00	G	250.00	ML	White HDPE bottle, blue lid			
LGLM-6301-93	Glutamine (L) Biotech, 500 G	In development	500.00	G	1000.00	ML	White HDPE bottle, blue lid			
LHMM-6250-01	Histidine (L) MonoHCl Monohydrate Biotech, 1 KG	Active	1.00	KG	1500.00	ML	White HDPE bottle, blue lid			
LHMM-6250-05	Histidine (L) MonoHCl Monohydrate Biotech, 5 KG	Active	5.00	KG	3.50	GAL	White HDPE round pail, blue lid	24 x 36	LDPE	
LHMM-6250-10	Histidine (L) MonoHCl Monohydrate Biotech, 10 KG	Active	10.00	KG	5.00	GAL	White HDPE round pail, blue lid	24 x 36	LDPE	
LHMM-6250-25	Histidine (L) MonoHCl Monohydrate Biotech, 25 KG	Active	25.00	KG	12.00	GAL	White HDPE round pail, blue lid	30 x 48	LDPE	
LHMM-6250-50	Histidine (L) MonoHCl Monohydrate Biotech, 50 KG	Active	50.00	KG	30.00	GAL	Dark blue HDPE drum	34 x 60	LDPE	
LHMM-6250-92	Histidine (L) MonoHCl Monohydrate Biotech, 100 G	Active	100.00	G	250.00	ML	White HDPE bottle, blue lid			
LHMM-6250-93	Histidine (L) MonoHCl Monohydrate Biotech, 500 G	Active	500.00	G	1000.00	ML	White HDPE bottle, blue lid			
LHMM-6350-01	Histidine (L) MonoHCl Monohydrate Biotech, 1 KG	Active	1.00	KG	1500.00	ML	White HDPE bottle, blue lid			
LHMM-6350-05	Histidine (L) MonoHCl Monohydrate Biotech, 5 KG	Active	5.00	KG	3.50	GAL	White HDPE round pail, blue lid	24 x 36	LDPE	
LHMM-6350-10	Histidine (L) MonoHCl Monohydrate Biotech, 10 KG	Active	10.00	KG	5.00	GAL	White HDPE round pail, blue lid	24 x 36	LDPE	
LHMM-6350-25	Histidine (L) MonoHCl Monohydrate Biotech, 25 KG	Active	25.00	KG	12.00	GAL	White HDPE round pail, blue lid	30 x 48	LDPE	
LHMM-6350-50	Histidine (L) MonoHCl Monohydrate Biotech, 50 KG	Active	50.00	KG	30.00	GAL	Dark blue HDPE drum	34 x 60	LDPE	
LHMM-6350-92	Histidine (L) MonoHCl Monohydrate Biotech, 100 G	Active	100.00	G	250.00	ML	White HDPE bottle, blue lid			
LHMM-6350-93	Histidine (L) MonoHCl Monohydrate Biotech, 500 G	Active	500.00	G	1000.00	ML	White HDPE bottle, blue lid			
MESH-3201-50	MES Hydrate, 50KG	In development	50.00	KG	30.00	GAL	Dark blue HDPE drum	34 x 60	LDPE	
MESH-3201-92	MES Hydrate, 100 G	Active	100.00	G	250.00	ML	White HDPE bottle, blue lid			
MESH-3250-01	MES Hydrate, 1KG	Active	1.00	KG	1500.00	ML	White HDPE bottle, blue lid			
MESH-3250-25	MES Hydrate, 25KG	Active	25.00	KG	12.00	GAL	White HDPE round pail, blue lid	30 x 48	LDPE	
MESH-3350-01	MES Hydrate, 1KG	In development	1.00	KG	1500.00	ML	White HDPE bottle, blue lid			
MESH-4250-01	MES Hydrate Biotech, 1 KG	Active	1.00	KG	1500.00	ML	White HDPE bottle, blue lid			
MESH-4250-05	MES Hydrate Biotech, 5 KG	Active	5.00	KG	3.50	GAL	White HDPE round pail, blue lid	24 x 36	LDPE	
MESH-4250-10	MES Hydrate Biotech, 10 KG	Active	10.00	KG	5.00	GAL	White HDPE round pail, blue lid	24 x 36	LDPE	
MESH-4250-25	MES Hydrate Biotech, 25 KG	Active	25.00	KG	12.00	GAL	White HDPE round pail, blue lid	30 x 48	LDPE	
MESH-4250-50	MES Hydrate Biotech, 50 KG	Active	50.00	KG	30.00	GAL	Dark blue HDPE drum	34 x 60	LDPE	
MESH-4250-92	MES Hydrate Biotech, 100 G	Active	100.00	G	250.00	ML	White HDPE bottle, blue lid			
MESH-4250-93	MES Hydrate Biotech, 500 G	Active	500.00	G	1000.00	ML	White HDPE bottle, blue lid			
MESH-9250-91	MES Hydrate, 50 G	In development	50.00	G	250.00	ML	White HDPE bottle, blue lid			
MESM-3220-01	MES Monohydrate, 1 KG	Active	1.00	KG	1500.00	ML	White HDPE bottle, blue lid			
MESM-3220-05	MES Monohydrate, 5 KG	Active	5.00	KG	3.50	GAL	White HDPE round pail, blue lid	24 x 36	LDPE	
MESM-3220-10	MES Monohydrate, 10 KG	Active	10.00	KG	5.00	GAL	White HDPE round pail, blue lid	24 x 36	LDPE	
MESM-3220-12	MES, Monohydrate, 12 KG	Active	12.00	KG	5.00	GAL	White HDPE round pail, blue lid	24 x 36	LDPE	
MESM-3220-25	MES, Monohydrate, 25 KG	Active	25.00	KG	12.00	GAL	White HDPE round pail, blue lid	30 x 48	LDPE	
MESM-3220-50	MES, Monohydrate, 50 KG	Active	50.00	KG	30.00	GAL	Dark blue HDPE drum	34 x 60	LDPE	
MESM-3220-91	MES Monohydrate, 50G	Active	50.00	G	250.00	ML	White HDPE bottle, blue lid			
MESM-3220-92	MES, Monohydrate, 100G	Active	100.00	G	250.00	ML	White HDPE bottle, blue lid			
MESM-3220-93	MES, Monohydrate, 500G	Active	500.00	G	1000.00	ML	White HDPE bottle, blue lid			
MESM-3221-10	MES Monohydrate, 10 KG	Active	10.00	KG	5.00	GAL	White HDPE round pail, blue lid	24 x 36	LDPE	
MESM-3221-93	MES, Monohydrate, 500G	Active	500.00	G	1000.00	ML	White HDPE bottle, blue lid			
MESM-3222-01	MES Monohydrate, 1 KG	Active	1.00	KG	1500.00	ML	White HDPE bottle, blue lid			
MESM-3222-03	MES Monohydrate, 3 KG	Active	3.00	KG	3.50	GAL	White HDPE round pail, blue lid	24 x 36		
MESM-3222-25	MES, Monohydrate, 25 KG	Active	25.00	KG	12.00	GAL	White HDPE round pail, blue lid	30 x 48	LDPE	

Product	Product Description	Product Status	Product size	Product UOM	Package size	Package UOM	Package type	Liner size	Liner type	Desiccant
MESM-3222-50	MES, Monohydrate, 50 KG	Active	50.00	KG	30.00	GAL	Dark blue HDPE drum	34 x 60	LDPE	
MESM-3222-91	MES Monohydrate, 50G	Active	50.00	G	250.00	ML	White HDPE bottle, blue lid			
MESM-3222-92	MES, Monohydrate, 100G	Active	100.00	G	250.00	ML	White HDPE bottle, blue lid			
MESM-3222-93	MES, Monohydrate, 500G	In development	500.00	G	1000.00	ML	White HDPE bottle, blue lid			
MESM-3223-05	MES, Monohydrate, 5 KG	Active	5.00	KG	3.50	GAL	White HDPE round pail, blue lid	24 x 36	LDPE	
MESM-3250-06	MES, Monohydrate, 1 KG	Active	1.00	KG	0.60	GAL	White HDPE round pail, blue lid	(2x) 9 x 15	LDPE	
MESM-3250-10	MES Monohydrate, 10 KG	Active	10.00	KG	5.00	GAL	White HDPE round pail, blue lid	24 x 36	LDPE	
MESM-3250-12	MES, Monohydrate, 12 KG	Active	12.00	KG	5.00	GAL	White HDPE round pail, blue lid	24 x 36	LDPE	
MESM-3250-25	MES, Monohydrate, 25 KG	Active	25.00	KG	12.00	GAL	White HDPE round pail, blue lid	30 x 48	LDPE	
MESM-3250-50	MES, Monohydrate, 50 KG	Active	50.00	KG	30.00	GAL	Dark blue HDPE drum	34 x 60	LDPE	
MESM-3250-91	MES, Monohydrate, 250G	Active	250.00	G	1000.00	ML	White HDPE bottle, blue lid			
MESM-3250-92	MES, Monohydrate, 100G	Active	100.00	G	250.00	ML	White HDPE bottle, blue lid			
MESM-3250-93	MES, Monohydrate, 500G	Active	500.00	G	1000.00	ML	White HDPE bottle, blue lid			
MESM-3250-95	MES, Monohydrate, 10G	In development	10.00	G	250.00	ML	White HDPE bottle, blue lid			
MESM-3251-25	MES, Monohydrate, 25 KG	Active	25.00	KG	12.00	GAL	White HDPE round pail, blue lid	30 x 48	LDPE	
MESM-3251-92	MES, Monohydrate, 100G	Active	100.00	G	250.00	ML	White HDPE bottle, blue lid			
MESM-3320-25	MES, Monohydrate, 25 KG	In development	25.00	KG	12.00	GAL	White HDPE round pail, blue lid	30 x 48	LDPE	
MESM-3320-92	MES, Monohydrate, 100G	In development	100.00	G	250.00	ML	White HDPE bottle, blue lid			
MESM-3321-10	MES Monohydrate, 10 KG	In development	10.00	KG	5.00	GAL	White HDPE round pail, blue lid	24 x 36	LDPE	
MESM-3321-93	MES, Monohydrate, 500G	In development	500.00	G	1000.00	ML	White HDPE bottle, blue lid			
MESM-3350-25	MES, Monohydrate, 25 KG	In development	25.00	KG	12.00	GAL	White HDPE round pail, blue lid	30 x 48	LDPE	
MESM-3350-92	MES, Monohydrate, 100G	In development	100.00	G	250.00	ML	White HDPE bottle, blue lid			
MESM-4220-01	MES, Monohydrate, 1 KG	Active	1.00	KG	1500.00	ML	White HDPE bottle, blue lid			
MESM-4220-10	MES, Monohydrate, 10 KG	Active	10.00	KG	5.00	GAL	White HDPE round pail, blue lid	24 x 36	LDPE	
MESM-4220-25	MES, Monohydrate, 25 KG	Active	25.00	KG	12.00	GAL	White HDPE round pail, blue lid	30 x 48	LDPE	
MESM-4220-26	MES, Monohydrate, 25 KG	Active	25.00	KG	26.00	GAL	Fiber drum	(2x) 34 x 60	LDPE	
MESM-4221-01	MES, Monohydrate Biotech, 1 KG	Active	1.00	KG	1500.00	ML	White HDPE bottle, blue lid			
MESM-4221-05	MES, Monohydrate Biotech, 5 KG	Active	5.00	KG	3.50	GAL	White HDPE round pail, blue lid	24 x 36	LDPE	
MESM-4221-10	MES, Monohydrate Biotech, 10 KG	Active	10.00	KG	5.00	GAL	White HDPE round pail, blue lid	24 x 36	LDPE	
MESM-4221-25	MES, Monohydrate Biotech, 25 KG	Active	25.00	KG	12.00	GAL	White HDPE round pail, blue lid	30 x 48	LDPE	
MESM-4221-50	MES, Monohydrate Biotech, 50 KG	Active	50.00	KG	30.00	GAL	Dark blue HDPE drum	34 x 60	LDPE	
MESM-4221-92	MES, Monohydrate Biotech, 100 G	Active	100.00	G	250.00	ML	White HDPE bottle, blue lid			
MESM-4221-93	MES, Monohydrate Biotech, 500 G	Active	500.00	G	1000.00	ML	White HDPE bottle, blue lid			
MESM-5201-26	MES, Monohydrate, 25 KG	Active	25.00	KG	26.00	GAL	Fiber drum	(2x) 34 x 60	LDPE	
MESM-5202-10	MES, Monohydrate, 10 KG	Active	10.00	KG	5.00	GAL	White HDPE round pail, blue lid	24 x 36	LDPE	
MESM-5220-10	MES, Monohydrate, 10 KG	Active	10.00	KG	5.00	GAL	White HDPE round pail, blue lid	24 x 36	LDPE	
MESM-5220-92	MES, Monohydrate, 100 G	Active	100.00	G	250.00	ML	White HDPE bottle, blue lid			
MESM-5220-93	MES, Monohydrate, 500 G	Active	500.00	G	1000.00	ML	White HDPE bottle, blue lid			
MESN-5220-10	MES Sodium, 10 KG	Active	10.00	KG	5.00	GAL	White HDPE round pail, blue lid	24 x 36	LDPE	
MESN-5220-26	MES Sodium, 25 KG	Active	25.00	KG	26.00	GAL	Fiber drum	(2x) 34 x 60	LDPE	
MESN-5220-50	MES Sodium, 50 KG	Active	50.00	KG	30.00	GAL	Dark blue HDPE drum	34 x 60	LDPE	
MESN-5220-93	MES Sodium, 300 G	Active	300.00	G	1000.00	ML	White HDPE bottle, blue lid			
MOPS-3201-10	MOPS, 10 KG	Active	10.00	KG	5.00	GAL	White HDPE round pail, blue lid	24 x 36	LDPE	
MOPS-3201-51	MOPS, 50 KG	Active	50.00	KG	30.00	GAL	Dark blue HDPE drum	34 x 60	LDPE	
MOPS-3201-93	MOPS, 500 G	Active	500.00	G	1000.00	ML	White HDPE bottle, blue lid			

Product	Product Description	Product Status	Product size	Product UOM	Package size	Package UOM	Package type	Liner size	Liner type	Desiccant
MOPS-3220-01	MOPS, 1 KG	Active	1.00 KG		1500.00 ML		White HDPE bottle, blue lid			
MOPS-3220-25	MOPS, 25 KG	Active	25.00 KG		12.00 GAL		White HDPE round pail, blue lid	30 x 48	LDPE	
MOPS-3220-50	MOPS, 50 KG	Active	50.00 KG		30.00 GAL		Dark blue HDPE drum	34 x 60	LDPE	
MOPS-3220-51	MOPS, 50 KG	Active	50.00 KG		30.00 GAL		Dark blue HDPE drum	34 x 60	LDPE	
MOPS-3220-92	MOPS, 100 G	Active	100.00 G		250.00 ML		White HDPE bottle, blue lid			
MOPS-3220-93	MOPS, 500 G	Active	500.00 G		1000.00 ML		White HDPE bottle, blue lid			
MOPS-4220-05	MOPS, 5 KG	Active	5.00 KG		3.50 GAL		White HDPE round pail, blue lid	24 x 36	LDPE	
MOPS-4221-01	MOPS Biotech, 1 KG	Active	1.00 KG		1500.00 ML		White HDPE bottle, blue lid			
MOPS-4221-05	MOPS Biotech, 5 KG	Active	5.00 KG		3.50 GAL		White HDPE round pail, blue lid	24 x 36	LDPE	
MOPS-4221-10	MOPS Biotech, 10 KG	Active	10.00 KG		5.00 GAL		White HDPE round pail, blue lid	24 x 36	LDPE	
MOPS-4221-25	MOPS Biotech, 25 KG	Active	25.00 KG		12.00 GAL		White HDPE round pail, blue lid	30 x 48	LDPE	
MOPS-4221-50	MOPS Biotech, 50 KG	Active	50.00 KG		30.00 GAL		Dark blue HDPE drum	34 x 60	LDPE	
MOPS-4221-92	MOPS Biotech, 100 G	Active	100.00 G		250.00 ML		White HDPE bottle, blue lid			
MOPS-4221-93	MOPS Biotech, 500 G	Active	500.00 G		1000.00 ML		White HDPE bottle, blue lid			
MOPS-5201-26	MOPS, 50 KG	Active	50.00 KG		26.00 GAL		Fiber drum	(2x) 34 x 60	LDPE	
MOPS-5202-25	MOPS, 25 KG	Active	25.00 KG		12.00 GAL		White HDPE round pail, blue lid	30 x 48	LDPE	
MOPS-5203-10	MOPS, 10 KG	Active	10.00 KG		5.00 GAL		White HDPE round pail, blue lid	24 x 36	LDPE	
MOPS-5220-51	MOPS, 50 KG	Active	50.00 KG		30.00 GAL		Dark blue HDPE drum	34 x 60	LDPE	
MOPS-5220-92	MOPS, 100 G	Active	100.00 G		250.00 ML		White HDPE round pail, blue lid			
MOPS-5221-92	MOPS, 100 G	In development	100.00 G		250.00 ML		White HDPE round pail, blue lid			
MOPS-5222-25	MOPS, 25 KG	Active	25.00 KG		12.00 GAL		White HDPE round pail, blue lid	30 x 48	LDPE	
NACL-3150-02	Sodium Chloride 5M Solution, 1 L	Active	1.00 L		1.00 L		Tepoflex sample			
NACL-3150-69	Sodium Chloride 5M Solution, 200 L	Active	200.00 L		58.00 GAL		Blue HDPE drum with metal ring for BPC		BPC bag	
NACL-3150-76	Sodium Chloride 5M Solution, 1025 L	Active	1025.00 L		1000.00 L		BPC shell	1000L	BPC bag	
NAHY-3150-76	Sodium Hydroxide 0.5N Solution, 1005 L	In development	1005.00 L		1000.00 L		BPC shell	1000L	BPC bag	
NAHY-3153-75	Sodium Hydroxide 1N Solution, 1005 L	In development	1005.00 L		275.00 GAL		Tote			
NAHY-3153-76	Sodium Hydroxide 1N Solution, 1005 L	In development	1005.00 L		1000.00 L		BPC shell	1000L	BPC bag	
NAHY-4122-70	Sodium Hydroxide 10N Solution, 1 L	Active	1.00 L		1.00 L		Non-fluorinated bottle, black lid			
NAHY-4122-72	Sodium Hydroxide 10N Solution, 200 L	Active	200.00 L		55.00 GAL		Natural HDPE drum			
NAHY-4122-75	Sodium Hydroxide 10N Solution, 940 L	Active	940.00 L		275.00 GAL		Tote			
NAHY-4122-96	Sodium Hydroxide 10N Solution, Bulk	Active	13059.00 L		13059.00 L		Bulk tanker			
NAHY-4123-70	Sodium Hydroxide 2N Solution, 1 L	Active	1.00 L		1.00 L		Non-fluorinated bottle, black lid			
NAHY-4123-71	Sodium Hydroxide 2N Solution, 1 L X 6 combo pack	Active	6.00 L		6.00 L		HDPE combo pack			
NAHY-4123-73	Sodium Hydroxide 2N Solution, 4 x 3.78 L	In development	15.12 L		15.12 L		4x1 Gallon combo pack			
NAHY-4124-78	Sodium Hydroxide 5N Solution, 10 L	Active	10.00 L		2.60 GAL		HDPE carboy			
NAHY-4126-71	Sodium Hydroxide 10N Solution Biotech, 6 X 1 L	Active	6.00 L		6.00 L		HDPE combo pack			
NAHY-4126-72	Sodium Hydroxide 10N Solution Biotech, 200 L	Active	200.00 L		55.00 GAL		Black HDPE drum			
NAHY-4126-73	Sodium Hydroxide 10N Solution Biotech, 4 X 3.78 L	Active	15.12 L		15.12 L		4x1 Gallon combo pack			
NAHY-4126-75	Sodium Hydroxide 10N Solution Biotech, 1000 L	Active	1000.00 L		275.00 GAL		Tote			
NAHY-4126-78	Sodium Hydroxide 10N Solution Biotech, 10 L	Active	10.00 L		2.60 GAL		HDPE carboy			
NAHY-4150-70	Sodium Hydroxide 10N Solution, 1 L	Active	1.00 L		1.00 L		Non-fluorinated bottle, black lid			
NAHY-4150-71	Sodium Hydroxide 10N Solution, 1 L X 6	Active	6.00 L		6.00 L		HDPE combo pack			
NAHY-4150-73	Sodium Hydroxide 10N Solution, 4 x 3.78 L	Active	15.12 L		15.12 L		4x1 Gallon combo pack			
NAHY-4150-78	Sodium Hydroxide 10N Solution, 10 L	Active	10.00 L		2.60 GAL		HDPE carboy			
NAHY-4150-79	Sodium Hydroxide 10N Solution, 942.5 L	Active	942.50 L		275.00 GAL		Tote			

Product	Product Description	Product Status	Product size	Product UOM	Package size	Package UOM	Package type	Liner size	Liner type	Desiccant
NAHY-4150-81	Sodium Hydroxide 10N Solution, 4 x 3.78 L	Active	15.12 L		15.12 L		4x1 Gallon combo pack			
NAHY-4150-83	Sodium Hydroxide 10N Solution, 4 x 4 L	Active	16.00 L		16.00 L		4 x 4L combo pack			
NAHY-4151-70	Sodium Hydroxide 8N Solution, 1 L	Active	1.00 L		1.00 L		Non-fluorinated bottle, black lid			
NAHY-7122-74	Sodium Hydroxide 10N Solution, 4 L	Active	4.00 L		4.00 L		HDPE combo pack			
NDEC-3201-01	Sodium Decanoate, 1kg	Active	1.00 KG		3.50 GAL		White HDPE round pail, blue lid	(2x) 24 x 36	LDPE	
NDEC-3201-03	Sodium Decanoate, 3 KG	Active	3.00 KG		3.50 GAL		White HDPE round pail, blue lid	24 x 36	LDPE	
NDEC-3201-05	Sodium Decanoate, 5 KG	Active	5.00 KG		5.00 GAL		White HDPE round pail, blue lid	24 x 36	LDPE	
NDEC-3201-10	Sodium Decanoate, 10 KG	Active	10.00 KG		12.00 GAL		White HDPE round pail, blue lid	30 x 48	LDPE	
NDEC-3201-25	Sodium Decanoate, 25 KG	Active	25.00 KG		30.00 GAL		Dark blue HDPE drum	34 x 60	LDPE	
NDEC-3201-92	Sodium Decanoate, 100 gr	Active	100.00 G		1000.00 ML		White HDPE bottle, blue lid			
NDEC-3201-94	Sodium Decanoate, 250 gr	Active	250.00 G		1500.00 ML		White HDPE bottle, blue lid			
THCL-3203-01	Tris Hydrochloride, 1 KG	Active	1.00 KG		1500.00 ML		White HDPE bottle, blue lid			
THCL-3203-05	Tris Hydrochloride, 5 KG	Active	5.00 KG		3.50 GAL		White HDPE round pail, blue lid	24 x 36	LDPE	
THCL-3203-10	Tris Hydrochloride, 10 KG	Active	10.00 KG		5.00 GAL		White HDPE round pail, blue lid	24 x 36	LDPE	
THCL-3203-25	Tris Hydrochloride, 25 KG	Active	25.00 KG		12.00 GAL		White HDPE round pail, blue lid	30 x 48	LDPE	
THCL-3203-92	Tris Hydrochloride, 100 G	Active	100.00 G		250.00 ML		White HDPE bottle, blue lid			
THCL-3203-93	Tris Hydrochloride, 500 G	Active	500.00 G		1000.00 ML		White HDPE bottle, blue lid			
THCL-3220-01	Tris Hydrochloride, 1 KG	Active	1.00 KG		1500.00 ML		White HDPE bottle, blue lid			
THCL-3220-05	Tris Hydrochloride, 5 KG	Active	5.00 KG		3.50 GAL		White HDPE round pail, blue lid	24 x 36	LDPE	
THCL-3220-10	Tris Hydrochloride, 10 KG	Active	10.00 KG		5.00 GAL		White HDPE round pail, blue lid	24 x 36	LDPE	
THCL-3220-12	Tris Hydrochloride, 12 KG	Active	12.00 KG		5.00 GAL		White HDPE round pail, blue lid	24 x 36	LDPE	
THCL-3220-13	Tris Hydrochloride, 12 KG	Active	12.00 KG		5.00 GAL		White HDPE round pail, blue lid	(2x) 30 x 48	LDPE	
THCL-3220-25	Tris Hydrochloride, 25 KG	Active	25.00 KG		12.00 GAL		White HDPE round pail, blue lid	30 x 48	LDPE	
THCL-3220-50	Tris Hydrochloride, 50 KG	Active	50.00 KG		30.00 GAL		Dark blue HDPE drum	(2x) 34 x 60	LDPE	4 Unit
THCL-3220-58	Tris Hydrochloride, 50 KG	Active	50.00 KG		30.00 GAL		Dark blue HDPE drum	(1x) 34 x 60, (1x) 3 LDPE/Tyvek		5 x 8 Unit
THCL-3220-92	Tris Hydrochloride, 100 G	Active	100.00 G		250.00 ML		White HDPE bottle, blue lid			
THCL-3220-93	Tris Hydrochloride, 500 G	Active	500.00 G		1000.00 ML		White HDPE bottle, blue lid			
THCL-3221-25	Tris Hydrochloride, 25 KG	Active	25.00 KG		12.00 GAL		White HDPE round pail, blue lid	30 x 48	LDPE	
THCL-3221-92	Tris Hydrochloride, 100 G	Active	100.00 G		250.00 ML		White HDPE bottle, blue lid			
THCL-3250-01	Tris Hydrochloride, 1 KG	Active	1.00 KG		1500.00 ML		White HDPE bottle, blue lid			
THCL-3250-05	Tris Hydrochloride, 5 KG	Active	5.00 KG		3.50 GAL		White HDPE round pail, blue lid	24 x 36	LDPE	
THCL-3250-10	Tris Hydrochloride, 10 KG	Active	10.00 KG		5.00 GAL		White HDPE round pail, blue lid	24 x 36	LDPE	
THCL-3250-12	Tris Hydrochloride, 12 KG	Active	12.00 KG		5.00 GAL		White HDPE round pail, blue lid	24 x 36	LDPE	
THCL-3250-13	Tris Hydrochloride, 12 KG	Active	12.00 KG		5.00 GAL		White HDPE round pail, blue lid	(2x) 30 x 48	LDPE	
THCL-3250-25	Tris Hydrochloride, 25 KG	Active	25.00 KG		12.00 GAL		White HDPE round pail, blue lid	(2x) 30 x 48	LDPE	
THCL-3250-28	Tris Hydrochloride, 25 KG	Active	25.00 KG		12.00 GAL		White HDPE round pail, blue lid	30 x 48	LDPE	
THCL-3250-41	Tris Hydrochloride, 10 KG	Active	10.00 KG		5.00 GAL		White HDPE round pail, blue lid	24 x 36	LDPE	
THCL-3250-42	Tris Hydrochloride, 12 KG	Active	12.00 KG		5.00 GAL		White HDPE round pail, blue lid	24 x 36	LDPE	
THCL-3250-58	Tris Hydrochloride, 50 KG	Active	50.00 KG		30.00 GAL		Dark blue HDPE drum	(1x) 34 x 60, (1x) 3 LDPE/Tyvek		5 x 8 Unit
THCL-3250-91	Tris Hydrochloride, 50 G	Active	50.00 G		250.00 ML		White HDPE bottle, blue lid			
THCL-3250-92	Tris Hydrochloride, 100 G	Active	100.00 G		250.00 ML		White HDPE bottle, blue lid			
THCL-3250-93	Tris Hydrochloride, 500 G	Active	500.00 G		1000.00 ML		White HDPE bottle, blue lid			
THCL-3253-10	Tris Hydrochloride, 10 KG	Active	10.00 KG		5.00 GAL		White HDPE round pail, blue lid	24 x 36	LDPE	
THCL-3253-12	Tris Hydrochloride, 12 KG	Active	12.00 KG		5.00 GAL		White HDPE round pail, blue lid	24 x 36	LDPE	
THCL-3253-13	Tris Hydrochloride, 12 KG	Active	12.00 KG		5.00 GAL		White HDPE round pail, blue lid	(2x) 30 x 48	LDPE	

Product	Product Description	Product Status	Product size	Product UOM	Package size	Package UOM	Package type	Liner size	Liner type	Desiccant
THCL-3253-92	Tris Hydrochloride, 100 G	Active	100.00 G		250.00 ML		White HDPE bottle, blue lid			
THCL-3253-93	Tris Hydrochloride, 500 G	Active	500.00 G		1000.00 ML		White HDPE bottle, blue lid			
THCL-3254-01	Tris Hydrochloride, 1 KG	Active	1.00 KG		1500.00 ML		White HDPE bottle, blue lid			
THCL-3254-12	Tris Hydrochloride, 12 KG	Active	12.00 KG		5.00 GAL		White HDPE round pail, blue lid	24 x 36	LDPE	
THCL-3254-13	Tris Hydrochloride, 12 KG	Active	12.00 KG		5.00 GAL		White HDPE round pail, blue lid	30 x 48	LDPE	
THCL-3254-25	Tris Hydrochloride, 25 KG	Active	25.00 KG		12.00 GAL		White HDPE round pail, blue lid	(2x) 30 x 48	LDPE	
THCL-3254-50	Tris Hydrochloride, 50 KG	Active	50.00 KG		30.00 GAL		Dark blue HDPE drum	(2x) 34 x 60	LDPE	4 Unit
THCL-3254-92	Tris Hydrochloride, 100 G	Active	100.00 G		250.00 ML		White HDPE bottle, blue lid			
THCL-3254-93	Tris Hydrochloride, 500 G	Active	500.00 G		1000.00 ML		White HDPE bottle, blue lid			
THCL-3257-10	Tris Hydrochloride, 10 KG	Active	10.00 KG		5.00 GAL		White HDPE round pail, blue lid	24 x 36	LDPE	
THCL-3257-93	Tris Hydrochloride, 500 G	Active	500.00 G		1000.00 ML		White HDPE bottle, blue lid			
THCL-3259-01	Tris Hydrochloride, 1 KG	Active	1.00 KG		1500.00 ML		White HDPE bottle, blue lid			
THCL-3259-05	Tris Hydrochloride, 5 KG	Active	5.00 KG		3.50 GAL		White HDPE round pail, blue lid	24 x 36	LDPE	
THCL-3259-10	Tris Hydrochloride, 10 KG	Active	10.00 KG		5.00 GAL		White HDPE round pail, blue lid	24 x 36	LDPE	
THCL-3259-12	Tris Hydrochloride, 12 KG	Active	12.00 KG		5.00 GAL		White HDPE round pail, blue lid	24 x 36	LDPE	
THCL-3259-25	Tris Hydrochloride, 25 KG	Active	25.00 KG		12.00 GAL		White HDPE round pail, blue lid	30 x 48	LDPE	
THCL-3259-58	Tris Hydrochloride, 50 KG	Active	50.00 KG		30.00 GAL		Dark blue HDPE drum	(1x) 34 x 60, (1x) 3 LDPE/Tyvek		5 x 8 Unit
THCL-3259-91	Tris Hydrochloride, 250 G	Active	250.00 G		1000.00 ML		White HDPE bottle, blue lid			
THCL-3259-92	Tris Hydrochloride, 100 G	Active	100.00 G		250.00 ML		White HDPE bottle, blue lid			
THCL-3259-93	Tris Hydrochloride, 500 G	Active	500.00 G		1000.00 ML		White HDPE bottle, blue lid			
THCL-3260-01	Tris Hydrochloride, 1 KG	Active	1.00 KG		1500.00 ML		White HDPE bottle, blue lid			
THCL-3260-92	Tris Hydrochloride, 100 G	Active	100.00 G		250.00 ML		White HDPE bottle, blue lid			
THCL-4251-01	Tris Hydrochloride Biotech, 1 KG	Active	1.00 KG		1500.00 ML		White HDPE bottle, blue lid			
THCL-4251-05	Tris Hydrochloride Biotech, 5 KG	Active	5.00 KG		3.50 GAL		White HDPE round pail, blue lid	24 x 36	LDPE	
THCL-4251-10	Tris Hydrochloride Biotech, 10 KG	Active	10.00 KG		5.00 GAL		White HDPE round pail, blue lid	24 x 36	LDPE	
THCL-4251-25	Tris Hydrochloride Biotech, 25 KG	Active	25.00 KG		12.00 GAL		White HDPE round pail, blue lid	30 x 48	LDPE	
THCL-4251-50	Tris Hydrochloride Biotech, 50 KG	In development	50.00 KG		30.00 GAL		Dark blue HDPE drum	(2x) 34 x 60	LDPE	4 Unit
THCL-4251-92	Tris Hydrochloride Biotech, 100 G	Active	100.00 G		250.00 ML		White HDPE bottle, blue lid			
THCL-4251-93	Tris Hydrochloride Biotech, 500 G	Active	500.00 G		1000.00 ML		White HDPE bottle, blue lid			
THCL-6201-01	Tris Hydrochloride, 1KG	Active	1.00 KG		1500.00 ML		White HDPE bottle, blue lid			
THCL-6201-92	Tris Hydrochloride, 100 G	In development	100.00 G		250.00 ML		White HDPE bottle, blue lid			
TRED-3250-03	Trehalose Dihydrate, 3 KG	Active	3.00 KG		3.50 GAL		White HDPE round pail, blue lid	24 x 36	LDPE	
TRED-3250-05	Trehalose Dihydrate, 5 KG	Active	5.00 KG		3.50 GAL		White HDPE round pail, blue lid	24 x 36	LDPE	
TRED-3250-06	Trehalose Dihydrate, 1 KG	Active	1.00 KG		0.60 GAL		White HDPE round pail, blue lid	(2x) 9 x 15	LDPE	
TRED-3250-10	Trehalose Dihydrate, 10 KG	Active	10.00 KG		5.00 GAL		White HDPE round pail, blue lid	24 x 36	LDPE	
TRED-3250-25	Trehalose Dihydrate, 25 KG	Active	25.00 KG		12.00 GAL		White HDPE round pail, blue lid	30 x 48	LDPE	
TRED-3250-50	Trehalose Dihydrate, 50 KG	Active	50.00 KG		30.00 GAL		Dark blue HDPE drum	34 x 60	LDPE	
TRED-3250-92	Trehalose Dihydrate, 100g	Active	100.00 G		250.00 ML		White HDPE bottle, blue lid			
TRED-3252-06	Trehalose Dihydrate, 1 KG	Active	1.00 KG		0.60 GAL		White HDPE round pail, blue lid	(2x) 9 x 15	LDPE	
TRED-3252-10	Trehalose Dihydrate, 10 KG	Active	10.00 KG		5.00 GAL		White HDPE round pail, blue lid	24 x 36	LDPE	
TRED-3252-25	Trehalose Dihydrate, 25 KG	Active	25.00 KG		12.00 GAL		White HDPE round pail, blue lid	30 x 48	LDPE	
TRED-3252-92	Trehalose Dihydrate, 100g	Active	100.00 G		250.00 ML		White HDPE bottle, blue lid			
TRED-3253-06	Trehalose Dihydrate, 1 KG	In development	1.00 KG		0.60 GAL		White HDPE round pail, blue lid	(2x) 9 x 15	LDPE	
TRED-3253-10	Trehalose Dihydrate, 10 KG	Active	10.00 KG		5.00 GAL		White HDPE round pail, blue lid	24 x 36	LDPE	
TRED-4250-05	Trehalose Dihydrate Biotech, 5 KG	Active	5.00 KG		3.50 GAL		White HDPE round pail, blue lid	24 x 36	LDPE	

Product	Product Description	Product Status	Product size	Product UOM	Package size	Package UOM	Package type	Liner size	Liner type	Desiccant
TRED-4250-06	Trehalose Dihydrate Biotech, 1 KG	Active	1.00 KG		0.60 GAL		White HDPE round pail, blue lid	(2x) 9 x 15	LDPE	
TRED-4250-10	Trehalose Dihydrate Biotech, 10 KG	Active	10.00 KG		5.00 GAL		White HDPE round pail, blue lid	24 x 36	LDPE	
TRED-4250-25	Trehalose Dihydrate Biotech, 25 KG	Active	25.00 KG		12.00 GAL		White HDPE round pail, blue lid	30 x 48	LDPE	
TRED-4250-50	Trehalose Dihydrate Biotech, 50 KG	Active	50.00 KG		30.00 GAL		Dark blue HDPE drum	34 x 60	LDPE	
TRIS-1201-01	Tris/Tromethamine, 1 KG	Active	1.00 KG		1500.00 ML		White HDPE bottle, blue lid			
TRIS-1201-10	Tris/Tromethamine, 10 KG	Active	10.00 KG		5.00 GAL		White HDPE round pail, blue lid	24 x 36	LDPE	
TRIS-1201-25	Tris/Tromethamine, 25 KG	Active	25.00 KG		12.00 GAL		White HDPE round pail, blue lid	30 x 48	LDPE	
TRIS-1201-91	Tris/Tromethamine, 250 G	Active	250.00 G		1000.00 ML		White HDPE bottle, blue lid			
TRIS-1201-92	Tris/Tromethamine, 100 G	Active	100.00 G		250.00 ML		White HDPE bottle, blue lid			
TRIS-1201-93	Tris/Tromethamine, 500 G	Active	500.00 G		1000.00 ML		White HDPE bottle, blue lid			
TRIS-1220-01	Tris/Tromethamine, 1 KG	Active	1.00 KG		1500.00 ML		White HDPE bottle, blue lid			
TRIS-1220-25	Tris/Tromethamine, 25 KG	Active	25.00 KG		12.00 GAL		White HDPE round pail, blue lid	30 x 48	LDPE	
TRIS-1220-91	Tris/Tromethamine, 250 G	Active	250.00 G		1000.00 ML		White HDPE bottle, blue lid			
TRIS-1220-92	Tris/Tromethamine, 100 G	Active	100.00 G		250.00 ML		White HDPE bottle, blue lid			
TRIS-1220-93	Tris/Tromethamine, 500 G	Active	500.00 G		1000.00 ML		White HDPE bottle, blue lid			
TRIS-2255-01	Tris/Tromethamine, 1 KG	Active	1.00 KG		1500.00 ML		White HDPE bottle, blue lid			
TRIS-2255-10	Tris/Tromethamine, 10 KG	Active	10.00 KG		5.00 GAL		White HDPE round pail, blue lid	24 x 36	LDPE	
TRIS-2255-91	Tris/Tromethamine, 250 G	Active	250.00 G		1000.00 ML		White HDPE bottle, blue lid			
TRIS-2255-92	Tris/Tromethamine, 100 G	Active	100.00 G		250.00 ML		White HDPE bottle, blue lid			
TRIS-2255-93	Tris/Tromethamine, 500 G	Active	500.00 G		1000.00 ML		White HDPE bottle, blue lid			
TRIS-2257-01	Tris/Tromethamine, 1 KG	Active	1.00 KG		1500.00 ML		White HDPE bottle, blue lid			
TRIS-2257-25	Tris/Tromethamine, 25 KG	Active	25.00 KG		12.00 GAL		White HDPE round pail, blue lid	30 x 48	LDPE	
TRIS-2257-91	Tris/Tromethamine, 250 G	Active	250.00 G		1000.00 ML		White HDPE bottle, blue lid			
TRIS-2257-92	Tris/Tromethamine, 100 G	Active	100.00 G		250.00 ML		White HDPE bottle, blue lid			
TRIS-2257-93	Tris/Tromethamine, 500 G	Active	500.00 G		1000.00 ML		White HDPE bottle, blue lid			
TRIS-3201-01	Tris/Tromethamine, 1 KG	Active	1.00 KG		1500.00 ML		White HDPE bottle, blue lid			
TRIS-3201-92	Tris/Tromethamine, 100 G	Active	100.00 G		250.00 ML		White HDPE bottle, blue lid			
TRIS-3201-93	Tris/Tromethamine, 500 G	Active	500.00 G		1000.00 ML		White HDPE bottle, blue lid			
TRIS-3220-01	Tris/Tromethamine, 1 KG	Active	1.00 KG		1500.00 ML		White HDPE bottle, blue lid			
TRIS-3220-10	Tris/Tromethamine, 10 KG	Active	10.00 KG		5.00 GAL		White HDPE round pail, blue lid	24 x 36	LDPE	
TRIS-3220-25	Tris/Tromethamine, 25 KG	Active	25.00 KG		12.00 GAL		White HDPE round pail, blue lid	30 x 48	LDPE	
TRIS-3220-50	Tris/Tromethamine, 50 KG	Active	50.00 KG		30.00 GAL		Dark blue HDPE drum	(2x) 34 x 60	LDPE	4 Unit
TRIS-3220-92	Tris/Tromethamine, 100 G	Active	100.00 G		250.00 ML		White HDPE bottle, blue lid			
TRIS-3220-93	Tris/Tromethamine, 500 G	Active	500.00 G		1000.00 ML		White HDPE bottle, blue lid			
TRIS-3251-01	Tris/Tromethamine, 1 KG	Active	1.00 KG		1500.00 ML		White HDPE bottle, blue lid			
TRIS-3251-05	Tris/Tromethamine, 5 KG	Active	5.00 KG		3.50 GAL		White HDPE round pail, blue lid	24 x 36	LDPE	
TRIS-3251-10	Tris/Tromethamine, 10 KG	Active	10.00 KG		5.00 GAL		White HDPE round pail, blue lid	24 x 36	LDPE	
TRIS-3251-92	Tris/Tromethamine, 100 G	Active	100.00 G		250.00 ML		White HDPE bottle, blue lid			
TRIS-3251-95	Tris/Tromethamine, 80 G	Active	80.00 G		1000.00 ML		White HDPE bottle, blue lid	9 x 15	LDPE	
TRIS-3252-25	Tris/Tromethamine, 25 KG	Active	25.00 KG		12.00 GAL		White HDPE round pail, blue lid	30 x 48	LDPE	
TRIS-3252-60	Tris/Tromethamine, 50 KG	Active	50.00 KG		30.00 GAL		Dark blue HDPE drum	(1x) 34 x 60, (2x) 3 LDPE/(2x) Tyvek		(2x) 5 x 8 Unit
TRIS-3254-01	Tris/Tromethamine, 1KG	Active	1.00 KG		1500.00 ML		White HDPE bottle, blue lid			
TRIS-3254-05	Tris/Tromethamine, 5 KG	Active	5.00 KG		3.50 GAL		White HDPE round pail, blue lid	24 x 36	LDPE	
TRIS-3254-10	Tris/Tromethamine, 10 KG	Active	10.00 KG		5.00 GAL		White HDPE round pail, blue lid	24 x 36	LDPE	
TRIS-3254-25	Tris/Tromethamine, 25 KG	Active	25.00 KG		12.00 GAL		White HDPE round pail, blue lid	30 x 48	LDPE	

Product	Product Description	Product Status	Product size	Product UOM	Package size	Package UOM	Package type	Liner size	Liner type	Desiccant
TRIS-3254-50	Tris/Tromethamine, 50 KG	Active	50.00 KG		30.00 GAL		Dark blue HDPE drum	(2x) 34 x 60	LDPE	4 Unit
TRIS-3254-58	Tris/Tromethamine, 50 KG	Active	50.00 KG		30.00 GAL		Dark blue HDPE drum	(1x) 34 x 60, (1x) 3 LDPE/Tyvek		5 x 8 Unit
TRIS-3254-92	Tris/Tromethamine, 100 G	Active	100.00 G		250.00 ML		White HDPE bottle, blue lid			
TRIS-3254-93	Tris/Tromethamine, 500 G	Active	500.00 G		1000.00 ML		White HDPE bottle, blue lid			
TRIS-3255-01	Tris/Tromethamine, 1 KG	Active	1.00 KG		1500.00 ML		White HDPE bottle, blue lid			
TRIS-3255-05	Tris/Tromethamine, 5 KG	Active	5.00 KG		3.50 GAL		White HDPE round pail, blue lid	24 x 36	LDPE	
TRIS-3255-07	Tris/Tromethamine, 2 KG	Active	2.00 KG		1.25 GAL		White HDPE round pail, blue lid	24 x 36	LDPE	
TRIS-3255-10	Tris/Tromethamine, 10 KG	Active	10.00 KG		5.00 GAL		White HDPE round pail, blue lid	24 x 36	LDPE	
TRIS-3255-12	Tris/Tromethamine, 12 KG	Active	12.00 KG		5.00 GAL		White HDPE round pail, blue lid	24 x 36	LDPE	
TRIS-3255-25	Tris/Tromethamine, 25 KG	Active	25.00 KG		12.00 GAL		White HDPE round pail, blue lid	30 x 48	LDPE	
TRIS-3255-41	Tris/Tromethamine, 10 KG	Active	10.00 KG		5.00 GAL		White HDPE round pail, blue lid	24 x 36	LDPE	
TRIS-3255-42	Tris/Tromethamine, 12 KG	Active	12.00 KG		5.00 GAL		White HDPE round pail, blue lid	24 x 36	LDPE	
TRIS-3255-50	Tris/Tromethamine, 50 KG	Active	50.00 KG		30.00 GAL		Dark blue HDPE drum	(2x) 34 x 60	LDPE	4 Unit
TRIS-3255-58	Tris/Tromethamine, 50 KG	Active	50.00 KG		30.00 GAL		Dark blue HDPE drum	(1x) 34 x 60, (1x) 3 LDPE/Tyvek		5 x 8 Unit
TRIS-3255-89	Tris/Tromethamine, 200 G	Active	200.00 G		1000.00 ML		White HDPE bottle, blue lid			
TRIS-3255-91	Tris/Tromethamine, 50 G	Active	50.00 G		250.00 ML		White HDPE bottle, blue lid			
TRIS-3255-92	Tris/Tromethamine, 100 G	Active	100.00 G		250.00 ML		White HDPE bottle, blue lid			
TRIS-3255-93	Tris/Tromethamine, 500 G	Active	500.00 G		1000.00 ML		White HDPE bottle, blue lid			
TRIS-3255-94	Tris/Tromethamine, 250 G	Active	250.00 G		1000.00 ML		White HDPE bottle, blue lid			
TRIS-3257-10	Tris/Tromethamine, 10 KG	Active	10.00 KG		5.00 GAL		White HDPE round pail, blue lid	24 x 36	LDPE	
TRIS-4223-01	Tris (Tromethamine) Biotech, 1 KG	Active	1.00 KG		1500.00 ML		White HDPE bottle, blue lid			
TRIS-4223-05	Tris (Tromethamine) Biotech, 5 KG	Active	5.00 KG		3.50 GAL		White HDPE round pail, blue lid	24 x 36	LDPE	
TRIS-4223-10	Tris (Tromethamine) Biotech, 10 KG	Active	10.00 KG		5.00 GAL		White HDPE round pail, blue lid	24 x 36	LDPE	
TRIS-4223-25	Tris (Tromethamine) Biotech, 25 KG	Active	25.00 KG		12.00 GAL		White HDPE round pail, blue lid	30 x 48	LDPE	
TRIS-4223-50	Tris (Tromethamine) Biotech, 50 KG	In development	50.00 KG		30.00 GAL		Dark blue HDPE drum	(2x) 34 x 60	LDPE	4 Unit
TRIS-4223-92	Tris (Tromethamine) Biotech, 100 G	Active	100.00 G		250.00 ML		White HDPE bottle, blue lid			
TRIS-4223-93	Tris (Tromethamine) Biotech, 500 G	Active	500.00 G		1000.00 ML		White HDPE bottle, blue lid			
TRIS-4224-92	Tris (Tromethamine), 100 G	Active	100.00 G		250.00 ML		White HDPE bottle, blue lid			
TRIS-5201-26	Tris/Tromethamine, 25 KG	Active	25.00 KG		26.00 GAL		Fiber drum	(2x) 34 x 60	LDPE	4 Unit
TRIS-5203-59	Tris/Tromethamine, 50 KG	Active	50.00 KG		26.00 GAL		Fiber drum	(2x) 34 x 60	LDPE	4 Unit
TRIS-5204-26	Tris/Tromethamine, 25 KG	Active	25.00 KG		26.00 GAL		Fiber drum	(2x) 34 x 60	LDPE	
TRIS-5207-26	Tris/Tromethamine, 25 KG	Active	25.00 KG		26.00 GAL		Fiber drum	(2x) 34 x 60	LDPE	4 Unit
TRIS-5220-10	Tris/Tromethamine, 10 KG	Active	10.00 KG		5.00 GAL		White HDPE round pail, blue lid	24 x 36	LDPE	
TRIS-5222-10	Tris/Tromethamine, 10 KG	Active	10.00 KG		5.00 GAL		White HDPE round pail, blue lid	24 x 36	LDPE	
TRIS-5223-89	Tris/Tromethamine, 200 G	Active	200.00 G		1000.00 ML		White HDPE bottle, blue lid			
TRIS-6201-01	Tris/Tromethamine, 1 KG	Active	1.00 KG		1500.00 ML		White HDPE bottle, blue lid			
TRIS-6201-92	Tris/Tromethamine, 100 G	Active	100.00 G		250.00 ML		White HDPE bottle, blue lid			
URAC-4201-05	Uracil, 5 KG	Active	5.00 KG		3.50 GAL		White HDPE round pail, blue lid	24 x 36	LDPE	
URAC-4202-01	Uracil, 1 KG	Active	1.00 KG		1500.00 ML		White HDPE bottle, blue lid			
URAC-4202-05	Uracil, 5 KG	Active	5.00 KG		3.50 GAL		White HDPE round pail, blue lid	24 x 36	LDPE	
URAC-4202-91	Uracil, 50 G	Active	50.00 G		250.00 ML		White HDPE bottle, blue lid			
URAC-4202-92	Uracil, 100 G	Active	100.00 G		250.00 ML		White HDPE bottle, blue lid			
URAC-4202-95	Uracil, 10g	Active	10.00 G		60.00 ML		White polypropylene bottle			
URAC-4250-01	Uracil, 1 KG	Active	1.00 KG		1500.00 ML		White HDPE bottle, blue lid			
URAC-4250-95	Uracil, 10g	Active	10.00 G		60.00 ML		White polypropylene bottle			

Product	Product Description	Product Status	Product size	Product UOM	Package size	Package UOM	Package type	Liner size	Liner type	Desiccant
URAC-4301-05	Uracil, 5 KG	Active	5.00 KG		3.50 GAL		White HDPE round pail, blue lid	(2x) 24 x 36	LDPE	
URAC-4301-06	Uracil, 5 KG	Active	5.00 KG		3.50 GAL		White HDPE round pail, blue lid	(5x) 24 x 36	LDPE	
UREA-3120-02	Urea 6M Solution, 500 ML	Active	500.00 ML		17.00 OZ		Black HDPE UN leakproof bottle			
UREA-3120-76	Urea 6M Solution, 750 L	In development	750.00 L		1000.00 L		BPC shell	1000L	BPC bag	
UREA-3120-78	Urea 6M Solution, 1200 L	Active	1200.00 L		330.00 GAL		Tote		LDPE	
UREA-3220-10	Urea, 10 KG	Active	10.00 KG		5.00 GAL		White HDPE round pail, blue lid	24 x 36	LDPE	
UREA-3220-25	Urea, 25 KG	Active	25.00 KG		12.00 GAL		White HDPE round pail, blue lid	34 x 60	LDPE	
UREA-3220-50	Urea, 50 KG	Active	50.00 KG		30.00 GAL		Dark blue HDPE drum	(1x) 34 x 60, (1x) 3 LDPE/Tyvek		5 x 8 Unit
UREA-3220-89	Urea, 200g	Active	200.00 G		1000.00 ML		White HDPE bottle, blue lid			
UREA-3220-92	Urea, 100g	Active	100.00 G		250.00 ML		White HDPE bottle, blue lid			
UREA-3220-93	Urea, 500 G	Active	500.00 G		1000.00 ML		White HDPE bottle, blue lid			
UREA-3222-10	Urea, 10 KG	Active	10.00 KG		5.00 GAL		White HDPE round pail, blue lid	24 x 36	LDPE	
UREA-3222-25	Urea, 25 KG	Active	25.00 KG		12.00 GAL		White HDPE round pail, blue lid	34 x 60	LDPE	
UREA-3223-25	Urea, 25 KG	Active	25.00 KG		12.00 GAL		White HDPE round pail, blue lid	34 x 60	LDPE	
UREA-3223-58	Urea, 50 KG	Active	50.00 KG		30.00 GAL		Dark blue HDPE drum	(1x) 34 x 60, (1x) 3 LDPE/Tyvek		5 x 8 Unit
UREA-3250-01	Urea, 1KG	Active	1.00 KG		1500.00 ML		White HDPE bottle, blue lid			
UREA-3250-12	Urea, 12 KG	Active	12.00 KG		5.00 GAL		White HDPE round pail, blue lid	24 x 36	LDPE	4 Unit
UREA-3250-25	Urea, 25 KG	Active	25.00 KG		12.00 GAL		White HDPE round pail, blue lid	34 x 60	LDPE	
UREA-3250-92	Urea, 100g	Active	100.00 G		250.00 ML		White HDPE bottle, blue lid			
UREA-3250-93	Urea, 500 G	Active	500.00 G		1000.00 ML		White HDPE bottle, blue lid			
UREA-4221-01	Urea Biotech, 1 KG	Active	1.00 KG		1500.00 ML		White HDPE bottle, blue lid			
UREA-4221-05	Urea Biotech, 5 KG	Active	5.00 KG		3.50 GAL		White HDPE round pail, blue lid	24 x 36	LDPE	
UREA-4221-10	Urea Biotech, 10 KG	Active	10.00 KG		5.00 GAL		White HDPE round pail, blue lid	24 x 36	LDPE	
UREA-4221-25	Urea Biotech, 25 KG	Active	25.00 KG		12.00 GAL		White HDPE round pail, blue lid	30 x 48	LDPE	
UREA-4221-50	Urea Biotech, 50 KG	Active	50.00 KG		30.00 GAL		Dark blue HDPE drum	34 x 60	LDPE	
UREA-4221-92	Urea Biotech, 100 G	Active	100.00 G		250.00 ML		White HDPE bottle, blue lid			
UREA-4221-93	Urea Biotech, 500 G	Active	500.00 G		1000.00 ML		White HDPE bottle, blue lid			
UREA-5203-05	Urea, 5 KG	Active	5.00 KG		3.50 GAL		White HDPE round pail, blue lid	24 x 36	LDPE	
UREA-5203-15	Urea, 15 KG	Active	15.00 KG		12.00 GAL		White HDPE round pail, blue lid	34 x 60	LDPE	
UREA-5220-25	Urea, 25 KG	Active	25.00 KG		12.00 GAL		White HDPE round pail, blue lid	30 x 48	LDPE	
URID-3250-01	Uridine, 1 KG	Active	1.00 KG		1500.00 ML		White HDPE bottle, blue lid			
URID-3250-10	Uridine, 10 KG	Active	10.00 KG		5.00 GAL		White HDPE round pail, blue lid	24 x 36	LDPE	
URID-3250-25	Uridine, 25 KG	Active	25.00 KG		12.00 GAL		White HDPE round pail, blue lid	30 x 48	LDPE	
URID-3250-91	Uridine, 250 G	Active	250.00 G		1000.00 ML		White HDPE bottle, blue lid			
URID-3250-92	Uridine, 100 G	Active	100.00 G		250.00 ML		White HDPE bottle, blue lid			
URID-3250-93	Uridine, 500 G	Active	500.00 G		1000.00 ML		White HDPE bottle, blue lid			
URID-4250-01	Uridine Biotech, 1 KG	Active	1.00 KG		1500.00 ML		White HDPE bottle, blue lid			
URID-4250-05	Uridine Biotech, 5KG	Active	5.00 KG		3.50 GAL		White HDPE round pail, blue lid	24 x 36	LDPE	
URID-4250-10	Uridine Biotech, 10 KG	Active	10.00 KG		5.00 GAL		White HDPE round pail, blue lid	24 x 36	LDPE	
URID-4250-25	Uridine Biotech, 25 KG	Active	25.00 KG		12.00 GAL		White HDPE round pail, blue lid	30 x 48	LDPE	
URID-4250-50	Uridine Biotech, 50 KG	Active	50.00 KG		30.00 GAL		Dark blue HDPE drum	34 x 60	LDPE	
URID-4250-92	Uridine Biotech, 100 G	Active	100.00 G		250.00 ML		White HDPE bottle, blue lid			
URID-4250-93	Uridine Biotech, 500 G	Active	500.00 G		1000.00 ML		White HDPE bottle, blue lid			
URID-4250-95	Uridine Biotech, 50 G	Active	100.00 G		250.00 ML		White HDPE bottle, blue lid			
WAFI-3150-68	Water for Injection (WFI), 200L	In development	200.00 L		200.00 L		BPC shell	200L	Typoflex BPC bag	

Product	Product Description	Product Status	Product size	Product UOM	Package size	Package UOM	Package type	Liner size	Liner type	Desiccant
WAFI-3150-74	Water for Injection (WFI), 3.78L	Active	3.78 L		3.78 L		Non-fluorinated HDPE bottle			