

19L CARBOY SPECIFICATIONS



5 Gallon / 19L HDPE Carboy	
Test	Specifications
Capacity	5 Gallons
Material	Natural HDPE
Appearance and Color	White, rectangular with built-in handle, vent lid and spigot adapter
Outside Length (Internal Tolerance)	10.00 inches \pm 0.5 inch
Outside Width (Internal Tolerance)	9.75 inches \pm 0.5 inch
Outside Height (Internal Tolerance)	15 inches \pm 0.5 inch including lid
Mouth Opening (Internal Tolerance)	2 and 3/8 inches \pm 0.5 inches

For further information, please contact marketing@biospectra.us.



**Connie Augustine | Director
Business Integration & Retention
BioSpectra Inc.**

This is not considered a controlled document. We are not responsible for any errors or omissions, and the user is responsible for any decisions based on the information herein.